

Product datasheet for **RG206184**

KHDC1 (NM_030568) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: KHDC1 (NM_030568) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: KHDC1
Synonyms: bA257K9.4; C6orf147; C6orf148; Em:AC019205.8; KHDC1L; NDG1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG206184 representing NM_030568
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGACATGGGAACGAGTGCTCTCAGCAAGAAGCCGTGGTGGACCCTGCCTCAAACCTTTCATGCACCAA
 TGGTGTTCACATGGAAGAGGACCAGGAGGAGCTCATCTTCGGGCATGGTGACACATACCTTCGCTGCAT
 TGAGGTGCACAGCCACACCCTTATTCAGCTGGAGAGTTGGTTCACAGCTACAGGCCAGACTCGTGTGACT
 GTTGGACCACACAGGCAAGGCAGTGGCTGCTGCACATGTTTTGTGTGGGGAGCCAGGACTCCTATC
 ATCATGCTCGAGGCTGGAGATGCTGGAACGTGTCCGAAGCCAGCCTCTGACCAACGATGACCTGGTCAC
 CTCCATTAGTGTGCCACCATACTGGAGACTGTCTTTGGCTCCAGGATAAGTGAACCGTGTGTCTT
 AGCGTTCTCAGCCTTCCCTTACCAAGTGATTGGTTGTCGGGTTTCCATTTGAGTTCACTGTATCCG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG206184 representing NM_030568
 Red=Cloning site Green=Tags(s)

MDMGTSALSKKPWWTL PQNFHAPMV FHMEDQEELIFGHGDTYLRCEVHSHTLIQLESWFTATGQTRVT
 VGPHRARQWLLHMFCCVGSQDSYHHARGLEMLERVRSQPLTNDLVTLSISVPPYTGDLSLAPRISGTVCL
 SVPQPSYQVIGCSGFHLSLYP

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_030568

ORF Size: 489 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030568.3](#), [NP_085045.3](#)

RefSeq Size: 1147 bp

RefSeq ORF: 495 bp


Locus ID: 80759

UniProt ID: [Q4VXA5](#)

Cytogenetics: 6q13

Protein Families: Transmembrane

Product images:



Circular map for RG206184