

Product datasheet for **RG206202**

C10orf82 (NM_144661) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C10orf82 (NM_144661) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: C10orf82
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG206202 representing NM_144661
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGCCTTCCAAGACCTTCATGAGAACTGCCAATCACACCAGGCTATAGCGGCTTTGTGCCATTCC
 TCAGCTGCCAAGGAATGTCCAAGGAGGATGACATGAACCACTGTGTGAAAACCTTCCAGGAGAAAACACA
 GCGCTATAAAGAACAGCTGCGGGAATTGTGCTGCGCAGTGGCCACTGCCCGAACTGAAACCTGTCAAC
 TCCGAGGAGACGGTCTGCAGGCCCTGCACCAGTACAATCTGCAGTACCACCCCTGATCCTGGAATGCA
 AATATGTAAGAAACCTCTCCAGGAGCCCCGATCCCTGGCTGGCAGGCTACCTGCCGAGCCAAAGGT
 CACTGAATTTGGCTGTGGCACGAGATACACTGTCATGGCCAAAACCTGCTACAAGGACTTCTGGAGATC
 ACGGAGAGGGCCAAGAAGGCACATCTGAAACCATGAAGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG206202 representing NM_144661
 Red=Cloning site Green=Tags(s)

MEPSKTFMRNLPITPGYSGFVFPFLSCQGMKEDDMNHCVKTFQEKTRRYKEQLRELCCAATAPKLPVN
 SEETVLQALHQYNLQYHPLILECKYVKKPLQEPPIPGWAGYLPRAKVTEFGCGTRYVMKNCYKDFLEI
 TERAHKAHLKPYEE

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_144661

ORF Size: 462 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_144661.4](#)

RefSeq Size: 1009 bp


RefSeq ORF: 465 bp

Locus ID: 143379

UniProt ID: [Q8WW14](#)

Cytogenetics: 10q25.3

Product images:



Circular map for RG206202