

## Product datasheet for **RG206220**

### **TNRC6B (NM\_001024843) Human Tagged ORF Clone**

#### **Product data:**

|                                  |                         |
|----------------------------------|-------------------------|
| <b>Product Type:</b>             | Expression Plasmids     |
| <b>Tag:</b>                      | TurboGFP                |
| <b>Symbol:</b>                   | TNRC6B                  |
| <b>Mammalian Cell Selection:</b> | Neomycin                |
| <b>Vector:</b>                   | pCMV6-AC-GFP (PS100010) |
| <b>E. coli Selection:</b>        | Ampicillin (100 ug/mL)  |



**ORF Nucleotide Sequence:** >RG206220 representing NM\_001024843  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGCAACCAATGATGAGGGAAGATATCGGAAGAAAGCAGTGTCCAAGGTGGAAACAGGAGGATTTGTATGTG  
AAGGGCATGGCAAGACTCCACCTCCTGGTGAAGAAAGCAAACAAGAGAAGGAGCAAGAAAGGGAAGAACA  
GTTAATGGAAGACAAGAAAAGGAAGAAAGAGGATAAAAAAGAAAAAGAGCCACTCAGAAGGTACGCGAA  
CAAAAAACCAAGTGCCCGAAGTGACGAAACCAAGTTTAAGCCAACCAACGGCCGCCAGCCAATTGGCA  
GCTCTCCATCGCCACCAGTCAATGGTGGCAACAATGCCAAAAGGGTGGCAGTGCCGAACGGACAACCGCC  
AAGCGCGCCGCTACATGCCTCGGAGAGTGCCGCCGCGATTCCGTTGCCAGCAGGACCACAAAGTTGA  
CTAAAACGTTGGGCAAGCCCTCCACCTGCTCGTACGCTTCTGGGGTGGGCGAGGCCCTCTCCCTGCA  
CAGCACTGTGAGCAAAACCAACACGCACAAGTGACGAGGCGCTGCTGCAAGTGAAGTGGGAGTGC  
GCCAGTCTGGAGCAAAAGCACACCACCTGCTCCAGATAATGGTACTTCCGCTTGGGGTGAAGCAATGAA  
AGCAGTCTGGTGGGCGAGATGGATGATACGGGAGCATCGACCACAGGCTGGGGGAACAGCCCGCCA  
ACGCTCCCAATGCCATGAAGCCTAATTCCAATCTATGCAAGACGGCTGGGGGAGAGTGACGGGCCAGT  
CACAGGAGCTCGCCATCCAGCTGGGAAGAGGAGGAGGATGGAGGAGTCTGGAACACCACTGGCTCTCAG  
GGCAGTGCTTCTCCACAACCTCAGCAAGCTGGGGAACAGGAGGAAGAAACAAATGAAGTGCTCACTCA  
AAGGAGGAACAATGATTCATGGATGAATCCTCTGCCAAACAGTTTTCAAATATGGGATGCTGAGTCA  
GACTGAAGATAATCCAAGCAGCAAAATGGATTTGTCTGTAGGAAGCCTTTCAGATAAAAAATTTGATGTG  
GACAAGCGAGCGATGAATCTCGGGGATTTTAAATGATATCATGAGGAAGGATCGATCTGGGTTCCGTCAC  
CTAATTCCAAGACATGGGAACCACAGATAGTGGGCCTTATTTGAGAAGGGCGGTAGTCATGGTTTGT  
TGAAACAGCACAGCACAATCGAGAGGTCTGCACACACCCGTGCAGCCACTAAATCTTCTCCAGTCTC  
CGGGCGCAAGTGCCCTCCCAAGTTTATTTCCCCCAAGTTTCTGCCTCAATGCTCAAGCAGTTTCCAACA  
GTGGCCTGAGTCCAGGTCTTTTCAATGTGGGGCCCCAGTTATCTCTCAACAAATTGCCATGCTGAGCCA  
GCTTCCACAAATTCCCAAGTTTTCAGTTGGCATGTCAGCTTCTCTTGACGACGACGAACAGCAGCAGTTG  
TACAGAAACAGAGAAGATTCTCAAGCTGATCGCCAACAGCAAGACAGCAGCTGGCTCGAATGGTGA  
GTGCTACGACGACGACGACGACGACGACGACGACGACGACGACGACGACGACGACGACGACGACGACG  
TGGGCCCAAGCCGCATCTGGACAACATGGTACCAACGCAATTGAATGTGGGGCTCCAGACCTTCAAACC  
AAAGGGCAATACCTGGATATGGTTCTGGCTTTCAGCTCTGGCGGCATGGACTATGGCATGGTGGTGGGA  
AGGAGGCTGGAACCGAGTCTCGCTTTAAACAGTGGACCTCCATGATGGAGGGCTGCCCTCTGTAGCCAC  
ACAGGAAGCCAATATGCACAAAATGGCGCTATAGTGGCCCTGGTAAACCCGGGAGGGTCAACCGTAC  
AACCAGTTTGATATCATCCCTGGTGACACATGGGTGGCCATACGGGTCTGCTGGTGATAGCTGGTTAC  
CTGCCAAATCTCCACCAAAAATAAAATCGGAAGTAAATCCAGCAATGCCAGTTGGCTCCGAAGATCCA  
ACAGGAGTGCCATGAAAGGATATCCAAAACATTGACCTGAATGCCCTTATGTCACCCAGGAAGT  
GTGCTGGGGGTACAGCCACATCTCCATTGTAGATACTGACCACCAACTGCTGCGGGATAACACCACAG  
GGTCTAATTCTCCCTCAACACCTCGCTGCCTTCACCTGGTGCCTGGCCCTACAGTGCCTCTGACAACTC  
CTTTACCAACGTTTCATAGCACTTCAGCAAAGTTCCCTGATTACAAATCAACATGGTCCCAGATCCATA  
GGACACAACCCCACTCATCTCTCCAACAAGATGTGAAAAAACCATATTTCTCCAGGAACACTACACCGC  
TGCCCGCCCACTCTGGTCTGACCAACCCCAACCATCATCTCCCTGGAGCAGCAGCAGACCCCGATC  
AGTCAGGGGGTGGGGACACGAGGATACGGCTCGCCTCGGCCCTCACTTGGAGTGATGGTGGCTCAGTT  
CGTCTAGTTACTGGCTGGTTTCTCAACATCTACCCCAAGCATGATGGGTCAACCTTGAAGAACGATCT  
GCATGCAAGATGCGCCACTGCTGACATTCCATCTGAATCTAACCCAGGCGATGCCCTGATCCGATACAG  
CACCAAACAGGAGGCGGCCAAGGCCCAAACTGCACTGCACATGTGTGTGTTGGGAAACACTACCATCCTT  
GCTGAGTTTGCCACTGATGATGAAGTCAGCCGCTTTCTGGCACAAGCTCAGCCCCCTACACCTGCAGCAA  
CCCCAAGTGCGCCAGCTGCGGGGTGGCAGTCGCTGGAGACCGGCCAGAACCAGTCAGATCCCGTGGGACC  
TGCTCTGAATCTTTTTGGTGGGTCCACTGGGCTCGGGCAGTGGAGCAGCAGCGCTGGTGGCAGACGGGG  
GCCGATCTTGTGGCGCTTCATTGTGGGGGCCCCAACTATTCTTCTAGCTTATGGGGAGTCCCAACGG  
TGAAGATCCCATAGGATGGGCAGCCCTGCTCCTTTACTACCTGGTGACCTTCTGGGAGGAGGGTCGGA  
TTCAATC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG206220 representing NM\_001024843

Red=Cloning site Green=Tags(s)


MQTNEGEVSEESSKVEQEDFVMEGHGKTPPPGEE SKQEKEQEREEQLMEDKKRKKEDKKKKEATQKVTE  
QKTKVPEVTKPSLSQPTAASPIGSSPSPVNGGNNAKRVAVPNGQPPSAARYMPREVPPRFRCQQDHKVL  
LKRQPPPPSCMLLGGGAGPPCTAPGANPNNAQVTGALLQSESGTAPVWSKSTPPAPDNGTSAWGEPNE  
SSPGWGEMDDTGASTTGWGNTPANAPNAMKPNKSMQDGWGESDGPVTGARHPSWEEEDGGVWNTTGSQ  
GSASSHNSASWGGGKKQMKCSLKGNNDSWMNPLAKQFSNMGLLSQTEDNPSSKMDLSVGSLSDKKFDV  
DKRAMNLGDFNDIMRKDRSGFRPPNSKDMGTTDSGPYFEKGSGHGLFGNSTAQSRGLHTPVQPLNSSPSL  
RAQVPPQFISPQVSASMLKQFPNSGLSPGLFNVGPQLSPQQTAMLSQLPQIPQFLACQLLLQQQQQQQL  
LQNQRKISQAVRQQEQQLARMYSALQQQQQQQRQPGMKHSPSHPVGPKPHLDNMVPNALNVGLPDLQT  
KGPPIPGYGSFSSGGMDYGMVGGKEAGTESRFKQWTSMMEGLPSVATQEANMHKNGAIVAPGKTRGGSPY  
NQFDIIPGDTLGGHTGPAGDSWLPKSPPTNKIGSKSSNASWPPEFQPGVPWKGIQNIDPESDPYVTPGS  
VLGGTATSPIVDTHQLLRDNTTGSNSSLNTSLSPGAWPYSASDNSFTNVHSTSAKFPDYKSTWSPDPI  
GHNPTHLNKMWNHNISSRNTTLP RPPPGLTNPKPSSPWSSTAPRSVRGWTQDSRLASASTWSDGGSV  
RPSYWLVLHNLTPQIDGSTLRTICMQHGPLLTFHLNLTQGTALIRYSTKQEAQAQTALHMCVLGNTTIL  
AEFATDDEVSRLAQAQPPTPAATPSAPAAGWQSLETGQNQSDPVGPALNLFGGSTGLGQWSSAGGSSG  
ADLAGASLWGPPNYSSSLWGVPTVEDPHRMGSPAPLLPGDLLGGGSDSI

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:




ACCN:

NM\_001024843

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ORF Size:</b>              | 3087 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>                                                                      |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <a href="#">NM_001024843.1</a> , <a href="#">NP_001020014.1</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>RefSeq Size:</b>           | 15976 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq ORF:</b>            | 3090 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 23112                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <a href="#">Q9UPQ9</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 22q13.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Gene Summary:</b>          | <p>Plays a role in RNA-mediated gene silencing by both micro-RNAs (miRNAs) and short interfering RNAs (siRNAs). Required for miRNA-dependent translational repression and siRNA-dependent endonucleolytic cleavage of complementary mRNAs by argonaute family proteins. As scaffolding protein associates with argonaute proteins bound to partially complementary mRNAs and simultaneously can recruit CCR4-NOT and PAN deadenylase complexes.</p> <p>[UniProtKB/Swiss-Prot Function]</p>                              |

Product images:



Circular map for RG206220