

Product datasheet for **RG206227**

SPATA46 (NM_182581) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	SPATA46
Synonyms:	Clorf11l; HSD20
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG206227 representing NM_182581
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAGAACTTCTCACTCCTCAGCATCTCTGGACCTCCAATCTCTCCTCCGCCCTGAGTGCTTTTCCCG
ACATTATGTTCTCTCGTGCCACCAGCCTGCCAGACATTGCAAAGACAGCAGTACCCACTGAGGCATCCAG
CCCAGCTCAGGCCCTGCCACCCAGTACCAAAGCATCATTGTGAGGCAAGGGATACAGAACACAGTGCTC
TCACCAGACTGCAGCTTGGGGACACCCAGCACGGAGAGAAGCTGAGGCGGAAGTGCATATCTACCGGC
CCTGGTTCTCCCCCTACAGTACTTCGTGTGTGCAGACAAAGAGAGCCAGCTGGAGGCCTATGACTTCCC
AGAGGTGCAGCAGGATGAGGGCAAGTGGGACAAGTGCCTTTCTGAGGACATGGCTGAGAACATCTGTTTCG
TCCTCTTCTCCCCAGAGAACACTTGCCCTCGAGAAGCCACCAAGAAATCCAGGCATGGCCTGGACTCCA
TCACATCCCAGGACATCCTAATGGCTTCCAGATGGCACCCAGCACAGCAGAATGGCTACAAGTGGTGGC
CTGCTGCCGATGTACCCACCCCTGGACTTCTCAAGAGCCACATCAAGAGGGGCTTCAGGGAGGGCTTC
AGCTGCAAGGTGTACTACCACAAGCTCAAAGCCCTCTGGAGCAAGGAGCAGAAGGCCCGGCTGGGAGACA
GGCTCTCTCCGGCAGCTGCCAGGCCTCAATAGTCTGCTGAACACCTTAGGCAAATTGGCGGTGAAGC
CTACTTATGTCTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG206227 representing NM_182581
Red=Cloning site Green=Tags(s)

MENFSLLSISGPPISSSALSAPFDIMFSRATSLPDIKTAVPTEASSPAQALPPQYQSIIVRQGIQNTVL
SPDCSLGDTQHGEKLRRNCTIYRPWFSPYSYFVCADKESQLEAYDFPEVQQDEGKWDNCLSEDMAENICS
SSSSPENTCPREATKKSRHGLDSITSQDILMASRWHPAQNGYKCVACCRMYPTLDFLKSHIKRGFREGF
SCKVYVYHKLKALWSKEQKARLGDRLLSSGSCQAFNSPAEHLRQIGGEAYLCL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_182581

ORF Size: 783 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_182581.2](#)

RefSeq Size: 907 bp

RefSeq ORF: 786 bp


Locus ID: 284680

UniProt ID: [Q5T0L3](#)

Cytogenetics: 1q23.3

Gene Summary: Play a role in spermiogenesis and fertilization.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG206227