

## Product datasheet for **RG206238**

### PCDHA6 (NM\_018909) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PCDHA6 (NM_018909) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PCDHA6
Synonyms:	CNR2; CNRN2; CNRS2; CRNR2; PCDH-ALPHA6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG206238 representing NM\_018909  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**


ATGGTGTTTACCCCGGAGGATAGATTGGGAAAGCAATGTCTGCTCCTCCCGCTTCTGCTCCTCGCAGCCT  
 GGAAGGTGGGAGCGGCCAGCTCCACTACTCCGTACCCGAGGAGGCCAAACACGGCACCTTCGTGGGCCG  
 GATCGCGCAGGACCTGGGGCTGGAGCTGGCGGAGCTGGTGCCGCGCTGTTTCAGGATGGCCTCCAAGAC  
 CGCGAGGACCTTCTGGAGGTAATCTGCAGAATGGCATTGTTGTTGTGAATTCTCGGATCGACCGGAGG  
 AGCTGTGCGGGCGGAGCGGAGTGCAGCATCCACCTGGAGGTGATCGTGGACAGGCCGCTGCAGGTTTT  
 CCATGTGGACGTGGAGGTGAGGGACATTAACGACAACCCGCCCTTGTCCCGGTAGAGGAACAAAGAGTG  
 CTGATTTACGAATCTAGGCTGCCAGATTCTGTGTTTCCACTGGAGGGCGCGTCCGATGCAGATGTTGGCT  
 CAAATTCATCTTAACCTATAAACTCAGTTCTAGCGAATACTTCGGGCTAGATGTAAAAATAACAGTGA  
 TGACAATAAACAAATTGGGCTCTTATTAAGAAATCCTTGGACAGAGAGGAAGCTCTGCACACAACCTTA  
 TTCTGACAGCCACAGATGGGGGCAAACCTGAGCTCACAGGCACTGTTTCAGCTGCTGGTACAGTGTGG  
 ATGTGAATGATAATGCTCCCACTTTGAAACAGTCTGAATACGAAGTAAGAATATTCGAAAATGCAGACAA  
 CGGAACAACAGTTATCAGACTGAATGCTTCTGATCGGGATGAAGGAGCGAATGGGGCAATTCATATTCT  
 TTTAATAGCCTTGTTCAGCCATGGTTATTGACCACTTATGATAGATCGAAATACGGGAGAAAATAGTGA  
 TTCGGGGTAATTTGGATTTTGAACAAGAAAACCTATACAAAATCCTCATTGACGCCACGGACAAAGGCCA  
 TCCTCCATGGCGGGTATTGCACCGTTTTAGTGAGAAATTTGGATAAAAATGATAACGTCCCTGAGATA  
 GCACTGACTTCTTATCCTTGCCTGTACGTGAAGACGCTCAATTTGGTACTGTATCGCCCTAATTAGCG  
 TGAACGACCTCGATTAGGTGCCAACGGCAGGTGAAGTCTCGTGCAGCCTCAGTCCCTTTCAAGCT  
 GGTGTCCACCTTCAAGAATTACTACTCGTTGGTGTGGACAGTGCCTGGACCGGAGAGCGTTCGGCC  
 TATGAGTTGGTGTGACCGCGGGACGGGGCTCGCCTTCGCTGTGGGCCACAGCCAGCTTGTCTGTGG  
 AGGTGGCCGACATGAATGACAATGCTCCGGCGTTCGCGCAGCCGAGTACACAGTGTTCGTGAAGGAGAA  
 CAACCCGCGGGCTGCCACATCTTACGGTGTCTGCGGAGACGCGGACGCGCAGGAGAACGCGTGGTG  
 TCCTACTCGTGGTGGAGCGGGTGGGCGAGCGCGTTCGAGCTACATTTGGTGCACGCGGAGA  
 GCGGCAAGGTGTACGCGTGCAGCCGCTGGACCAGGAGCTAGAGCTGCTGCAGTTTCAGGTGAGCGC  
 GCGCGACGCGGGCTGCCCTCTGGGACGCAACGTGACGCTGCAGGTGTTCTGCTGGACGAGAACGAC  
 AACGCGCCGCGCTGCTGGCGCTCGGGTGGTGGTACTGGTGGTGCAGTGCAGGAGCTGGTCCCGGT  
 CACTGGGTGCAGCCAAGTGGTGGCGAAGTGCAGCGAGTTGACGCCGACTCAGGCTACAACGCGTGGCT  
 TTCGATGAGCTGCAGCCCCGGCAAGCAGCGCTCGCTTCCCGTTTCGCGTGGGGCTGTACACGGGCGAG  
 ATCAGCACCCTCGTGTCTGGACGAAGCGGACTCTCCGCGCCACCGGCTGCTGGTGTGGTGAAGACC  
 ACGGTGAGCCGGCGTGACAGCGACGGCCACGGTCTGGTGTGCTGGTGGAGAGTGGCCAGGCTCCAAA  
 GGCGTCATCACGGGCGTGGTGGCGCCGCGGGCCAGAGGCGGCGTGGTGGATGTCAACGTGTACCTG  
 ATCATCGCATCTGCGCGGTATCCAGCCTGCTGGTCTCAGCTACTGCTGTACACAGCGCTGCGGTGCT  
 CGGCGCCACCCACCGAGGGCGCGTGCACGGCGGACAAGCCACGCTGGTGTGCTCCAGCGCAGTGGGGAG  
 TGGTCTACTCGCAGCAGAGGGCGCAGAGGGTGTGCTCCGGGAGGGCCACCCAAGATGGATCTCATG  
 GCCTTTAGCCCCAGCCTTTCACCTTGTCTATTATGATGGGTAAGCGGAGAAATCAGATTTAAATGAAG  
 ATCATGATGCCAAACCAGCAGACCCAAACCTGACTGGCGTACTCTGCCTCCCTGAGAGCAGGCATGCA  
 CAGCTCTGTGCACCTAGAGGAGGCTGGCATTCTACGGGCTGGTCCAGGAGGGCTGATCAGCAGTGGCCA  
 ACAGTATCCAGTGAACACCAGAACCAGAGGAGGAGAGTGTCCCTCCAGTCCGTGCGGGTGTCAACA  
 GCAACAGCTGGACCTTTAAATACGGACAGGCAACCCCAAACAATCCGGTCCCGGTGAGTTGCCGACAA  
 ATTCATTATCCAGGATCTCTGCAATCATCTCCATCCGGCAGGAGCCTACTAACAGCCAAATTGACAAA  
 AGTGACTTCATAACCTTCGGCAAAAAGGAGGAGACCAAGAAAAAGAAGAAAAAGAAGGGTAACAAGA  
 CCCAGGAGAAAAAGAGAAAGGGAACAGCAGACTGACAACAGTGACCAG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_018909.2</a> , <a href="#">NP_061732.1</a>
<b>RefSeq Size:</b>	5287 bp
<b>RefSeq ORF:</b>	2853 bp
<b>Locus ID:</b>	56142
<b>UniProt ID:</b>	<a href="#">Q9UN73</a>
<b>Cytogenetics:</b>	5q31.3
<b>Domains:</b>	CA
<b>Protein Families:</b>	Secreted Protein, Transmembrane
<b>Gene Summary:</b>	This gene is a member of the protocadherin alpha gene cluster, one of three related gene clusters tandemly linked on chromosome five that demonstrate an unusual genomic organization similar to that of B-cell and T-cell receptor gene clusters. The alpha gene cluster is composed of 15 cadherin superfamily genes related to the mouse CNR genes and consists of 13 highly similar and 2 more distantly related coding sequences. The tandem array of 15 N-terminal exons, or variable exons, are followed by downstream C-terminal exons, or constant exons, which are shared by all genes in the cluster. The large, uninterrupted N-terminal exons each encode six cadherin ectodomains while the C-terminal exons encode the cytoplasmic domain. These neural cadherin-like cell adhesion proteins are integral plasma membrane proteins that most likely play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been observed and additional variants have been suggested but their full-length nature has yet to be determined. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG206238