

Product datasheet for **RG206313**

RALYL (NM_173848) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RALYL (NM_173848) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	RALYL
Synonyms:	HNRPCL3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG206313 representing NM_173848 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACTGGCAAAACAGACCAGCAACGTACCAATAAGAATGACCCCAAGTCCATCAACTCCCGTGTTT
TCATCGGCAATCTAAATACGGCAATTGTCAAGAAAGTTGACATTGAAGCCATTTTTCAAAGTATGGAAA
AATAGTTGGATGTTCCGTTACAAAGGTTATGCATTTGTACAGTACATGAGTGAGCGACATGCAAGAGCT
GCAGTGGCTGGAGAAAATGCCAGAGTCATCGCCGGCCAACCACTTGATATCAACATGGCAGGAGAGCCCA
AACCATACAGACAAAACCTGGAACAAGAGGCCCTTTCTGCACTTTACAGACTTGAATCAAAGGAACC
TTTCCTGTCTGTTGGCGTTATGTCTTTGACTATGATTACTACAGAGATGATTTCTACAATCGGTTATTT
GATTACCACGGGCGTGTGCCTCCACCTCCCGTGCAAGTAAATCCGCTGAAGCGTCCAGAGTGGCAGTCA
CAACGACTCGCAGGGGAAAGGAGTCTTTTCCATGAAAGGTGGATCGAGATCTACTGCCAGTGGGTCAAC
AGGTTCTAAATTGAAATCAGATGAGTTACAGACCATCAAGAAAGAATTAACCCAGATCAAACTAAAATT
GACTCCTTGCTAGGGCGCTGGAGAAGATTGAGAAACAGCAGAAGGCGGAGGCAGAAGCTCAGAAGAAGC
AATTGGAAGAGAGTCTAGTGCTGATCCAAGAGGAATGTGTGTGAGAGATTGCAGATCACTCTACAGAGGA
GCCTGCTGAAGGAGGCCAGATGCCGATGGAGAAGAGATGACAGATGGGATAGAGGAGGACTTCGATGAA
GATGGGGTTCATGAGCTGTTTCTACAGATAAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173848.7](#)

RefSeq Size: 2903 bp


RefSeq ORF: 876 bp

Locus ID: 138046

UniProt ID: [Q86SE5](#)

Cytogenetics: 8q21.2

Product images:



Circular map for RG206313