

Product datasheet for **RG206348**

ZNF595 (NM_182524) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF595
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG206348 representing NM_182524
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAACGTAACATTAGGGATGTGCCATAGAATTCTCCCTGAAGAGTGGAATGTCTGGACCTG
 CCCAGCAGAAATTTGTATAGAGATGTGATGTTGGAGAACTACAGGAACCTGGTCTCCCTGGGTTTTGTGAT
 CTCTAACCCAGACCTGGTCACCTGTCTGGAGCAAAATAAGAGCCCTGCAATTTGAAGATACATGAGACA
 GCAGCCAAACCCAGCTATATGTTCTCCTTCAGCCAAGACCTTTCACAGTGCAGGGGATAGAAGATT
 CATTCCACAACTTATACTGAAAAGATACGAGAAATGTGGACATGAGAATTTACAATTAAGAAAAGTTG
 TAAACGTGTGAATGAGTGAAGGTGCAGAAAGGAGTTAATAATGGAGTTTACCAGTGTGTTCAACTACC
 CAGAGCAAAATATTTCAATGTAATACATGTGTTAAAGTTTTAGTAAATTTCAAATTTCAAACAACATA
 AGATAAGACATACTGGAGAGAAACCTTTAAATGTACAGAATGTGCAGATCGTTTTACATGTCACACCT
 AACTCAACATACAGGAATTCATGCTGGAGAGAAACCTACAAATGTGAAAAATGTGGCAAAGCCTTTAAT
 AGGTCCACATCACTTAGTAAACATAAGAGAATTCATACTGGAGAGAAACCTACACATGTGAAGAATGTG
 GCAAAGCCTTTAGACGGTCCACAGTTCTGAACGAACATAAGAAAATTCATACTGGAGAGAAACCTACAA
 ATGTGAAGAATGTGGCAAAGCCTTTACAAGGTCCACAACACTGAATGAACACAAGAAAATTCATACTGGA
 GAGAAACCTTACAAATGTAAGAATGTGGCAAAGCCTTTAGATGGTCCACAAGCATGAATGAACATAAGA
 ATATTCATACTGGAGAGAAACCTACAAATGTAAGAATGTGGCAAAGCCTTTAGACAGTCCAGGAGCCT
 GAATGAACATAAAAAATTTACTACTGGCGAAAAACCTACACATGTGAAAAATGTGGCAAAGCCTTTAAC
 CAATCCTCAAGTCTTATTATACACAGGAGCATTCTGTAACAAAACTTTACAAATGTGAAGAATGTG
 GCAAAGCCTTTACTTGGTCTCATCCCTTAATAACATAAGAGAATTCATACTGGAGAGAAACCTACGC
 ATGTGAAGAATGTGGCAAAGCCTTTTATAGGTCCTCACACCTTGCTAAACATAAGAGAATTCATACTGGA
 GAGAAACCTACACGTGCGAAGAATGTGGCAAAGCCTTTAACCAATCCTCAACTCTTATATTACACAAGA
 GAATCCATTCTGGGCAAAAACCTTACAAATGTGAAGAATGTGGCAAAGCCTTTACACGGTCCACAACACT
 GAACGAACATAAGAAAATTCATACTGGCGAGAAACCTACAAATGTGAAGAATGTGGCAAAGCCTTTTATA
 TGGTCCGCAAGCCTGAATGAACATAAGAAATTTACTACTGGAGAGAAACCTACAAATGTAAGAATGTG
 GCAAAGCCTTTAACCAATCCTCAGGCCTTATTATACACAGGAGCATTCTGTAACAAAACTTTACAA
 ATGTGAAGAATGTGGCAAAGCCTTTACTCGGTCCACAGCCCTGAATGAACATAAGAAAATTCATTCTGGA
 GAGAAACCTTACAAATGCAAGAATGTGGCAAAGCCTTAACTTATCCTCAACCCTTACTAACATAAGA
 GAATTCATACTGGAGAGAAACCTTCACATGTGAAGAATGTGGCAAAGCCTTCAATTGGTCTCATCCCT
 TACTAACATAAGATAATTCATACTGGAGAGAAATCCTACAAATGTGAAGAATGTGGCAAAGCCTTTAAT
 CGGCCCTCAACCTTACTGTACACAAGCGAATTCATACTGGCAAGGAACATAGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG206348 representing NM_182524
 Red=Cloning site Green=Tags(s)

MELVTFRDVAIEFSPEEWKCLDPAQQNLYRDVMLENYRNLVSLGFVISNPDLVTCLEQIKEPCNLKIHET
 AAKPPAICSPFSQDLSPVQGIEDSFHKLILKRYEKCCHENLQLRKGCGRVNECKVQKGVNNGVYQCLSTT
 QSKIFQCNTCVKVFSKFSNSNKHKIRHTGEKPFKCTECGRSFMYSHTLQHTGIHAGEKPYKCEKCGKAFN
 RSTSLSKHKRIHTGEKPYTCEECGKAFRRSTVLNEHKKIHTGEKPYKCEECGKAFTRSTTLNEHKKIHTG
 EKPYKCEKCGKAFRWSTSMNEHKNIHTGEKPYKCEKCGKAFRQSRSLNEHKNIHTGEKPYTCEKCGKAFN
 QSSSLIIHRSIHSEQKLYKCEECGKAFWSSSLNKHKRIHTGEKPYACEECGKAFYRSSHAKHKRIHTG
 EKPYTCEECGKAFNQSSSLILHKRIHSGQKPYKCEECGKAFTRSTTLNEHKKIHTGEKPYKCEECGKAFI
 WSASLNEHKNIHTGEKPYKCEKCGKAFNQSSSLIIHRSIHSEQKLYKCEECGKAFTRSTALNEHKKIHTG
 EKPYKCEKCGKAYNLSSTLTAKHKRIHTGEKPFTEECGKAFNWSSSLTKHKIHTGEKSYKCEECGKAFN
 RPSTLTVHKRIHTGKEHS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_182524

ORF Size: 1944 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_182524.1, NP_872330.1</u>
RefSeq Size:	2895 bp
RefSeq ORF:	1947 bp
Locus ID:	152687
UniProt ID:	<u>Q8IYB9</u>
Cytogenetics:	4p16.3
Protein Families:	Transcription Factors
Gene Summary:	This gene encodes a protein belonging to the Cys2His2 zinc finger protein family, whose members function as transcription factors that can regulate a broad variety of developmental and cellular processes. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Oct 2013]

Product images:



Circular map for RG206348