

Product datasheet for **RG206392**

FGFBP2 (NM_031950) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FGFBP2 (NM_031950) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FGFBP2
Synonyms:	HBP17RP; KSP37
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG206392 representing NM_031950 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGTTCGTCCCCTGCCTCCTGCTGGTGACCTTGTCTGCCTGGGGACTTTGGGTGAGCCCCGAGGC
AAAAGCAAGGAAGCACTGGGGAGGAATTCATTTCCAGACTGGAGGGAGAGATTCCTGCACTATGCGTCC
CAGCAGCTTGGGGCAAGGTGCTGGAGAAGTCTGGCTTCGCGTCGACTGCCGCAACACAGACCAGACCTAC
TGGTGTGAGTACAGGGGGCAGCCAGCATGTGCCAGGCTTTTGTCTGCTGACCCAAACCTTACTGGAATC
AAGCCCTGCAGGAGCTGAGGCGCCTTACCATGCGTGCCAGGGGGCCCCGGTGTAGGCCATCCGTGTG
CAGGGAGGCTGGACCCAGGCCATATGCAGCAGGTGACTTCCAGCCTCAAGGGCAGCCAGAGCCCAAC
CAGCAGCCTGAGGCTGGGACGCCATCTCTGAGGCCAAGGCCACAGTGAACTCACAGAAGCAACACAGC
TGGGAAAGGACTCGATGGAAGAGCTGGGAAAAGCCAAACCCACCCGACCCACAGCCAAACCTACCCA
GCCTGGACCCAGGCCCGGAGGGAATGAGGAAGCAAGAAGAAGGCTGGGAACATTGTTGGAAACCTTCC
CAGGCCCTGTGCGCCTTCTCATCAGCTTCTCCGAGGG

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG206392 representing NM_031950
Red=Cloning site Green=Tags(s)

MKFVPCLLLVTLSCLGTLGQAPRQKQGSGTEEFHFQTGGRDSCMRPSSLGQGAGEVWLRVDCRNTDQTY
 WCEYRQPSPMCQAF AADPKPYWNQALQELRRLHHACQGAPVLRPSVCREAGPQAHMQQVTSLLKGSPEPN
 QQPEAGTPSLRPKATVKL TEATQLGKDSMEELGKAKPTTRPTAKPTQPGPRPGGNEEAKKKAWEHCWKPF
 QALCAFLISFFRG

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_031950

ORF Size: 669 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_031950.2</u> , <u>NP_114156.1</u>
RefSeq Size:	1199 bp
RefSeq ORF:	672 bp
Locus ID:	83888
UniProt ID:	<u>Q9BYJ0</u>
Cytogenetics:	4p15.32
Protein Families:	Secreted Protein
Gene Summary:	This gene encodes a member of the fibroblast growth factor binding protein family. The encoded protein is a serum protein that is selectively secreted by cytotoxic lymphocytes and may be involved in cytotoxic lymphocyte-mediated immunity. An increase in the amount of gene product may be associated with atopic asthma and mild extrinsic asthma.[provided by RefSeq Staff, Oct 2008]

Product images:



Circular map for RG206392