

Product datasheet for **RG206437**

MX2 (NM_002463) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MX2 (NM_002463) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MX2
Synonyms:	MXB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG206437 representing NM_002463
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCTAAGGCCACAAGCCTTGCCCTACCGGAGGAGAAGTCAATTTTCTTCTCGAAAATACCTGAAAA
AAGAAATGAATTCCTTCCAGCAACAGCCACCGCCATTCGGCACAGTGCCACCACAAATGATGTTTCTCCTC
AAACTGGCAGGGGCGAGAGAAGGACGCTGCTTTCCTCGCCAAGGACTTCAACTTTCTCACTTTGAACAAT
CAGCCACCACCAGGAAACAGGAGCCAACCAAGGCAATGGGGCCCGAGAACAACCTGTACAGCCAGTACG
AGCAGAAGGTGCGCCCTGCATTGACCTCATCGACTCCCTGCGGGCTCTGGGTGTGGAGCAGGACCTGGC
CCTGCCAGCCATCGCGTCATCGGGACCAGAGCTCGGGCAAGAGCTCTGTGCTGGAGGCACTGTCAGGA
GTCGCGCTTCCCAGAGGCAGCGGAATCGTAACCAGGTGTCGCTGGTGTGAACTGAAAAAGCAGCCCT
GTGAGGCATGGGCGGAAGGATCAGCTACCGGAACACCGAGCTAGAGCTTCAGGACCCTGGCCAGGTGGA
GAAAGAGATACACAAAGCCAGAACGTCATGGCCGGGAATGGCCGGGCATCAGCCATGAGCTCATCAGC
CTGGAGATCACCTCCCCTGAGGTTCCAGACCTGACCATCATTGACCTTCCCGGCATCACCGAGGTGGCTG
TGGACAACCAGCCCCGAGACATCGGACTGCAGATCAAGGCTCTCATCAAGAAGTACATCCAGAGGACGCA
GACGATCAACTTGGTGGTGGTTCCTGTAACTGGACATTGCCACCACGGAGGCGCTGAGCATGGCCAT
GAGGTGGACCCGGAAGGGGACAGGACCATCGGTATCCTGACCAAACAGATCTAATGGACAGGGGCACTG
AGAAAAGCGTCATGAATGTGGTGCAGAACCTCACGTACCCCTCAAGAAGGGCTACATGATTGTGAAGTG
CCGGGGCCAGCAGGAGATCACAACAGGCTGAGCTTGGCAGAGGCAACCAAGAAAGAAATTACATTTCTT
CAAACACATCCATATTTAGAGTTCTCCTGGAGGAGGGTTCAGCCACGGTTCCTGACTGGCAGAAAAGAC
TTACCCTGAACATCATGCATATCAAAAAATCGCTCCGTTGTTAGAAGGACAATAAGGGAGAGCCCA
CCAGAAGGCGACCGAGGAGCTGCGGCTTGGGGGCTGACATCCCCAGCCAGGAGGCGACAAGATGTTT
TTTCTAATTGAGAAAATCAAGATGTTTAACTCAGGACATCGAAAAGTTAGTAGAAGGAGAAGAAGTTGTAA
GGGAGAATGAGACCCGTTTATACAACAAAATCAGAGAGGATTTTAAAAACTGGGTAGGCATACTTGAAC
TAATACCAAAAAAGTTAAAAATATTATCCACGAAGAAGTTGAAAAATATGAAAAGCAGTATCGAGGCAAG
GAGCTTCTGGGATTTGTCAACTACAAGACATTTGAGATCATCGTGCATCAGTACATCCAGCAGCTGGTGG
AGCCCGCCCTTAGCATGCTCCAGAAAGCCATGAAATTTCCAGCAAGCTTTCATTAACGTGGCCAAAAA
ACATTTTGGCGAATTTTCAACCTTAACCAAACCTGTTCCAGACACGATTGAAGACATAAAAGTGAACAC
ACAGCAAAGGCAGAAAACATGATCCAATTCAGTTCAGAAATGGAGCAGATGGTTTTTTGTCAAGATCAGA
TTTACAGTGTGTTCTGAAGAAAGTCCGAGAAGAGATTTTAAACCCTCTGGGGACGCCTTACAGAAAT
GAAGTTGAACTCTCATTTTCCCAGTAATGAGTCTTCGTTTTCTCCTTTACTGAAATAGGCATCCACCTG
AATGCCTACTTCTTGAAACCAGCAAACGTCTCGCCAACCAGATCCCATTATAATTCAATTTATGTC
TCCGAGAGAATGGTGACTCCTTGCAGAAAGCCATGATGCAGATACTACAGGAAAAAAATCGCTATTCCTG
GCTGCTCAAGAGCAGAGTGAGACCGCTACCAAGAGAAGAATCCTTAAGGAGAGAATTTACCGGCTCACT
CAGGCGGACACGCACTCTGTCAATTCTCCAGCAAAGAGATCCAC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG206437 representing NM_002463
 Red=Cloning site Green=Tags(s)


MSKAHKPWYRRRSQFSSRKYLKEMNSFQQPPPGFTVPPQMMFPPNWQGAEKDAAFLAKDFNFLTNN
 QPPPGNRSQPRAMGPENNLYSQYEQKVRPCIDLIDSLRALGVEQDLALPAIAVIGDQSSGKSSVLEALSG
 VALPRGSGIVTRCPLVLKLLKQPCAWAGRISYRNTELELQDPGQVEKEIHKAQNVMAGNGRGISHELIS
 LEITSPEVPDLTIIDLPGITRVAVDNQPRDIGLQIKALIKKYIQRQQTINLVVPCNVDIATTEALSMAH
 EVDPEGDRITIGILTKPDLMDRGTEKSMNVVRNLTYPKKGMYIVKCRGQQEITNRLSLAEATKKEITFF
 QTHPYFRVLLLEEGSATVPRLAERLTTEELIMHIQKSLPLLEGQIREHQKATEELRRCGADIPSQEADKMF
 FLIEKIKMFNQDIEKLEVEGEEVVRENETRLYNKIREDFKNWVGILATNTQKVKNIIHEEVEKYEKQYRGK
 ELLGFVNYKTFEIIVHQYIQQLEVPALSMKAMEIIQAFINVAKKHFGFEFNLNQTQVSTIEDIKVKH
 TAKAENMIQLQFRMEQMVFCQDQIYSVVLKVVREEIFNPLGTPSQNMKLNHFPSSNESSVSSFTEIGIHL
 NAYFLETSKRLANQIPFIIQYFMLRENGDSLQKAMMQILQEKNRYSWLLQEQSETATKRRILKERIYRLT
 QARHALCQFSSKEIH


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_002463

ORF Size: 2145 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_002463.2](#)

RefSeq Size: 2961 bp

RefSeq ORF: 2148 bp

Locus ID: 4600

UniProt ID: [P20592](#)

Cytogenetics: 21q22.3

Domains: dynamin_2, dynamin, GED

Protein Families: Druggable Genome

Gene Summary: The protein encoded by this gene has a nuclear and a cytoplasmic form and is a member of both the dynamin family and the family of large GTPases. The nuclear form is localized in a granular pattern in the heterochromatin region beneath the nuclear envelope. A nuclear localization signal (NLS) is present at the amino terminal end of the nuclear form but is lacking in the cytoplasmic form due to use of an alternate translation start codon. This protein is upregulated by interferon-alpha but does not contain the antiviral activity of a similar myxovirus resistance protein 1. [provided by RefSeq, Jul 2008]