

Product datasheet for **RG206438**

TSBP1 (NM_006781) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	TSBP1
Synonyms:	C6orf10; TSBP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RG206438 representing NM_006781
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGACAGTCTTGGAAATACTTTGGCTGTATCCTGACTCTACTGGGACTTGCCATCCTGGCTATTTTGT
 TAACAAGATGGGCACGATGTAAGCAAAGTGAATGTATATCTCCAGATACAGTTCAGAACAAAGTGCTAG
 ACTTCTGGACTATGAGGATGGTAGAGGATCCCGACATGCATATTCAACACAAAGTGACATTCATATGAT
 AACCGAGAGAGATCCAAAGAGATTACACACCATCAACCAACTCTCTAGTGTCTATGGCATCTAAGTTCT
 CCCTGGGACAAACAGAACTCATTCTCTTTTGTGTGTTTTATTTAGCACTGTCTCGATCAAGTATTGG
 TAGTATAAAATGTTTACAAACAACCTGAAGAACCTCCTTCCAGAACTGCAGGAGCCATGATGCAATTCACA
 GCCCCTATTCGCGAGCTACAGGACCTATCAAGCTCTCTCAAAAAACCATTGTGCAAACTCCAGGACCTA
 TTGTACAATATCCTGGATCCAATGCTGGTCCACCTTCAGCACCCCGCGGTCCACCCATGGCACCCATAAT
 AATTTACAGAGAACCGAAGAATACCTCAAGTTCACACTATGGACAGTTCTGGAAAAATCACACTGACT
 CCTGTGGTTATTAACAGGTTACATGGATGAAGAACTTGCAAAAAATCTTGTTCAAAAATCCAGATTC
 TAAATGTGGAGGCACTGCAAGGTCTCAGAATAGCCGAGAAGAAAAACAAGGAAGCACTAAAGAAATGACAT
 CATATTTACGAATCTGTAGAATCCTTGAAATCAGCACACATAAAGGAGCCAGAAAGAGAAGGAAAAGGC
 ACTGATTTAGAGAAAGACAAAATAGGAATGGAGGTCAAGGTAGACAGTGACGCTGGAATACCAAAAAGAC
 AGGAAACCCAACTAAAAATCAGTGAGATGAGTATACCACAAGGACAGGGAGCCCAATAAAGAAAAGTGT
 GTCAGATGTACCAAGAGGACAGGAGTCCCAAGTAAAGAAGAGTGAGTCAGGTGTCCCAAAAGGACAAGAA
 GCCCAAGTAACGAAGAGTGGGTTGGTTGTACTGAAAGGACAGGAAGCCAGGTAGAGAAGAGTGAGATGG
 GTGTGCCAAGAAGACAGGAATCCCAAGTAAAGAAGAGTCAGTCTGGTGTCTCAAAGGACAGGAAGCCCA
 GGTAAAGAAGAGGAGTCAATTGTACTGAAAGGACAGGAAGCCAGGTAGAGAAGAGTGAGTTGAAGGTA
 CCAAAAGGACAAGAAGGCCAAGTAGAGAAGACTGAGGCAGATGTGCCAAAGGAACAAGAGGTCCAAGAAA
 AGAAGAGTGAGGCAAGGTGTACTGAAAGGACCAGAAATCCCAAGTAAAGAACACTGAGGTGAGTGTACCAGA
 AACACTGGAATCCCAAGTAAAGAAGAGTGAGTCAGGTGTACTAAAAGGACAGGAAGCCCAAGAAAAGAG
 GAGAGTTTTGAGGATAAAGGAAATATGATAAAGAAAAGGAGAGATGCAGAGAAAAGATCCAAATAAAA
 AAGAAAAGGTGACAAAAACAAAAAGGTGACAAAGGAAAGGACAAAGTTAAAGGAAAAGAGAGAATCAGA
 AATCAATGGTGAATAATCAAAAGGCTCGAAAGGGCGAAGGCAATACAGGAAGGAAGTACAACAAAAAA
 GTGGAAGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG206438 representing NM_006781
 Red=Cloning site Green=Tags(s)

MTVLEITLAVILTLLGLAILLILLTRWARCKQSEMYISRYSSSEQSARLLDYEDGRGSRHAYSTQSDTSYD
 NRERSKRDYPTNSLSVMSKFSLGQTELILLMCFILALSRSSIGSIKCLQTTEPPSRTAGAMMQFT
 APIPGATGPIKLSQKTIIVQTPGPIVQYPGSNAGPPSAPRGPPMAPIIISQRTARIPQVHTMDSSGKITLT
 PVVILTGYMDEELAKKCSKIQLKCGGTARSQNSREENKEALKNDIIFTNSVESLKSAAHIKEPEREGKG
 TDLEKDKIGMEVKVDSAGIPKRQETQLKISEMSIPQGGQAQIKKSVSDVPRGQESQVKKSESGVPKGQE
 AQVTKSGLVVLKGQEAQVEKSEMGVPRRQESQVKKSQSGVSKGQEAQVKKRESVVLKGQEAQVEKSELKV
 PKGQEGQVEKTEADVPKEQVEKSEAGVLKGPESQVKNTESVVPETLESQVKKSESGVLKGQEAQVEK
 ESFEDKGNNDKEKERDAEKDPNKKKGDKNKTKGDKGDKVKGKRESEINGEKSCKSKRAKANTGRKYNKK
 VEE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_006781

ORF Size: 1689 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_006781.5</u>
RefSeq Size:	2148 bp
RefSeq ORF:	1692 bp
Locus ID:	10665
UniProt ID:	<u>Q5SRN2</u>
Cytogenetics:	6p21.32
Protein Families:	Transmembrane

Product images:



Circular map for RG206438