

Product datasheet for **RG206718**

FAM81B (NM_152548) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM81B (NM_152548) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FAM81B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG206718 representing NM_152548 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCAATTACAATTCCTTGGTACATTGGCTTCCTCAGAAAAAGAAAAAATCACAGAGATTGTTTTTCA
AAAATATCAAATCTACAAAAATAAGCTGGAAGCAAGCATCATGAGTTCAGATACAAATGTAACAA
AAGTGCCTCTCCAAGTGCAGCTGCAGAGGAACAGCCAGTTGAACCTGATGGCCCCCTTCTGGCTCAGAC
AATAACCAAGAAAAAGTAAGATTATCTCCAGCCAAAATGTCAACCAAGAATTCTACAGATCTAGTTG
AATATGTTGACAAGAGTCATGCTTTTCTCCCATCATTCCAACACCCAGAGAGGTCAGCTAGAAGACAG
ACTGAACCAACCAGGCGCGTACCATAGCTTTCTTCTTGAACAAGCCTTCCGCATCAAGGAGGACATCTCT
GCTTGCCTGCAGGGGACCCATGGCTTTCGAAAAGAGGAATCGCTCGCCAGGAAGTTACTGAAAAGCCACA
TCCAGACCATCACCAGCATCGTCAAAAACTCAGCCAAAATATTGAGATTTAGAAAGACCAATAAGAGC
TCGAGATCAGGCGGCCACAGGAATACTTTGCAGTACAGGATAAACATCAACACCTACAAGGAGTT
GGAGATCTTCGAGGAAGAGTAGCCAGATGTGATTCAAGCATTGTGAAGCTTCTGGAGACATTCACTTAT
TCAGGCAAGAGCACCGGCAAATTGAGAAAGCCATTCAAGAATTCGTGCCCGCCCTGGAACTCTTCCAA
GAACCTGGACATGAAGGTGATGCAGCTCTTAGGAAAGATAGAACTGCCAGTTCTGAGCAACCTCGAAT
TTAAAGATGGTCCAGGGGATTATCGCCACGAAATGAACCTTTTGAATTCAAATTTCACTCACTTTCAA
GTAATCTGTACGAAGAAGTTGAGAATAATAAAAAATGGACAGAAAACCAATTTCTCAAATATAGAAAAGA
CCACCTGGGCCATATAAATGAATGTCTGAAGGTCTACAGGAGAACTGGAAAAGTCTGAAAATAAATG
GAAGAAAACTGCTGCAGCTTTCAAGCAAAGTAGAGAATTTCAATTAACACACAGAAACAGGAAACACAAC
TAAGTAAAGTAAAGCATATGGAATAAATTGTCCAAAAAGATGGAACAAATGGAAGCAGATCTGGGG
TGAATTAGAGACAATGCAGAATGAATATCAATCAGGATTTAAATCAATTCATGACTCTCTCAGCTCCCTC
CAACAAATACAGAAAACAAAGATGGATTTAGAGAAATATAAAGTACAGAAAGACCTAAAGAAATTACAGC
GCAAGATAGTGAACTCCAGGAAGTA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG206718 representing NM_152548
 Red=Cloning site Green=Tags(s)


MQLQLGLTLASSEKRKKSQRLFFKNIKSTKNKAGKASIMSSDTNVNKSASPTATAEEQPVEPDGGLPGSD
 NNQEKVRLSPAKMSTKNSTDLVEYVDKSHAFLPIIPNTQRGLEDRLNNQARTIAFLLEQAFRIKEDIS
 ACLQGTGFRKEESLARKLLESHIQITITSIVKKLSQNIIELEDQIRARDQAATGTNFVHEINIKHLQGV
 GDLRGRVARCDSSIVKLSGDIHLFRQEHRIEKAIQEFVPALETLSKNLDMKVMQLLGKIETASSEQTSN
 LKMVQGDYRHEMNLLFHFHSLSSNLYEEVENNKKWTENQFLKYRKDHLGHINECLKVLQEKLEKSENKM
 EEKLLQLSSKVENFINTQKQETQLSKVKHMENKLSKKMEQMEKQIWGELETMQNEYQSGFKSIHDSLSSL
 QQIQKTKMDLEKYKVQKDLKKLQRKIVELQEV

TRTRPLE – GFP Tag – V

Chromatograms: https://cdn.origene.com/chromatograms/ja1776_c11.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_152548


ORF Size: 1356 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_152548.3</u>
RefSeq Size:	1563 bp
RefSeq ORF:	1359 bp
Locus ID:	153643
UniProt ID:	<u>Q96LP2</u>
Cytogenetics:	5q15

Product images:



Circular map for RG206718