

Product datasheet for **RG206793**

ABAT (NM_000663) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ABAT (NM_000663) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ABAT
Synonyms:	GABA-AT; GABAT; NPD009
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG206793 representing NM_000663
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCTCCATGTTGCTCGCCAGCGCCTGGCCTGCAGCTTCCAGCACAGCTACCGCCTGCTGGTGCCTG
 GATCCAGACACATTAGTCAAGCTGCAGCCAAAGTCGACGTTGAATTTGATTATGATGGGCCCTCTGATGAA
 GACGGAAGTCCCAGGGCCTAGATCTCGGGAGTTAATGAAACAGCTGAATATAAATTCAGAATGCAGAGGCT
 GTGCATTTTTTCTGCAATTACGAAGAGAGCCGAGGCAATTACCTGGTTGATGTGGACGGCAACCGAATGC
 TGGATCTTTATCCAGATCTCCTCTGTCCCCATAGGTTACAGCCACCCCGCCCTGCTGAAACTCATCCA
 ACAGCCTCAAATGCGAGCATGTTGTCAACAGACCCGCCCTCGGAATCCTGCCTCCGGAGAACTTTGTG
 GAGAAGCTCCGGCAGTCTTGTCTCGGTGGCTCCCAAAGGGATGCCAGCTCATCACCATGGCCTGCC
 GCTCCTGTCCAATGAAAACGCCTAAAGACCATCTTCATGTGGTACCGAGCAAGGAAAGAGGGCAGAG
 GGGCTTCTCCAGGAGGAGCTGGAGACGTGCATGATTAACAGGCCCTGGCTGCCCGACTACAGCATC
 CTCTCCTTCATGGCGCGTTCCATGGGAGGACCATGGGTTGCTTAGCGACCACGCACTCTAAAGCCATTC
 ACAAGATCGACATCCCTTCTTTGACTGGCCATCGCACCGTTCCACAGGCTGAAATACCCCTCTGGAAGA
 GTTTGTGAAAGAGAACCAACAGGAGGAGGCCCGCTGTCTGGAAGAGGTGGAGGATCTGATTGTGAAATAT
 CGGAAAAAGAAGAAGACGGTGGCCGGGATCATCGTGGAGCCCATCCAGTCCGAGGGTGGAGACAACCACG
 CATCCGATGACTTCTTTTCGGAAGCTGAGAGACATCGCCAGGAAGCATGGCTGCGCCTTCTTGGTGGACGA
 GGTACAGACCGGAGGAGGCTGCACGGCAAGTTCTGGGCCATGAGCACTGGGGCCTGGATGACCCAGCA
 GACGTGATGACCTTCAGCAAGAAAATGATGACTGGGGCTTCTTCCACAAGGAGGAGTTCAGGCCTAATG
 CTCCTACCGGATCTTCAACACCTGGCTGGGGACCCGTCGAAGAACCTGTTGCTGGCTGAGGTCATCAA
 CATCATCAAGCGGGAGGACCTGCTAAATAATGCAGCCCATGCCGGGAAGGCCCTGCACAGGACTGCTG
 GACCTCCAGGCCCGGTACCCCCAGTTCATCAGCAGGGTGAAGGACGAGGCACCTTTTGCTCCTTCGATA
 CTCCCGATGATTCATACGGAATAAGCTCATTTTAATTGCCAGAAAACAAAGGTGTGGTGTGGGTGGCTG
 TGGTGACAAATCCATTTCGTTCCGTCCCACGCTGGTCTTCAGGGATCACACGCTCACCTGTTCTCAAT
 ATTTTCAGTGACATCTTAGCAGACTTCAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:


>RG206793 representing NM_000663
 Red=Cloning site Green=Tags(s)

MASMLLAQRLACSFQHSYRLLVPGSRHISQAAAKVDVEFDYDGPLMKTEVPGPRSRELMKQLNIIQNAEA
 VHFFCNYEESRGNLVDVDGNRMLDLYSQISSVPIGYSHPLLKLIQQPQNAMFVNRPALGILPPENFV
 EKLRLQSLLSVAPKGMSQLITMACGSCSNENALKTIFMWYRSKERGQRFSGEELETMINQAPGCPDYSI
 LSFMGAFHGRMGCLATTHSKAIHKIDIPFDWPIAPFPRLKYPLEEFVKENQQEEARCLEEVEDLIVKY
 RKKKKTIVAGIIVEPIQSEGGDNHASDDFFRKLRIARKHGCAFLVDEVQTTGGGCTGKFWAHEHWGLDDPA
 DVMTFSKMMTGGFFHKEEFRPNAPYRIFNTWLGDPKLNLLAEVINIIKREDLLNNAHAGKALLTGLL
 DLQARYPQFISRVGRGTFCSFDTPDSDIRNKLILIARNKGVVLGGCGDKSIRFRPTLVFRDHHHLFLN
 IFSIDILADFK

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Plasmid Map:


ACCN: NM_000663

ORF Size: 1500 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_000663.3 , NP_000654.2
RefSeq Size:	5463 bp
RefSeq ORF:	1503 bp
Locus ID:	18
UniProt ID:	P80404
Cytogenetics:	16p13.2
Domains:	aminotran_3
Protein Families:	Druggable Genome
Protein Pathways:	Alanine, aspartate and glutamate metabolism, beta-Alanine metabolism, Butanoate metabolism, Metabolic pathways, Propanoate metabolism, Valine, leucine and isoleucine degradation
Gene Summary:	<p>4-aminobutyrate aminotransferase (ABAT) is responsible for catabolism of gamma-aminobutyric acid (GABA), an important, mostly inhibitory neurotransmitter in the central nervous system, into succinic semialdehyde. The active enzyme is a homodimer of 50-kD subunits complexed to pyridoxal-5-phosphate. The protein sequence is over 95% similar to the pig protein. GABA is estimated to be present in nearly one-third of human synapses. ABAT in liver and brain is controlled by 2 codominant alleles with a frequency in a Caucasian population of 0.56 and 0.44. The ABAT deficiency phenotype includes psychomotor retardation, hypotonia, hyperreflexia, lethargy, refractory seizures, and EEG abnormalities. Multiple alternatively spliced transcript variants encoding the same protein isoform have been found for this gene. [provided by RefSeq, Jul 2008]</p>