

Product datasheet for **RG206801**

C19orf47 (NM_178830) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	C19orf47
Synonyms:	DKFZp686P05129; FLJ36888
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG206801 representing NM_178830
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTTGTGGATAATAGGATTGAGAGCATGCTGCTGGATCTCAATAAGGAGATAATGAATGAGCTGG
GCGTGACCGTGGTGGGTGACATCATCGCCATTCTCAAGCATGCCAAAGTGGTGACCGTCAGGACATGTG
CAAAGCTGCCACTGAGTCAGTACCCCTGCAGCCCTAGCCCCCTTGACGGCGAAATTCGCCGTGGCACCAGT
GCTGCCCTCCGAATGATCACCACAGCCTGAACCATGACTCTCCACCCAGCACACCCCCAGGCGCCCGG
ACACCAGCACCTCCAAGATCTCGGTCACTGTGTCCAACAAGATGGCAGCAAAGAGTGCCAAGGCCACTGC
AGCCCTGGCCCGCCGGGAGGAGAGAGCCTGGCTGTTCTGCAAGCGCGCCGGGTCACTGCTGAGATG
GAGGGGAAGTACGTCAATCAACATGCCAAAGGCACACCCCGCACCCGCAAGATCCTGGAGCAGCAGC
AGGCTGCAAAAGGTCTCATAGGACGTCTGTGTTGACCGCCTCGGCGCCGAGACCAAGGCAGACACCAC
GACAGGGAGTAAACCCACAGGAGTCTTCAGCCGCCTGGGGGCCACCCAGAAACGGACGAGGATCTGGCT
TGGGACAGCGACAATGACAGCAGCAGCTCTGTCTTGCAATATGCCGGGGTCTGAAGAAGCTAGGACGGG
GCCCAGCAAGGCCAGTCCCAGCCAGCACTGACTGTCAAAGCCAAGGCCACAAGCTCAGCGACAACGGC
TGCTGCCCCGACTGCGGCGCTGGCGCTTTCCTCAGGTCTGGGCTTGAGAGGAAGCCGAGTCTTG
TCTAAAGTCAGCATCATCAAGAGACTGGGCGCAGCTGCCCTTGTCGCCGAGGCCAGGACAGCCAGGTCA
CCAGCACCAAGAGTAAGTCTCAGCCGAGGTCAAGGTACCATTAAGAGGACTCTGGTGGGGCCCCGGGG
GAGCAGCTCCAGCGAGGGCCTTGGTGCCAGATGGACCACGCGGGCACTGTGAGCGTGTCAAAAGACTG
GGCCGACGACCTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG206801 representing NM_178830
Red=Cloning site Green=Tags(s)

MFVDNRIQKSMILLDLNKEIMNELGVTVVGDIIAILKHAKVVHRQDMCKAATESVPCSPSPLAGEIRRGTS
AASRMITNSLNHDSPPSTPPRRPDTSTKISVTVSNKMAAKSAKATAALARREESLAVPAKRRRVTAEM
EGKYVINMPKGTTPRTRKILEQQQAAKGLHRTSVFDRLGAETKADTTTGSKPTGVFSRLGATPETDEDLA
WDSNDNDSSSVLQYAGVLKKLGRGPAKASQPALTVKAKATSSATTAAPTLLRLALSSRGLERKPESL
SKVSIIRKLGAAALVPEAQDSQVTSTKSSSAEVKVTIKRTLVGPRGSSSEGLGAQMDHAGTVSVFKRL
GRRTF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_178830

ORF Size: 1065 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_178830.2, NP_849152.1


RefSeq Size: 2140 bp

RefSeq ORF: 1067 bp

Locus ID: 126526

Cytogenetics: 19q13.2

Product images:



Circular map for RG206801