

Product datasheet for **RG206829**

TBX3 (NM_016569) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TBX3 (NM_016569) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TBX3
Synonyms:	TBX3-ISO; UMS; XHL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG206829 representing NM_016569
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGAGCCTCTCCATGAGAGATCCGGTCATTCTGGGACAAGCATGGCCTACCATCCGTTCTACCTCACC
 GGGCGCCGGACTTCGCCATGAGCGCGGTGCTGGGTCAACCAGCCGCGTTCTTCCCGCGCTGACGCTGCC
 TCCCAACGGCGCGCGCGCTCTCGCTGCCGGCGCCCTGGCCAAGCCGATCATGGATCAATTGGTGGGG
 GCGGCCGAGACCGGCATCCCGTTCTCTCCCTGGGGCCCCAGGCGCATCTGAGGCCTTTGAAGACCATGG
 AGCCCCAAGAAGAGGTGGAGGACGACCCCAAGGTGCACCTGGAGGCTAAAGAAGTTGGGATCAGTTTCA
 CAAGCGGGGCACCGAGATGGTCATTACCAAGTCGGGAAGCGAATGTTTCTCCATTTAAAGTGAGATGT
 TCTGGGCTGGATAAAAAAGCCAAATACATTTTATTGATGGACATTATAGCTGCTGATGACTGTCTGTATA
 AATTTACAATTCTCGGTGGATGGTGGTGGTAAGGCCGACCCCGAAATGCCAAAGAGGATGTACATTCA
 CCCGGACAGCCCGCTACTGGGGAACAGTGGATGTCAAAGTCGTCACCTTCCACAAACTGAAACTCACC
 AACACATTTTCAGACAAACATGGATTTACTTTGGCCTTCCCAAGTGATCAGCTACGTGGCAGGGGAATT
 ATAGTTTTGGTACTCAGACTATATTGAATCCATGCACAATACCAGCCCCGGTTCACATTGTAAGAGC
 CAATGACATCTTGAAACTCCCTTATAGTACATTTCCGACATACTTGTCCCGAAACTGAATTCATCGCT
 GTGACTGCATACCAGAATGATAAGATAACCCAGTTAAAAATAGACAACAACCCCTTTTGCAAAAGGTTTCC
 GGGACTGGAATGGCCGAAGAGAAAAAGAAAAAGAGCTCACCCTGCAGTCCATGAGGGTGTGTTGATGA
 AAGACACAAAAAGGAGAATGGGACCTCTGATGAGTCTCCAGTGAACAAGCAGCTTTCAACTGCTTCGCC
 CAGGCTTCTTCTCCAGCCGCTCCACTGTAGGGACATCGAACCTCAAAGATTTATGTCCCAGCGAGGGTG
 AGAGCGACGCGGAGGCCGAGAGCAAAGAGGAGCATGGCCCCGAGGCTGCGACGCGGCCAAGATCTCCAC
 CACCACGTCGGAGGAGCCCTGCCGTGACAAGGGCAGCCCCGCGGTCAAGGCTCACCTTTTCGCTGTGAG
 CGGCCCGGGACAGCGGGCGGCTGGACAAGCGTCCGCCGACTCACGCCATAGCCCCGCCACCATCTCGT
 CCAGCACTCGCGGCTGGGCGGGAGGAGCGCAGGAGCCCGTTCGCGAGGGCACAGCGCCGGCCAAAGT
 GGAAGAGGCGCGCGCTCCCGGCAAGGAGGCCCTTCGCGCCGCTCACGGTGCAGACGGACGCGGCCGCC
 GCGCACCTGGCCAGGGCCCCCTGCCTGGCCTCGGCTTCGCCCGGGCTGGCGGGCCAACAGTTCTTCA
 ACGGGCACCCGCTCTTCTGCACCCAGCCAGTTTGCATGGGGGGCGCCTTCTCCAGCATGGCGGCCGC
 TGGCATGGTCCCCTCTGGCCACGTTTCTGGGGCTCCACCGTGTCTCGGGCTGGATTCCACGGCC
 ATGGCCTCTGCCGCTGCGGCGCAGGACTGTCCGGGGCTCCGCGGCCACCTGCCCTTCCACCTCCAGC
 AGCACGTCCTGGCCTCTCAGGGCTGGCCATGTCCCTTTCGGAAGCCTGTTCCCTTACCCTACACGTA
 CATGGCCGCAGCGCGGCCGCTCCTCTGCGGCAGCCTCCAGCTCGGTGCACCGCCACCCCTTCTCAAT
 CTGAACACCATGCGCCCGCGGCTGCGCTACAGCCCTACTCCATCCCGGTGCCGGTCCCGGACGGCAGCA
 GTCTGCTCACCACCGCCCTGCCCTCCATGGCGGCGGCCGCGGGGCCCTGGACGGCAAAGTCGCCGCCCT
 GGCCGCCAGCCCGCCTCGGTGGCAGTGGACTCGGGCTCTGAATCAACAGCCGCTCCTCCACGCTCTCC
 TCCAGCTCCATGTCCTTGTGCCAAACTCTGCGCGGAGAAAGAGGGGCCACCAGCGAATGCAGAGCA
 TCCAGCGTTGGTTAGCGGCTTGAAGCCAAGCCGGACAGTCCCGCAGCGGTCCCCG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG206829 representing NM_016569
 Red=Cloning site Green=Tags(s)

```


MSLSMRDPVIPGTSMAYHPFLPHRAPDFAMSAVLGHQPPFFPALTLPPNGAAALSLPGALAKPIMDQLVG
AAETGIPFSSSLGPQHLRPLKTMPEEEEEVEDDPKVHLEAKELWDQFHKRGTEMVITKSGRRMFPPFKVRC
SGLDKKAKYILLMDIIAADDCRYKFHNSRWMVAGKADPEMPKRMYPHDPSPATGEQWMSKVVTFHKLKLT
NNISDKHGFTLAFPSDHATWQGNYSFGTQITILNSMHKYQPRFHIVRANDILKLPYSTFRTYLFPETEFIA
VTAYQNDKITQLKIDNPNFAKGFRTDGNRRREKRKQLTLQSMRVFDERHKKENGTSESSSEQAAFNCFA
QASSPAASTVGTSNLKDLCPEGESDAEAEESKEEHGPEACDAAKISTTTSEEPCKDKGSPAVKAHLFAAE
RPRDSGRLDKASPDSPSRHSPATISSSTRGLGAEERRSPVREGTAPAKVEEARALPGKEAFAPLTVQTDAAA
AHLAQGPLPGLGFAPLAGQQFFNGHPLFLHPSQFAMGGAFSSMAAAGMGPLLATVSGASTGVSGLDSTA
MASAAAQGLSGASAATLPFHLQQHVLASQGLAMSPFGSLFPYPYTYMAAAAAASSAAASSSVHRHPFLN
LNTMRPRLRYPSPYIPVPVPDGSLLTTALPSMAAAGPLDGKVAALAASPASVAVDSGSELNSRSSTLS
SSSMSLSPKLCAEKEAATSELQSIQRLVSGLEAKPDRSRASP
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



ACCN: NM_016569

ORF Size: 2229 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016569.4](#)

RefSeq Size: 4814 bp

RefSeq ORF: 2232 bp

Locus ID: 6926

UniProt ID: [O15119](#)


Cytogenetics: 12q24.21

Domains: T-box

Protein Families: Druggable Genome, Transcription Factors

Gene Summary: This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. This protein is a transcriptional repressor and is thought to play a role in the anterior/posterior axis of the tetrapod forelimb. Mutations in this gene cause ulnar-mammary syndrome, affecting limb, apocrine gland, tooth, hair, and genital development. Alternative splicing of this gene results in three transcript variants encoding different isoforms; however, the full length nature of one variant has not been determined. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG206829