

## Product datasheet for **RG206910**

### STK33 (NM\_030906) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	STK33 (NM_030906) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	STK33
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG206910 representing NM\_030906  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGCTGATAGTGGCTTAGATAAAAAATCCACAAAATGCCCCGACTGTTTCATCTGCTTCTCAGAAAGATG  
 TACTTTGTGTATGTTCCAGCAAAAACAAGGGTTCCTCCAGTTTTGGTGGTAAATGTCACAGACATCAAG  
 CATTGGTAGTGCAGAATCTTTAATTTCACTGGAGAGAAAAAAGAAAAAATATCAACAGAGATATAACC  
 TCCAGGAAAGATTTGCCCTCAAGAACCTCAAATGTAGAGAGAAAAGCATCTCAGCAACAATGGGGTCGGG  
 GCAACTTTACAGAAGGAAAAGTTCCTCACATAAGGATTGAGAAATGGAGCTGCTATTGAGGAAATCTATAC  
 CTTTGGAAAGAATATTGGGAAAAGGGAGCTTTGGAATAGTCATTGAAGCTACAGACAAGGAAAACAGAAACG  
 AAGTGGCAATTAAGGAAAGTGAACAAAAGAAAAGGCTGGAAGCTCCGCTGTGAAGTACTTGAACGAGAGG  
 TGAACATTTCTGAAAAGTGTAAAACATGAACACATCATACTGGAACAAGTATTTGAAACGCCAAAGAA  
 AATGTACCTTGTGATGGAGCTTTGTGAGGATGGAGAACTCAAAGAAATTTGGATAGGAAAAGGCATTTTC  
 TCAGAGAATGAGACAAGGTGGATCATTCAAAGTCTCGCATCAGCTATAGCATATCTTCAATAATGATA  
 TTGTACATAGGGATCTGAAACTGGAAAATATAATGGTTAAAAGCAGTCTTATTGATGATAACAATGAAAT  
 AAACCTTAAACATAAAGGTGACTGATTTTGGCTTAGCGGTGAAGAAGCAAAGTAGGAGTGAAGCCATGCTG  
 CAGGCCACATGTGGGACTCCTATCTATATGGCCCTGAAGTTATCAGTGCCACGACTATAGCCAGCAGT  
 GTGACATTTGGAGCATAGGCGTCGTAATGTACATGTTATTACGTGGAGAACCACCCTTTTTGGCAAGCTC  
 AGAAGAGAAGCTTTTTGAGTTAATAAGAAAAGGAGAACTACATTTGAAAATGCAGTCTGGAATTCATA  
 AGTGACTGTGCTAAAAGTGTTTGAAACAACCTTATGAAAGTAGATCCTGCTCACAGAATCACAGCTAAGG  
 AACTACTAGATAACCAAGTGGTTAACAGGCAATAAACTTTCTTCGGTGAGACCAACCAATGTATTAGAGAT  
 GATGAAGGAATGGAAAAATAACCCAGAAAAGTGTGAGGAAAACACAACAGAAGAGAAGAATAAGCCGTCC  
 ACTGAAGAAAAGTTGAAAAGTTACCAACCTGGGAAAATGTCCCTGAGACCAATTACACTTCAGATGAAG  
 AGGAGGAAAAACAGTCTACTACTTATGAAAAGCAATTTCTGCAACCAGTAAGGACAACCTTTGATATGTG  
 CAGTTCAAGTTTACATCTAGCAAACCTCTCCAGCTGAAATCAAGGGAGAAATGGAGAAAAACCCCTGTG  
 ACTCCAAGCCAAGGAACAGCAACCAAGTACCCTGCTAAATCCGGCGCCCTGTCCAGAACCAAAAAGAAAC  
 TC

**ACGCGTACGCGGCCGCTCGAG** – GFP Tag – GTTTAA

**Protein Sequence:**

>RG206910 representing NM\_030906  
 Red=Cloning site Green=Tags(s)

MADSGLDKKSTKCPDCSSASQKDVLCVCSKTRVPPVLVEMSQTSIGSAESLISLERKKEKNINRDIT  
 SRKDLPSRTSNVERKASQQWGRGNFTEGKVPHIRIENGAAIEEYTFGRILGKGSFGIVIEATDKETET  
 KWAICKVNKEKAGSSAVKLLEREVNILKSVKHEHIIHLEQVFETPKKMYLVMELCEDGELKEILDRKGHF  
 SENETRWIIQSLASAIAYLHNNDIVHRDLKLENIMVKSLLIDNNEINLNIKVDFGLAVKKQSRSEAML  
 QATCGTPIYMAPEVISAHDYSQQCDIWSIGVVMYMLLRGEPFLASSEEKLFELIRKGELHFENAVNSI  
 SDCAKSVLKQLMKVDPAHRITAKELLDNQWLTGNKLSVRPTNVLEMMKEWKNNPESVEENTTEEKNKPS  
 TEEKLSYQPWGNVPETNYTSDEEEEKQSTTYEKQFPATSKDNFDMCSSSFTSSKLLPAEIKGEMEKTPV  
 TPSQGTATKYPKSGALSRTKKKL

**TRTRPLE** – GFP Tag – V

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_030906

**ORF Size:** 1542 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_030906.2](#), [NP\\_112168.1](#)

**RefSeq Size:** 2707 bp

**RefSeq ORF:** 1545 bp

**Locus ID:** 65975


**UniProt ID:** [Q9BYT3](#)

**Cytogenetics:** 11p15.4

**Protein Families:** Druggable Genome, Protein Kinase

**Gene Summary:** Serine/threonine protein kinase which phosphorylates VIME. May play a specific role in the dynamic behavior of the intermediate filament cytoskeleton by phosphorylation of VIME (By similarity). Not essential for the survival of KRAS-dependent AML cell lines.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for RG206910