

## Product datasheet for **RG206952**

### MCU (NM\_138357) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	MCU (NM_138357) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MCU
Synonyms:	C10orf42; CCDC109A; HsMCU
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG206952 representing NM_138357 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGGCCGCCGAGGTAGATCGCTCCTGCTCCTCTCCTCTCGGGCGGGCGGGCGGGGGCGCCG  
GCGGCTGCGGGCGCTGACTGCCGGTCTTCCCTGGGCTGGGCGTCAGCCGCCACCGCAGCAGCAGCA  
CCACCGGACGGTACACCAGAGGATCGTTCCTGGCAGAATTTGGGAGCTGTTTATTGCAGCACTGTTGTG  
CCCTCTGATGATGTTACAGTGGTTATCAAATGGTTACCTGTGATATCTGTGAGGCTACCATCCCGGC  
GTGAACGCTGTCACTCACACTCAAGCCTATCTCTGACTCTGTTGGTGTATTTTTACGACAACACTGCAAGA  
AGAGGATCGGGGAATTGACAGAGTTGCTATCTATTACCAGATGGTGTTCGCGTTGCTGCTTCAACAGGA  
ATAGACCTCCTCCTCTTATGACTTTAAGCTGGTCATTAATGACTTAACATACCACGTACGACCACCAA  
AAAGAGACCTCTTAAGTCATGAAAATGCAGCAACGCTGAATGATGTAAGACATTGGTCCAGCAACTATA  
CACCACACTGTGCATTGAGCAGCACCAGTTAAACAAGGAAAGGGAGCTTATTGAAAGACTAGAGGATCTC  
AAAGAGCAGCTGGCTCCCTGGAAAAGGTACGAATTGAGATTAGCAGAAAAGCTGAGAAGAGGACCACTT  
TGGTGCATGGGGTGGCCTTGCCTACATGGCCACACAGTTTGGCATTGTTGGCCGGCTTACCTGGTGGGA  
ATATTCCTGGGACATCATGGAGCCAGTAACATACTTCATCACTTATGGAAGTGCCATGGCAATGTATGCA  
TATTTTGTAAATGACACGCCAGGAATATGTTTATCCAGAAGCCAGAGACAGACAATACTTACTATTTTCC  
ATAAAGGAGCCAAAAGTCACGTTTTGACCTAGAGAAAATACAATCAACTCAAGGATGCAATTGCTCAGGC  
AGAAAATGGACCTTAAGAGACTGAGAGACCCATTACAAGTACATCTGCCTCTCCGACAAATTGGTAAAAA  
GAT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG206952 representing NM\_138357  
 Red=Cloning site Green=Tags(s)

MAAAAGRSLLLLLSSRGGGGGAGGCGALTAGCFPGLGVSRRHQQHHRTVHQRIASWQNLGAVYCSTVV  
 PSDDVTVVYQNGLPVISVRLPSRRERCQFTLKPISDSVGVFLRQLQEEDRGIDRVAIYSPDGVVAASTG  
 IDLLLLDDFKLVINDLTYHVRPPKRDLLSHENAATLNDVKTLVQQLYTTLCIEQHQLNKERELIERLEDL  
 KEQLAPLEKVRIEISRKAEKRTTLVLWGGLAYMATQFGILARLTWWEYSWDIMEPVTYFITYGSAAMAYA  
 YFVMTRQEVVYPEARDRQYLLFFHKGAKKSRFDLEKYNQLKDAIAQAEMDLKRLRDPLQVHLPLRQIGEK  
 D

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_138357

**ORF Size:** 1053 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_138357.2](#)

**RefSeq Size:** 2978 bp

**RefSeq ORF:** 1056 bp

**Locus ID:** 90550

**UniProt ID:** [Q8NE86](#)


**Cytogenetics:** 10q22.1

**Domains:** DUF607

**Protein Families:** Transmembrane

**Gene Summary:** This gene encodes a calcium transporter that localizes to the mitochondrial inner membrane. The encoded protein interacts with mitochondrial calcium uptake 1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2012]

## Product images:



Circular map for RG206952