

Product datasheet for **RG207009**

TAF1C (NM_139353) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TAF1C (NM_139353) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TAF1C
Synonyms:	MGC:39976; SL1; TAFI95; TAFI110
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG207009 representing NM_139353
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGCCTCCCCTCATCGATCCCTGGGACCTGGCCTGACTGCCCGGGACCTGCTTTCCGCGGAGGGT
 ACCGGTATCGGAAGCGGCCCGAGTCGTGCTGGATGTGACTGAGCAGATCAGCCGTTCCCTTTGGATCA
 TGGAGACGTAGCCTTTGCGCCCTGGGGAAGCTGATGCTGGAGAATTTCAAGCTGGAGGGAGCGGGGAGC
 CGCACTAAGAAGAAGACAGTGGTCAGTGTGAAGAAGCTGCTCCAGGACCTCGGTGGACACCAGCCCTGGG
 GGTGTCCCTGGGCTTACCTCAGCAACCGACAGCGCCGCTTCTCTATCCTCGGGGGCCCCATCTGGGCAC
 GTCGGTGGCGAGCCACTTGGCAGAGCTGCTGCACGAGGAGCTGGTGTGCGGTGGGAGCAGCTGCTTCTG
 GATGAGGCCTGCACTGGGGGCGCGCTGGCTGGGTTCTGGAAGGACACCCAGTTCGGGCAGCTGGTCT
 ACCCTGCTGGAGGCGCCAGGACAGGCTGCATTTCCAAGAGTCTGTTCTGACCCAGGTGACAATCCCCA
 ATTCCTTTGGAAACCTGGACGCATCCAGCTCCAGGGACCTGTCCGGCAAGTGGTACATGCACCGTCCAG
 GGAGAAACTCTGCTGGCGTCCGCTCTGACTACCACTGTGCCGTGTGAAAGTTTGGTAAACAGTGGCAGC
 CAACCTTCTGCAGGCGATGCAGGTGGAGAAAGGGGCCACGGGGATCAGCCTCAGCCCTCACTGCCCGG
 GGAGCTGGCCATCTGCAGCCGCTCGGGAGCCGTCTGCCTGTGGAGCCCTGAGGATGGGCTGCGGCAATC
 TACAGGGACCCTGAGACCCTCGTGTCCGGGACTCCTCTCTGTGGCGTTGGGAGACTTCACTGCGCACC
 CTCGGGTGCTGACCGTGGGTGACCGCACCGGAGTGAAGATGCTGGACACTCAGGGCCCGCCGGGCTGTGG
 TCTGTTGCTTTTTCGTTTGGGGGAGAGGCTTCGTGCCAGAAAGGGGAACGTGCTCTGTTACCCAGTAC
 CTGGGGCACTCCAGCCCCAAATGCCTCCCCCTACTCTTCATCTCGTCTGTACCCAGTCTCTCTTACC
 TAGTGGACGAGCGCCTTCCCCTGGTGGCGATGCTGAAGTGAACCATGGCCTCCCCTCCCGCTCCTGCT
 GGCCCGACTGCTGCCTCCGCCCCGGCCAGCTGCGGTGCAGCCCTGCTCCTCGGAGGCCAGGGTGGGCAG
 CTGCAGCTGCTGCACCTGGCAGAAGGGCGTGGTGGCCCGCTGGCAGGCCCCCCCAGTCTTCTCCTT
 CCAGGATCGACTCCCTCCCTGCATTTCTCTGCTGGAGCCTAAGATCCAGTGGCGGCTGCAGGAGCGCT
 GAAAGCACCGACCATAGGTCTGGCTGCCGTCTCCCGCCATGCCCTCAGCGCCACACCAGGCCTGGTG
 CTCTCCAGCTCTCGGCGGCGGAGATGTCTTCTACCAGCAGCTCCGCCCCAGGTGGACTCCAGCCTCC
 GCAGAGATGCTGGGCTCCTGGCAGACCCAACCTGACTGCCATGCCCCACAGCTTCTGGACCTCCCA
 GGACTGCGCGGCTGCAGCCAGTGGTGAAGGCCCTGCTAAAAGTGGCCCTGGCTCCTCCTGTGTGGACA
 GCACCCACCTTACCCACCGCCAGATGCTGGGCAGCACAGAGCTGCGGAGGGAGGAAGAGGAAGGGCAGC
 GGCTGGGTGTACTCCGCAAGGCCATGGCCCGAGGGCAGCTCCTGCTGCAGAGAGACTGGGCTCCCTCCC
 TGGCGCAGAGCCACCCCTGCACCCGAGTCAGGCCTAGAGGACAAGCTCAGTGAAGCGCTGGGGGAGCC
 TGGGCAGGCCGAGGGGCTGCCTGGTGGGAGAGGCAGCAGGGCAGGACCTCGGAGCCCGGGAGACAGACCA
 GGCGGCCCAAGCGCCGACCCAGCTGTCCAGCAGCTTTTCGCTCAGTGGCCATGTGGATCCCTCAGAGGA
 CACCAGTCCCCTCATAGCCCTGAGTGGCCACCTGCTGATGCTCTGCCCTGCCCCACGACCCCGCCC
 TCCCAGGAGTTGACTCCGGATGCATGCGCCAGGGCGTCCCATCAGAGCAGCGGCAGATGCTCCGTGACT
 ACATGGCCAAGCTACCACCCAGAGGGACACCCAGGCTGTGCCACCACCTCCCCTCCCAGGCCTC
 CAGCGTCCGGGCCACTCGCTCCCAGCAGCACACCCCGTCTCTAGCTCTCAGCCCTCCGTAAGAAG
 CCTCGAATGGGCTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG207009 representing NM_139353
 Red=Cloning site Green=Tags(s)


MLPPLIDPWPGLTARDLLFRGGYRYRKRPRVLDVTEQISRFLLDHGDVAFAPLGKLMLENFKLEGAGS
 RTKKKTVVSVKLLQLDGGHPWGPCWAYLSNRQRRFSILGGPILGTSVASHLAELLHEELVLRWEQLLL
 DEACTGGALAWVPGRTPOFGQLVYPAGGAQDRLHFQEVVLTGPDNPQFLGKPGRIQLQGPVVRQVVTCTVQ
 GETLLAVRSDYHCAVWKF GKQWQPTLLQAMQVEKGATGISLSPHLPGELAIICSRSGAVCLWSPEDGLRQI
 YRDPETLVFRDSSSWRWADFTAHRVLTVDGRTGVKMLDTQGPPGCGLLLFRLGAEASCQKGERVLLTQY
 LGHSSPKCLPPTLHLVCTQFSLYLVDERLPLVPMKWNHGLPSPLLLARLLPPPRPSCVQPLLLGGQGGQ
 LQLLHLAEGASVPRLAGPPQSLPSRIDSLEPAFPLEPKIQWRLQERLKAPTIGLAAVVPPMPSAPTGLV
 LFQLSAAGDVFYQLRPQVDSLLRRDAGPPGDTQPDCHAPTASWTSQDTAGCSQWLKALLKVPLAPPVWT
 APTFTHRQMLGSTELRREEEEGQRLGVLKAMARGQLLQRDLGSLPAAEPPAPESGLEDKLSERLGEA
 WAGRGAAWWERQQGRTSEPGRQTRRPKRRQLSSSFSLSGHVDPS EDTSSPHSPEWPPADALPLPPTTP
 SQELTPDACAQGVPSERQMLRDYMAKLPQRDTPGCATTPPHSQASSVRATRSQOHTPVLSSSQPLRKK
 PRMGF


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:


ACCN: NM_139353

ORF Size: 2325 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_139353.1](#), [NP_647610.1](#)

RefSeq Size: 3741 bp

RefSeq ORF: 2328 bp

Locus ID: 9013

UniProt ID: [Q15572](#)

Cytogenetics: 16q24.1

Protein Families: Transcription Factors

Gene Summary: Initiation of transcription by RNA polymerase I requires the formation of a complex composed of the TATA-binding protein (TBP) and three TBP-associated factors (TAFs) specific for RNA polymerase I. This complex, known as SL1, binds to the core promoter of ribosomal RNA genes to position the polymerase properly and acts as a channel for regulatory signals. This gene encodes the largest SL1-specific TAF. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2011]