

## Product datasheet for **RG207099**

### **Epsin 1 (EPN1) (NM\_013333) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Epsin 1 (EPN1) (NM_013333) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Epsin 1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG207099 representing NM\_013333  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGTCGACCTCGTCTTGAGGCGCCAGATGAAGAACATCGTCCCAACTACTCAGAGCGGAGATCAAGG  
 TTCGAGAGGCCACGAGCAATGACCCCTGGGGCCATCCAGCTCCCTCATGTGAGAGATTGCCGACCTCAC  
 CTAACAACGTTGTGCGCTTCTCGGAGATCATGAGCATGATCTGGAAGCGGCTCAATGACCATGGCAAGAAC  
 TGGCGTCACGTTTACAAGGCCATGACGCTGATGGAGTACCTCATCAAGACCGGCTCGGAGCGCGTGTGCG  
 AGCAGTGAAGGAGAACATGTACGCCGTGACAGCCTGAAGGACTTCCAGTACGTGGACCGGACGGCAA  
 GGACCAGGGCGTGAACGTGCGTGAGAAAGCTAAGCAGCTGGTGGCCCTGCTGCGGACGAGGACCGGCTG  
 CGGAAGAGCGGGCGCACGCGCTCAAGACCAAGGAAAAGCTGGCAGACCGCCACGCGCTCATCAGCAG  
 CTGTGGGCTCAGGCCCTCCCGAGGCGGAGCAGGCGTGGCCGACAGCAGCGGGGAGGAGGAGCTGCA  
 GCTCCAGCTGGCCCTGGCCATGAGCAAGGAGGAGCCGACAGGAGGAGCGGATCCGTCGCGGGGATGAC  
 CTGCGGCTGCAGATGGCAATCGAGGAGAGCAAGAGGGAGACTGGGGGCAAGGAGGAGTCCCTCCATGG  
 ACCTTGCTGACGTCTTACGCGCCCAAGCTCTGCCCGACACAGACCCCTGGGGGGGCCAGCACCCAT  
 GGCTGCTGCCGTCCCACGGCTGCCCCACCTCGGACCCCTGGGGCGGCCCTGTCCCTCCAGCTGCT  
 GATCCCTGGGAGGTCCAGCCCCACGCGGCCTCTGGGACCCCTGGAGGCTGCTGCCCTGCAGGAC  
 CCTCAGTTGACCTTGGGGTGGGACCCAGCCCTGCAGCTGGGAGGGGGCCACGCCTGATCCATGGGG  
 AAGTCCGATGGTGGGGTCCCGTCACTGGGCCCTCAGCTCCGATCCCTGGACACCGGGCCCGGCTTC  
 TCAGATCCCTGGGAGGGTCACTGCCAAGCCAGCACAATGGCACAACAGCGGGGATTGACACGG  
 AGCCCGACGAGTTCTCTGACTTTGACCGACTCCGACGGCACTGCCGACCTCCGGGAGCAGCGCAGGAGA  
 GCTGGAGCTGCTGGCAGGAGAGGTGCCGCGCCGAAGCCCTGGGGGCTTGGACATGAGTGGGGTCAAGGGA  
 TCTCTGGCTGAGGCTGTGGGAGCCCCACCTGCAGCCACACCAACTCCACGCCCCCACCAGGAGA  
 CGCCGGAGTATTCTGGGGCCCAATGCAGCCCTCGTCGACTGGACTCGCTGGTGGAGCCGGCCGGGCC  
 CACGCGCCTGGAGCCAAGGCCTCAACCCCTTCTGCCAGGCGGAGGCCAGCCACTGGCCCTCCGTC  
 ACCAACCCCTTCCAGCCCGCCTCCCGGACGCTCACCTGAACCAGCTCCGTCTCAGTCTGTGCCTC  
 CCGTCCCTGGAGCGCCACCCACGTACATCTCTCCCTTGGCGGGGCCCTGGCTGCCCCCATGATGCC  
 CCCGGCCCCCGGCCCAACTAATCCCTTCTCTA

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:**

>RG207099 representing NM\_013333  
 Red=Cloning site Green=Tags(s)

MSTSSLRQMKNI VHNYSAEIKVREATSNDPWGPSSSLMSEIADLTYNVVFSEIMSMIWKRLNDHGKN  
 WRHVYKAMTLMELYIKTSERVSQCKENMYAVQTLKDFQYVDRDQKQGVNVREKAKQLVALLRDEDRL  
 REERAHALKTKKLAQTATASSA AVSGPPPEAEQAWPQSSGEEELQLQLALAMSKEEADQEERIRRGDD  
 LRLQMAIEESKRETGGKEESSLMDLADVFTAPAPPTTDPWGGPAPMAAAVPTAAPTSDPWGPPVPPAA  
 DPWGGPAPTPASGDPWRPAAPAGPSVDPWGGTPAPAAGEPTDPWSSDGGVPVSGPSASDPWTPAPAF  
 SDPWGSPAKPSTNGTTAGGFDTEPDEF SDFDRLRLALPTSGSSAGELELLAGEVPARSPGAFDMSGVRG  
 SLAEAVGSPPPAATPTPTPPTKTPESFLGPNAALVDLDSL VSRPGTPPGAKASNPFLPGGGPATGPSV  
 TNFPQPAPPATLTLNQLRLSPVPPVPGAPPTYISPLGGGPGLPMPMPGPPAPNTNPFLL

**TRTRPLE** - GFP Tag - V

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_013333

**ORF Size:** 1650 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_013333.3](#)

**RefSeq Size:** 2430 bp

**RefSeq ORF:** 1653 bp

**Locus ID:** 29924

**UniProt ID:** [Q9Y6I3](#)


**Cytogenetics:** 19q13.42

**Domains:** ENTH, UIM

**Protein Pathways:** Endocytosis

**Gene Summary:** This gene encodes a member of the epsin protein family. The encoded protein binds clathrin and is involved in the endocytosis of clathrin-coated vesicles. Loss of function of this gene is associated with reduced tumor growth and progression in certain cancer types. [provided by RefSeq, Mar 2016]

**Product images:**



Circular map for RG207099