

## Product datasheet for **RG207351**

### **R3HDM2 (NM\_014925) Human Tagged ORF Clone**

#### **Product data:**

|                                  |                         |
|----------------------------------|-------------------------|
| <b>Product Type:</b>             | Expression Plasmids     |
| <b>Tag:</b>                      | TurboGFP                |
| <b>Symbol:</b>                   | R3HDM2                  |
| <b>Synonyms:</b>                 | CAG6; PR01365           |
| <b>Mammalian Cell Selection:</b> | Neomycin                |
| <b>Vector:</b>                   | pCMV6-AC-GFP (PS100010) |
| <b>E. coli Selection:</b>        | Ampicillin (100 ug/mL)  |



**ORF Nucleotide Sequence:** >RG207351 representing NM\_014925  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTGCACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGTCTAACAGTAACACTACTCAAGAGACCCTGGAAATATGAAAGAATCAGAAAAAACTGGTGGAAG  
AATCTGTAAACAAAAACAAGTTTATCTAAGACTCCAAGTAAGGAAGAAATTGAGAAAGAATGTGAAGA  
TACCAGTTTTCGCTCAGGAGACACAGAGGCGGACATCTAACCATGGTCATGCCAGGAAAAGAGCCAAGTCT  
AATTCCAAGCTAAAGTTGGTGCCTAGCCTGGCAGTGTGTGAGGAGTCCCTCCACCCCATTTGCTGATGGGC  
CATTAGAAACCCAGGATATAATTCAATTGCACATCAGTTGCCCTTCTGACAAGGAGGAAGAAAGTCCAC  
AAAAGATGTCTCTGAAAAGGAAGACAAGGACAAAAACAAGAAAAGATCCCAAGGAAGATGCTGTCCAGA  
GACTCCAGCCAGGAATATACGGAAGTCCACTGGAATAGACCTACATGAATTTCTGTAAATACACTGAAAA  
AGAACCAAGGGACAGAATGATGCTGCTAAAATTAGAACAGGAGATTCTGGAATTTATTAATGACAACAA  
TAATCAGTTCAAGAAGTTCCCTCAGATGACCTCATATCACCGGATGCTATTACACCGGGTAGCTGCCTAT  
TTTGGGATGGACCACAATGTTGATCAAACTGGGAAAGCTGTATCATCAACAAAAGTAAACACAAAGAA  
TCCCTGAACAGAGGTTCTCAGAACATATAAAGGATGAGAAGAATACAGAATTTCAACAGAGGTTCAATTCT  
CAAGAGAGATGATGCCAGTATGGACCGAGATGATAACCAGACTGGCCAGAACGGATATCTAAATGACATC  
AGACTCTCCAAAGAAGCCTTTTCTTCTAGCTCTCACAAGAGAAGGCAGATTTTATAGGGGAACCGTGAAG  
GACTGAGCCGCACCTCAAGCAGCCGCCAGAGCAGCAGACAGCAGCAACTCAAAATCCCTGGAGCCACGCC  
TTGGAGCAGCAGACTCTGATGGCTCTGTCCGGAGCATGCGACCCCTGTACCAAAGCTAGCAGCTTC  
AGTGGAATCTCTATCCTTACCCGAGGTGACAGCATCGGCAGCAGTAAGGCGGCAGTGCAGGGAAGGATCT  
CCAGGCCAGGTATGGCACTAGGTGCCCGAGAAGTGTGAACAGGTACCTCATCCAGTCTGTCCGGGG  
GCTTCTCCCTTGTACTGCCAGCAGCAACAGCAGCAGCAGCAGCAACTTCTGTCTCTCCACCCACAG  
CCTCAGCAACAGCCACCCTTGAATAATCACATGATCTCACAGGCAGATGACCTCAGCAACCCCTTTGGAC  
AAATGAGCCTTAGTCGCCAAGGTTCTACTGAAGCAGTGACCCATCTGCAGTCTATTCCAGACCCCACT  
TATCTCCAGCACCCCTCAGCAGACTAGCTTCATCATGGCTTCCACGGGTGAGCCCTCCCACTTCCAAC  
TATTCACCTCTAGCCATGCACCACCTACTCAGCAAGTTCTGCCACCCAGGGGTACATGCAGCCCTCTC  
AACAGATCCAGGTTTCTTACTATCCCTGGACAATATCCTAACTCCAACCAGCAATATCGACCTCTCTC  
TCACCCGGTGGCTATAGCCCCAACGTGGTGCAGCAGTGCCTCAGCCATCCCAGCAGCCTGGTTTACAG  
CCCATGATGCCTAACAGCAGCAGGCGGCTTACCAAGGCATGATTGGGTCCAGCAGCCACAGAACCAGG  
GCCTGCTCAGCAGCCAGAGGAGCAGCATGGGGGGCCAGATGCAAGGCCTGGTGGTTCACTACACTCCACT  
GCCTTCTTACCAAGTTCAGTGGGTAGTGACTCGCAAAATGTGGTCCAGCCGCTTTCCAGCAACCCATG  
CTGGTCCCTGTGAGCCAGTCTGTGCAAGGAGGCCTCCAGCAGCGGGGTACCAAGTGTACTATAGCATGA  
TCCACCTGCTCAGCAGAACGGTACGAGCCCTTCTGTAGGGTTTCTGCAACCCCTGGCTCTGAGCAGTA  
CCAGATGCCTCAGTCTCCCTCTCCCTGCAGTCCACCACAGATGCCACAGCAGTACTCAGGAGTGTACCT  
TCTGGACCAGGTGTAGTGGTATGCAGCTGAATGTCCCTAATGGACCCAGCCCTCAGAACCCATCCA  
TGGTCCAGTGGAGTCATTGTAATACTACAGCATGGACCAGCGGGGCAGAAGCCTGGAGACCTGTACAG  
TCCTGACAGCAGCCCCAGGCCAACACAAAATGAGCAGCAGCCCTGTACATCTCTACCCAGTCTCCA  
GCACCTCTCCTGTACCAAGCCTCAGCAGTGTCTGCACAGGACTCAGTCCCTGCCTGTCTCACACAGT  
TCCCCCGGCTGGGGTCCAGCACAGGGTGTGGGCGTACTCCCTTTTGGGCCAGCCATTACAGTACAA  
TCTGTCCATCTGCCCTCCCTTGTCTCATGGCCAGTCACTTACACGGTGCACCAGGGACAGAGTGGACTG  
AAGCATGGAACCGGGGCAAGAGACAAGCACTCAAATCTGCCTCCACTGACCTGGGGACAGCAGATGTTG  
TCCTGGGGCGGGTGTGGAGGTGACAGATCTCCCTGAGGGCATCACCCGTAAGGAGCGGACAAACTCTT  
CACGCAGCTCGCCATGTCTGGCGCCAAGATCCAGTGGCTCAAGGATGCTCAGGGGCTGCCTGGAGGGGT  
GGGGGGGACAACAGTGGGACTGCTGAGAATGGCCGCCACTCGGACCTCGCTGCCTTGTACACCATTTG  
CTGTGTTCCCGAGCCCTGGCTGCCAAAATGCCTCCCTTGTCTCAACAACTCCGTGAGTCGCTTCAA  
ACTTCAATGGCAAAAAGAACTATGACCTGAGGATCTGGAGCGAGCCAGTCCCAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG207351 representing NM\_014925

Red=Cloning site Green=Tags(s)


MSNSNTTQETLEIMKESEKKLVEESVNKNKFISKTPSKEEIEKECEDTSLRQETQRRTSNHGHARKRAKS  
NSKLLVRS LAVCEESSTPFADGPLETQDIIQLHISCPDSEEEKSTKDVSEKEDKDNKEKIPRKMLSR  
DSSQEYTDSTGIDLHEFLVNTLKKNPRDRMMLLKLEQEILEFINDNNNQFKFPQMTSYHRMLLHRVAAY  
FGMDHNVDQTGKAVIINKTSNTRIPERFSEHIKDEKNTFQQRFILKRDDASMDRDNQTQNGYLNDI  
RLSKEAFSSSSHKRRQIFRGNREGLSRTSSSRQSTDSELKSLEPRWSSDSDGVSVMRPPVTKASSF  
SGISILTRGDSIGSSKGGSGAGISRPGMALGAPEVCNQVTSSQSVRGLLPCTAQQQQQQQQQLPALPPT  
PQQQPPLNNHMSQADDLSNPFQMSLSRQGSTEAADPSAALFQTLISQHPQQT SFIMASTGQPLPTS  
YSTSSHAPPTQQVLPQGYMQPPQIQVSYYPGQYPNSNQYRPLSHPVAYSPQRGQQLPQPSQQPGLQ  
PMMPNQQAAYQGMIGVQPQNQGLSSQSSSMGGMQGLVVQYTPLPSYQVPVGSQNVVQPPFQQPM  
LVPVSQSVQGLPAAGVPVYYSMIPPAQQNGTSPSVGFLQPPGSEQYQMPQSPSPCSPQMPQYSGVSP  
SGPGVVVMQLNVPNGPQPQNPSMVQWSHCKYYSMQDQKQPGDLYSPDSSPQANTQMSSPVTSTPTQSP  
APSPVTSLSVCTGLSPLPVLTFPRPGPAQGDGRYLLGQPLQYNLSICPPLLHGQSTYTVHQQSGL  
KHGNGRGRQALKSASTDLGTADVVLGRVLEVDLPEGITRTEADKLFTQLAMSGAKIQWLKDAQGLPGG  
GGDNSGTAENGRHSDLAALYTIIVAVFPSPLAAQNASRLNNSVSRFKLRMAKKNYDLRILERASSQ

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:




ACCN:

NM\_014925

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ORF Size:</b>              | 2928 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>                                                                      |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <a href="#">NM_014925.5</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>RefSeq Size:</b>           | 3974 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 2931 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 22864                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <a href="#">Q9Y2K5</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 12q13.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Domains:</b>               | R3H                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |

Product images:



Circular map for RG207351