

## Product datasheet for **RG207407**

### CABCOC01 (NM\_173554) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** CABCOC01 (NM\_173554) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** CABCOC01  
**Synonyms:** C10orf107  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG207407 representing NM\_173554  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGATTCTCTATTATTCAGTATTCAAAAATTTGACTTTACTAGCTATGTCACCTTCAAATCTTAAAA  
 CACTACATATGTCCTTAGAGGAAAGCATAAAATGGCTTGGAGAAGTTATGGCTGAAATAGGACCAACACA  
 TTCGAAAAGAGTGAGGACTGGAATATCTTTGATGTAAAAACAAGCCAATGCTATCATTGATTACTTAAAA  
 ATCAGCTATTTCAACTACAAGCTATACGAGTTTATGTTCTACTCTGCCAGGGAAGAAATTTGTGATAG  
 GAACTGAGCAAGTGATAGAGGTTGTCAAGTCTGCATGTGGCCCTTCCCAAATCCTCTGGAAGAAGGAAT  
 CTATTTGATATTTATCAACATTCATAGAGCCCCCACAATATTGGATACGGAAATGAAGAGGTTGGAT  
 CAAGAACAAGGCCCTGAGGAGTCTCAGCCAGAAACTGACACCTCAGACATGGATCCTTTAGTTGGTTTCA  
 CCATTGAAGATGTAAAATCAGTGCTGGATCAAGTCACAGATGACATTTTAAATCGGCATTCAGACCGAGAT  
 AAACGAAAAACTGCAAATACAGGAAGAGGCCTTTAATGCACGAATAGAAAAATTGAAAAAGGCC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >RG207407 representing NM\_173554  
 Red=Cloning site Green=Tags(s)

MDFSIIQYSKFMILLAMSLQNLKTLHMSLEESIKWLGEVMAEIGPHTSQKSEDWNIFDVKQANAIIDYLK  
 ISLFQHYKLYEFMFYSAREEIVIGTEQVIEVVKSSACGPFNPLEEIGISFDIYSTFIEPPTILDTEMKRLD  
 QEQGPESQPETDTSMDPLVGFITIEDVKSVLQVTDLILIGIQTEINEKLIQEEAFNARIEKLLKA

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_173554

**ORF Size:** 624 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_173554.1](#), [NP\\_775825.1](#)

**RefSeq Size:** 1301 bp

**RefSeq ORF:** 627 bp

**Locus ID:** 219621

**UniProt ID:** [Q8IVU9](#)

**Cytogenetics:** 10q21.2

**Gene Summary:** Calcium-binding protein. May be involved in the control of sperm flagellar movement. [UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for RG207407