

Product datasheet for **RG207478**

C5orf54 (ZBED8) (NM_022090) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	C5orf54
Synonyms:	Buster3; C5orf54
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG207478 representing NM_022090
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTCGAAGAAACGCAAAATGGGATGATGACTATGTTTCGCTACTGGTTCACCTGTACAACGGAAGTAGATG
 GAACTCAGCGTCCACAGTGTGTGTTGTAACTCAGTATTTTCAAATGCTGACCTCAGACCATCAAACT
 GTCAGACCATTTTAAACAGACAGCATGGTGGTGTAGCTGGGCATGATCTCAATAGCCTGAAGCATATGCCA
 GCACCATCTGATCAGAGTGAAACCTTGAAAGCATTGGAGTTGCATCTCATGAGGATACTTTATTACAAG
 CATCGTATCAATTTGCGTATTTATGTGCCAAGGAGAAGAATCCTCATACAGTAGCTGAAAAGTTAGTGAA
 ACCTTGTGCACTGGAAATAGCACAAATAGTTTTGGGACCAGATGCACAAAAAAGCTTCAGCAGGTACCC
 TTATCAGATGATGTGATCCATTCTAGAATTGATGAAATGAGCCAGGATATCTTACAGCAAGTTCTAGAAG
 ATATCAAAGCCAGTCTCTTAAAGTGGTATTGAGCTTGCTGAGACAACTGACATGGATGACTGCAGTCA
 GCTAATGGCATTGTGCGCTATATAAAGAAAGAGAGATCGTAGAAGAATTTCTTCTGTGAACCATTG
 CAGCTATCCATGAAAGGAATAGATGTGTTCAATCTCTTCAAGAGACTTCTTCTGAAGCATAAGATAGCAC
 TTGATGTATGTGGCTCTGTTTGTACTGATGGTGCCTCCTCTATGCTAGGAGAAAATTCGAGTTTGTTC
 CTACGTGAAAAAGAGATACCTCATATCGTAGTTACACATTGTTTATTGAATCCTCATGCACTTGTCTATA
 AAGACATTGCCTACAAAAGTGGGATGCTCTGTTTACTGTGGTGAGGGTAATAAATTTATCAAAAGGGA
 GAGCTCCAAATCATCGCCTATTTGAGGCTTTCTTTGAAGAAATTGGCATAGAATATAGTGTCTCCTCTTT
 CCATACTGAAATGAGGTGGCTTTCCCGAGGCCAAATACTTACTCACATTTTGAATGTATGAAGAAATA
 AATCAGTTTCTTACCACAAGAGCAGTAATTTAGTTGATGGCTTTGAAAATAAGAGTTTAAATTCACC
 TAGCATACCTTGCAGATTTATCAAACTTAAATGAAGTGCATCTATGCAGAGGACTGGGATGAA
 CACAGTATCAGCTAGAGAGAAGTTATCTGCTTTTGTAGGAAGTTTCCATTTTGGCAAAAACGAATTGAG
 AAAAGAAATTTTACCAATTTTCTTTTCTTGAAGAAATAATGTTTCAGATAATGAAGGCATATTCATTG
 CAGCTGAAATAACACTGCATCTGCAACAATTGAGCAACTTCTTCCATGGATATTTTCCATTGGAGATCT
 TAATGAGGCAAGTAAATGGATATTGGATCCATTCTTTTAAATATCGATTTTGTGATGATAGTTATTTA
 ATGAAAAATGATCTTGCTGAATTACGAGCTAGTGCCCAATCCTAATGGAATTTGAGACAATGAAGCTTG
 AGGATTTCTGGTGTGCCAATTCACAGCATTTCCAAACCTGGCAAAGACAGCTCTAGAAATCCTTATGCC
 ATTTGCAACTACATACCTTTGTGAGTTGGGATTTTATCACTTTTACATTTCAAACAAAGTCCAGAAGC
 TGCTTTAATCTGAGTGATGATATCCGTGTGGCTATTTCAAAAAAGTTCCTCGTTTCTCGGACATCATTG
 AACAAAGCTACAGCTACAGCAGAAGTCACTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG207478 representing NM_022090
 Red=Cloning site Green=Tags(s)

MSKKRWDDDYVRYWFTCTTEVDGTQRPQCVLCNSVFSNADLRPSKLSDFHNRQHGGVAGHDLNSLKHMP
 APSDQSETLKAFGVASHEDTLLQASYQFAYLCAKEKNPHTVAEKLKPCALEIAQIVLGPDAQKKLQQVP
 LDDVIHSRIDEMSQDILQQVLEDIKASPLKVGIQLAETTDMDDCQLMAFVRYIKEREIVEEFLFCEPL
 QLSMKGIDVFNLFDRDFLKHKIALDVCVSVCTDGASSMLGENSEFVAYVKKEIPHIVVTHCLLNPHALVI
 KTLPTKLRDALFTVVRVINFIKGRAPNHRLFQAFEEIGIEYSVLLFHEMRWLSRGQILTHIFEMYEI
 NQFLHHKSSNLVDGFENKEFKIHLAYLADLFKHLNELSASMQRGMNTVSAREKLSAFVRKFPFWQKRIE
 KRNFTNFPFLEEIIIVSDNEGIFIAAEITLHLQLSNFFHGYFSIGDLNEASKWILDPFLFNIDFVDDSYL
 MKNDLAE LRASGQILMEFETMKLEDFWCAQFTAFPNLAKTALEILMPFATTYLCELGFSSLLHFKTKSRS
 CFNLSDDIRVAISKVPRFSIIIEQKLQLQKSL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_022090

ORF Size: 1782 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_022090.5](#)

RefSeq Size: 2808 bp
 RefSeq ORF: 1785 bp
 Locus ID: 63920
 UniProt ID: [Q8IZ13](#)
 Cytogenetics: 5q33.3

Product images:



Circular map for RG207478