

Product datasheet for **RG207538**

EVI2A (NM_014210) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	EVI2A
Synonyms:	EVDA; EVI-2A; EVI2
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG207538 representing NM_014210
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCCACGGACATGGAACACACAGGACATTACCTACATCTTGCCTTTCTGATGACAACAGTTTTTCTT
TGTCTCCTGGAACAAAAGCAAACCTATACCCGTCTGTGGGCTAACAGTACTTCTCCTGGGATTCAAGTTAT
TCAAAACAAGACAGGCAGAAACCAAATGAAAACATTAAACAAACCTATAACTCCTGAAGTAGATTAT
AAAGGTAATTCTACAAACATGCCTGAAACATCTCACATCGTAGCTTAACTTCTAAATCTGAACAGGAGC
TTTATATACCTTCTGTCTGTCAGCAACAGTCCTTCAACAGTACAGAGCATTGAAAACACAAGCAAAAGTCA
TGGTGAATTTTCAAAAAGGATGTCTGTGCGGAAAACAACAACATGGCTATGCTAATTTGCTTAATT
ATAATTGCAGTGCTTTTCTTATCTGTACCTTTCTATTTCTATCAACTGTGGTTTGGCAAACAAAGTCT
CTTCTCTCAGACGATCAAAACAAGTAGGCAAGCGTCAGCCTAGAAGCAATGGCGATTTTCTGGCAAGCGG
TCTATGGCCCGCTGAATCAGACACTTGGAAAAGAACAAACAGCTCACAGGACCCAACCTAGTGATGCAA
TCTACTGGAGTGCTCACAGCTACAAGGAAAGAAAAGATGAAGAAGGAAGTGAAGAACTTACTAACAAC
AGATAGGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG207538 representing NM_014210

Red=Cloning site Green=Tags(s)


MPTDMEHTGHYHLAFLMTTVFSLSPGKANYTRLWANSTSSWDSVIQNKTRNQENINTNPITPEVDY
KGNSTNMPETSHIVALTSEQELYIPSVVSNPSTVQSIENTSKSHGEIFKKDVCAENNNNMAMLI
IIAVLFLICTFLFLSTVVLANKVSSLRRSKQVGKRQPRSNQDFLASGLWPAESDTWKRTKQLTGPNLVMQ
STGVLATATREKDEEGTEKLTNKQIG

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_014210

ORF Size:

708 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:


This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_014210.4</u>
RefSeq Size:	1622 bp
RefSeq ORF:	711 bp
Locus ID:	2123
UniProt ID:	<u>P22794</u>
Cytogenetics:	17q11.2
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	May complex with itself or/and other proteins within the membrane, to function as part of a cell-surface receptor.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG207538