

Product datasheet for **RG207585**

FCER1G (NM_004106) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: FCER1G (NM_004106) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: FCER1G
Synonyms: FCRG
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG207585 representing NM_004106
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATTCCAGCAGTGGTCTTGCTCTTACTCCTTTGGTTGAACAAGCAGCGGCCCTGGGAGAGCCTCAGC
TCTGCTATATCCTGGATGCCATCCTGTTTCTGTATGGAATTGTCCTCACCTCCTCTACTGTCGACTGAA
GGTAATCCAAGTGCAGAAAGGCAGCTATAACCAGCTATGAGAAATCAGATGGTGTACACGGCCTGAGC
ACCAGGAACCAGGAGACTTACGAGACTCTGAAGCATGAGAAACCACCACAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG207585 representing NM_004106
Red=Cloning site Green=Tags(s)
MIPAVVLLLLLLVEQAAALGEPQLCYILDAILFLYGIVLTLLYCRCLKVIQVRKAAITSYEKSDGVYVYGLS
TRNQETYETLKHEKPPQ

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_004106

ORF Size: 261 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_004106.1](#), [NP_004097.1](#)

RefSeq Size: 591 bp

RefSeq ORF:	261 bp
Locus ID:	2207
UniProt ID:	P30273
Cytogenetics:	1q23.3
Domains:	ITAM
Protein Families:	Transmembrane
Protein Pathways:	Asthma, Fc epsilon RI signaling pathway, Natural killer cell mediated cytotoxicity
Gene Summary:	The high affinity IgE receptor is a key molecule involved in allergic reactions. It is a tetramer composed of 1 alpha, 1 beta, and 2 gamma chains. The gamma chains are also subunits of other Fc receptors. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG207585