

## Product datasheet for **RG207656**

### VSTM2L (NM\_080607) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** VSTM2L (NM\_080607) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** VSTM2L  
**Synonyms:** C20orf102; dj1118M15.2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG207656 representing NM\_080607  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGGGGCCCGCTCGCGTAGCGCTGGGCGCCCTCCACTACCTGGCACTTTTCTGCAACTCGGCGGCG  
CCACGCGGCCCGCCGACGCGCCCTGGGACAACCACGTCTCCGGCCACGCCCTGTTACAGAGACACC  
CCATGACATGACAGCACGGACGGCGGAGGACGTGGAGATGGCCTGCTCCTCCGCGGCAGCGGCTCCCC  
TCCTACTCGTGGAGATCCAGTGGTGTATGTACGGAGCCACGGGACTGGACCGACAAGCAGGCGTGGG  
CCTCGAACCAGCTAAAAGCATCTCAGCAGGAAGACGCAGGAAGGAGGCAACCAAAATAAGTGTGGTCAA  
GGTGGTGGGACAGCAACATCTCCACAAGCTGCGCCTGTCCGGGTGAAGCCACGGACGAAGGCACCTAC  
GAGTGCCGCGTACGACTTCAGCGACGCAAGGCCCGGCACCACAAGGTCAAGGCTACCTGCGGGTGC  
AGCCAGGGGAGAACTCCGTCTGCATCTGCCCGAAGCCCTCCCGCCGCGCCCGCCCGCCCGCCCA  
GCCAGGCAAGGAGCTGAGGAAGCGCTCGGTGGACCAGGAGGCCTGCAGCCTC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >RG207656 representing NM\_080607  
Red=Cloning site Green=Tags(s)

MGAPLAVALGALHYLALFLQLGGATRPAGHAPWDNHVSGHALFTETPHDMTARTGEDVEMACSFRRGSGSP  
SYSLEIQWVYVRSHRDWTDKQAWASNQLKASQQEDAGKEATKISVVKVVGSNISHKLRSLRVKPTDEGT  
ECRVIDFSDGKARHHKVKAYLRVQPGENSVLHLPEAPPAAPPPPKPKGKELRKRSVDQEACSL

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_080607

**ORF Size:** 612 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_080607.3](#)

**RefSeq Size:** 1952 bp


**RefSeq ORF:** 615 bp

**Locus ID:** 128434

UniProt ID: Q96N03

Cytogenetics: 20q11.23

**Product images:**



Circular map for RG207656