

Product datasheet for **RG207818**

CBLN2 (NM_182511) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: CBLN2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG207818 representing NM_182511
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCAGGCGCCGGCCGGGGGCCACTCGGGCTGCGGCTGATGATGCCGGGCGCCGGGGGGCGCTGCGCG
AGCCTGGCGGCTGCGGATCCTGCCTGGGGTGGCGCTGGCCCTGCTGTTGCTGCTACTGCCCGCTGCTG
CCCCGTGCGGGCGCAGAACGACACGAGGCCATCGTGCTGGAGGGCAAGTGCCTGGTGGTGTGCGACTCC
AGCCCGTCGGCGGACGGCGCCGTACCTCCTCCCTAGGCATCTCCGTGCGCTCCGGCAGCGCCAAGGTGG
CCTTCTCCGCCACGCGGAGCACCAACCAGAGCCGTCCGAGATGAGCAACCGCACCATGACCATCTATTT
CGACCAGGTATTAGTAAATATTGGCAACCACTTTGATCTTGCTTCCAGTATATTGTAGCACCGAGAAAA
GGGATTTATAGCTTCAGCTTCCACGTGGTCAAAGTGTATAACAGACAAACCATCCAGGTCACTTTAATGC
AGAATGGCTACCCAGTGATCTCGGCCTTTGCAGGAGACCAGGATGTACCCAGAGAAGCTGTAGCAATGG
CGTGCTGCTGCTCATGAAAGGGAAGACAAAGTGCATCTCAAACCTGAGAGAGGCAACCTCATGGGGGGC
TGAAATACTCCACATTCTCGGGCTTCTTGGTGTTCCTCTA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG207818 representing NM_182511
Red=Cloning site Green=Tags(s)

MQAPGRGPLGLRLMMPGRRGALREPGGCGSCLGVALALLLLLLPACCPVRAQNDTEPIVLEGKCLVVCDS
SPSADGAVTSSLGISVRSKSAKVAFSATRSTNHEPSEMSNRTMTIYFDQVLVNIGNHFDLASSIFVAPRK
GIYSFSFHVKVYNRQTIQVSLMQNGYPVISAFAGDQDVTREASNGVLLLMEREDKVHLKLERGNLMGG
WKYSTFSGLVFPL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




Cloning Scheme:



ACCN:	NM_182511
ORF Size:	672 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_182511.2 , NP_872317.1
RefSeq Size:	2750 bp

RefSeq ORF:	675 bp
Locus ID:	147381
UniProt ID:	Q8IUUK8
Cytogenetics:	18q22.3
Gene Summary:	May play role in synaptogenesis induction.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for RG207818