

Product datasheet for **RG207836**

EBF1 (NM_024007) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EBF1 (NM_024007) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	EBF1
Synonyms:	COE1; EBF; O/E-1; OLF1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG207836 representing NM_024007
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTTTGGGATTCAGAAAGCATCCAACGGAGTGAAGCAGCATGAAGGAAGAGCCGCTGGCGACGGCA
 TGAAACGGGTGCGGACGTGGATGCAAGGCGCCGGGTGCTGGACGCCAACACGGCGCGCAGAGCGGGT
 GGGTCTGGCCCGGCTCACTTTGAGAAGCAGCCGCTTCCAATCTGCGGAAATCCAACTTCTTCCACTTC
 GTCCTGGCCCTACGACAGACAGGGCCAGCCGTGGAGATCGAGAGGACAGCGTTTGTGGGGTTCGTGG
 AGAAGGAAAAAGAAGCCAACAGCGAAAAGACCAATAACGGAACTCACTACCGGCTTCAGTTCTCTACAG
 CAATGGGATAAGGACGGAGCAGGATTTCTACGTGCGCCTCATTGACTCCATGACAAAACAAGCCATAGTG
 TATGAAGGCCAAGACAAGAACCAGAAATGTGCCGAGTCTTGCTCACACATGAGATCATGTGCAGCCGCT
 GTTGTGACAAGAAAAGCTGTGGCAACCGAAATGAGACTCCCTCAGATCCAGTGATAATTGACAGGTTCTT
 CTTGAAATTTTTCTCAAATGTAACAAAATTGCCTAAAGAATGCGGGAAACCCAGTGACATGCGGAGA
 TTCCAGGTCGTGGTGTCTACGACAGTCAATGTGGATGGCCATGTCTTGGCAGTCTCTGATAACATGTTTG
 TCCATAATAATTCCAAGCATGGGCGGAGGGCTCGGAGGCTTGACCCCTCGGAAGGTACGCCCTTATCT
 GGAACATGCTACTCCCTGTATCAAAGCCATCAGCCGAGTGAAGGATGGACGACGGGAGGTGCGACTGTG
 ATCATCATAGGGGACAATTTCTTTGATGGGTTACAGGTCAATTCGGTACCATGCTGGTCTGGAGTGAGT
 TGATCACTCCTCATGCCATCCGTGTGCAGACCCCTCCTCGGCACATCCCTGGTGTGTGGAAGTCACACT
 GTCCTACAAATCTAAGCAGTTCTGCAAAGGAACACCAGGCAGATTCAATTTACAGCGCTCAACGAACCC
 ACCATCGATTATGGTTTCCAGAGGTTACAGAAGGTCATTCTCGGCACCCCTGGTGACCTGAGCGTTTGC
 CAAAGGAAGTAATACTCAAAGGGCTCGGATCTGGTAGAAGCACTGTATGGGATGCCACACAACCA
 GGAAATCATTCTGAAGAGAGCGGCCGACATTGCCGAGGCCCTGTACAGTGTCCCCGCAACCACAACCAA
 CTTCCCGCCCTTCTAACACCTCGGTCCACGACGGATGATGGGCGTGAATTCGTTCACTGGACAACCTGG
 CCGTGAATGTCTCCGAGGCATCACAAGCCACCAATCAGGGTTTCAACCGCAACTCAAGCAGCGTATCACC
 ACACGGGTACGTGCCGAGCACCCTCCCGCAGACCAACTATAACTCCGTCACCACGAGCATGAACGGA
 TACGGCTCTGCCCAATGTCCAATTTGGGCGGCTCCCCACCTTCTCAACGGCTCAGCTGCCAATCCC
 CCTATGCCATAGTCCATCCAGCCCCACCATGGCTCCTCCACAAGCCTCCCTCCAACCTGCAGCAGCTC
 CTCGGGCATCTTCTCTCTCACCAGCCAACATGGTCTCAGCCGTGAAACAGAAGAGTGCTTTCGACCA
 GTCGTGACACCCAGACCTCCACCTCCACCTGCACCAGCACAACGGGAACAGCCTGCAAGCGATAT
 CTGGCATGATTGTTCTCTATG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG207836 representing NM_024007
 Red=Cloning site Green=Tags(s)

MFGIQESIQRSGSSMKKEPLGSGMNAVRTWMQGAGVLDANTAAQSGVGLARAHFEKQPPSNLRKSNFFHF
 VLALYDRQGQPVEIERTAFVGFVEKEKEANSEKTNNGIHYRLQLLYSNGIRTEQDFYVRLIDSMTKQAI
 YEGQDKNPEMCRVLLTHEIMCSRCCDKKSCGNRNETPSDPVIIDRFFLKFLLKCNQNLKNAGNPRDMRR
 FQVVVSTTVNVDGHVLAUSDNMFVHNSKHGRRARRLDPSEGTPSYLEHATPCIKAI SPSEGWTTGGATV
 I IIGDNFFDGLQVIFGMLVWSELITPHAIRVQTPPRHIPGVVEVTL SYKSKQFCKGTPGRFIYALNEP
 TIDYGFQRLQKVI PRHPGDPERLPKEVILKRAADLVEALYGMPHNNQEIILKRAADIAEALYSVPRNHQ
 LPALANTSVHAGMMGVNSFSGQLAVNVSEASQATNQGFTRNSSVSPHGYVPSTTPQQTNYNSVTTSMNG
 YGSAAMSNLGGSPFTFLNGSAANSPYAI VPSSTMASSTSLPNSCSSSSGIFSFSPANMVS AVKQKSAFAP
 VVRPQTSPPPTCTSTNGNSLQAI SGMIVPPM

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_024007

ORF Size: 1773 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_024007.5](#)


RefSeq Size: 2762 bp

RefSeq ORF: 1776 bp

Locus ID: 1879

UniProt ID: [Q9UH73](#)
Cytogenetics: 5q33.3
Protein Families: Transcription Factors
Gene Summary: Transcriptional activator which recognizes variations of the palindromic sequence 5'-ATTCCCNNGGGAATT-3'. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG207836