

Product datasheet for **RG207852**

CCNJL (NM_024565) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	CCNJL
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG207852 representing NM_024565
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGATGGATGAGCCGTGGTGGGAAGGGCGCGTCGCTCGGACGTCCACTGCACCCTGCGCGAGAAGGAAC
TGAAGTGTCCACCTTCGAGCCCACTCCCACTCCTGAAGAGCCGCGTTCTTCGTGGACATCCTGAC
CCTGCTGAGCAGCCACTGCCAGCTCTGCCCTGCAGCCCGGCACCTGGCCGTCTACCTGCTGGACCACTTC
ATGGATCGCTACAACGTCAACCACTCCAAGCAGCTCTACACCGTGGCCGTCTCCTGCCTCCTGCTTGCAA
ACGGAGTTTCGCTCTTATCGCCAGGCTGAAGTGCAGTGGCATGATCTCAGCTCATTGCAACCTCCACCT
CCCAGGTTCAAGCAATTCTCCTGCTTCAGCGCCCCACCCCCACCCACACCGCCCCAAGTAGCTGAGACT
ACAGGTAAGTTTCGAGGATCGGGAAGACCACTCCCAAGTTGGAGCAAAATAACAGCACGAGGATCCTGA
GCAGCCAGAACTTCACCTCACCAAGAAGGAGCTGCTGAGCACAGAGCTGCTGCTCCTGGAGGCCTTCAG
CTGGAACCTCTGCCTGCCACGCTGCCCACTTCCTGGACTACTACCTCTTGGCCTCCGTGAGCCAGAAG
GACCACCACTGCCACACCTGGCCACCACTGCCCCGCAAGACCAAGAGTGCCTCAAGGAGTATGCC
ATTACTTCCTAGAGGTACCCCTGCAAGATCACATATTCTACAAATTCAGCCTTCTGTGGTCTGCTGCGGC
CTGTGTTGGGCTCCAGGATTTGCCTGCAGCTTTCTCCCTACTGGACCAGAGACCTGCAGAGGATCTCA
AGCTATTCCTGGAGCACCTCAGCAGTGTATTGAAATCCTGCTGGTAGTGTATGACAACGTCCTCAAGG
ATGCCGTAGCCGTCAAGAGCCAGGCTTGGAATGGTGCCCGGCACACCCCCACCCCACTCAAGTGCT
GTTCCAGCCACGACCTACCCGGCCCTCGGCCAGCCAGCAGCACCCTGGCACAGTTCCAGACCCCGTG
CAGGACCTATGCTTGGCTATCGGGAATCTTGCAGGCCACCGTTGAGGAGCCTGCTCTCGGGGAGTA
CAGGCTCATCCCTCCACACCCGTAACAACCGCTGCAGCCCTTGGATATGTGTCCCGTGGCCGTCCCTGC
ATCCCTTAGCATGCATATGGCAATTCAGCTGAGCCAGGCACTGCCTCGCCACCACCTATGGAAGCAGC
TACTTCAGTGGGAGCCACATGTTCCCAACCGGCTGCTTTGACAGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG207852 representing NM_024565
Red=Cloning site Green=Tags(s)

MMDEPWWEGRVASDVHCTLREKELKLPTRAHSPLLKSRFFVDILTLLSSHQCLCPAARHLAVYLLDHF
MDRYNVTTSKQLYTVAVSCLLLANGVSLSPRLKCSGMISAHCNLHLPSSNSPASAPHPPTPPQVAET
TGKFEDREDHVPKLEQINSTRILSSQNFTLTKEKLLSTELLLEAFSWNLCLPTPAHFLDYLLASVSQK
DHHCHTWPTTCPRKTKECLKEYAHYFLEVLQDHIFYKFQPSVAAACVGASRICQLSPYWTRDLQRIS
SYSLEHLSTCIEILLVVDNVLKDAVAVKSQALAMVGPPTPTQVLFQPPAYPALGQPATTLAQFQTPV
QDLCLAYRDSLQAHRSGLLSGSTGSSLHTPYQPLQLDMCPVPVPASLSMHMAIAAEPRHCLATTYGSS
YFSGSHMFPTGCFDR

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_024565

ORF Size: 1305 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_024565.7](#)

RefSeq Size: 3274 bp


RefSeq ORF: 1308 bp

Locus ID: 79616

UniProt ID: [Q8IV13](#)

Cytogenetics: 5q33.3

Product images:



Circular map for RG207852