

## Product datasheet for **RG207855**

### ZCCHC7 (NM\_032226) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZCCHC7 (NM_032226) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZCCHC7
Synonyms:	AIR1; HSPC086
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RG207855 representing NM\_032226  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGTTTGGTGGCTATGAGACTATAGAAGCATACGAAGATGATCTTTATCGAGATGAGTCATCTAGTGAAC  
 TGAGTGTGATAGTGAGGTGGAATTTCAACTCTATAGCCAAATCATTATGCCAAAGATCTTGATGATGT  
 CATCAGGGAGGAAGAGCATGAAGAAAAGAACTCTGGGAATTCGGAATCTTCGAGTAGTAAACCAATCAG  
 AAGAAGCTAATCGTCCTTTAGATAGTGAGGTCATCCAGCTGTCAGATGGGTGAGAGGTCATCACTTTGT  
 CTGATGAAGACAGTATTTATAGATGTAAGGAAAGAATGTTAGAGTTCAAGCACAAAGAAAATGCCCATGG  
 TCTTTCTTCTTCTTCAATCTAATGAGCTGGTTGATAAGAAATGCAAGAGTGATATTGAGAAGCCTAAA  
 TCTGAAGAGAGATCAGGTGTAATCCGAGAGGTCATGATTATAGAGGTGAGTCAAGTGAAGAGGAAGAGA  
 GCACCATTTCAGAAGGTGATAATGTGAAAGCTGGATGCTACTGGGATGTGAAGTAGATGATAAAGATGA  
 TGATATCCTTCTCAACCTTGTGGGATGTGAAAACCTCTGTTACTGAAGGAGAAGATGGTATAAACTGGTCC  
 ATCAGTGACAAAGACATTGAGGCCAGATAGCTAATAACCGAACACCTGGAAGATGGACCCAGCGGTACT  
 ATTCAGCCAACAAAAACATTATCTGTAGAAATGTGACAAACGTGGTCATTTATCAAAAACTGCCCCTT  
 ACCACGAAAAGTTCGTCGCTGCTTCTGTGCTCCAGGAGAGGACATCTCCTGTATTCTGTCCAGCCCC  
 CTTTGGCAATACTGTCCTGTGCCTAAGATGTTGGACCACTCATGTCTTTTCAGACATTCCTGGGATAAAC  
 AGTGTGACCGATGTCATATGCTAGGCCACTATACAGATGCTTGCACAGAAATCTGGAGGCAGTATCACCT  
 AACGACCAACCTGGACCACCAAAAAGCCGAAGACCCCTTCAAGACCATCAGCCTTAGCATATTGCTAT  
 CACTCGCGCAAAAAGGCCATTATGGACACGAATGTCCAGAAAGAGAAGTGTATGACCCGTCTCCAGTAT  
 TCTCCATCTGCTACTATGATGACAAATATGAAATTCAGGAGAGAGAAAAGAGACTAAAACAAAAAT  
 AAAAGTACTCAAGAAAAATGGGGTTATCCAGAGCCATCCAAGTACCTTATATAAAAGCAGCAAAATGAG  
 AACCCCAACCATGATATAAGGAAGGGCCGTGCCTCATGAAAAGCAACAGGTGGCCTCAAGAAAAATAAG  
 AAACACAAAAAGAAATGAAGAACAAGAATAGAACTGGGAGAAGCACAGGAAGGCTGACAGACATCGTGA  
 AGTGGATGAGGATTTTCCAGGGGCCCAAAACCTACTTCTCCTGGCAGTTTTAAACCAGAAGCCT  
 TCTAAGCCCTTTCACCGTTCATCACATTACCACACGTCAAGAGAAGACAAGTCTCCAAGGAAGGCAAGA  
 GGGCAAGCAGAAGAAAAGGAGAGGTGCTGGGAAGATGATGACAATGATAACTATTTCTATTAAGCA  
 GAGAAAAAAAAGTCT

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:**

>RG207855 representing NM\_032226  
 Red=Cloning site Green=Tags(s)

MFGGYETIEAYEDDLRDESSSELVDSEVEFQLYSQIHYAQDLDDVIREEEHEEKNSGNSSESSSKPNQ  
 KKLIVLSDSEVIQLSDGSEVITLSDEDSIYRCKGNVRVQAQENAHGLSSSLQSNELVDKCKSDIEKPK  
 SEERSGVIREVMIIEVSSSEEEESTISEGDNVESWMLLGCEVDDKDDILLNLVGCENSVTGEDGINWS  
 ISDKDIEAQIANNRTPGRWTQRYYSANKNIICRNCDKRGHLSKNCPLPRKVRRFCFSRRGHLLYSCAP  
 LCEYCPVKMLDHSCLFRHSWDKQCDRCHMLGHYTDACTEIWRQYHLTTKPGPPKPKTPSRPSALAYCY  
 HCAQKGHYGHECPEREVYDPSVSPFICYDDKYEIQEREKRLKQKIKVLKKNGVPEPSKLPYIKAANE  
 NPHHDIRKGRASWKSNRWPQENKETQKEMKNKRNWEKHRKADRHREVEDDFPRGPKTYSSPGSFKTQKP  
 SKPFHRSSHYSREDKSPKEGKRGKQKKKERCWEDDDNDNLFLIKQRKKKS

**TRTRPLE** - GFP Tag - V

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_032226

**ORF Size:** 1626 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_032226.1](#), [NP\\_115602.1](#)

**RefSeq Size:** 2701 bp


**RefSeq ORF:** 1632 bp

**Locus ID:** 84186

**UniProt ID:** [Q8N3Z6](#)

Cytogenetics: 9p13.2  
 Domains: zf-CCHC  
 Protein Pathways: RNA degradation

**Product images:**



Circular map for RG207855