

Product datasheet for **RG207865**

LIX1L (NM_153713) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	LIX1L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG207865 representing NM_153713
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGACTATGCGAGCGCAGCGGCTGCAGCCTGGTGTGGGCACCAGCGGGAGGGGCACTCTCCGAGCGC
TGCGGCCCGGAGTGACTGGGGCCGCGGTGCCACCGCCACACCCCTGCGGGCCCCCGCTGCCCGCC
GCCTCCCGCACCGCCCCGCGCGCTGCTCCTGTCTGGGGCCCCAGGACTACCCCTGCCCGCCGCGCC
GCCGGCAGCCCGGCACTGCTGCGAGAGGCCGTGGAGGCCGTGGTGAAGGAGCTTCGCCAAGCACACGAGG
GCTATGGCCGAGTGAATGTGGTGGAGGCACTTCAGGAATTCGGCAGATGAAGCAGTCCCGTGGTGTGA
CTTAAGAATGGGGCTCTAGTGGTTATGAGATGGTCCCTCCAACAGCCCTCCTTATGCTGTATGTC
ACCCTGCCTGGGGGAAGCTGCTTTGGGAGTTTCCAGTTTGGCCCCACAAAGCTGAGGCCCGGAGGAGTG
CTGCAAAGATTGCGCTAATGAATTCTGTGTTAATGAACATCCTTCCCGAAGAATCACTGATGAGTTCAT
CGAGAAGAGTGTCTCTGAGGCCCTGGCATCTTTAATGGCAACAGGGAGGAAGCTGACAACCCAAATACA
GGGATTGGTGCCTTCGATTTCATGCTGGAATCCAACAAGGGCAAAATCAATGTTGGAGTTCAGGAGCTAA
TGACAGTTTTTCACTGCTACACTGGAATGGCAGCCTTAAGGCCATGAGGGAACGACAATGCTCTCGGCA
GGAGGTGTTGGCTCATTATTCGACCGGGCCCTGGATGATGATTCGCCACCAAAATGGCCTTGACTGG
GTGAGCCGGGAGCAGAGTGTGCCGGGGCACTGTCTAGAGAGCTGGCCTCTACTGAGCGGGAGCTGGATG
AAGCCCGACTGGCAGGCAAGGAGCTGCGCTTCCACAAGGAGAAGAAAGATATTCTGTGCTGGCTGCTGG
GCAGTTGGCAATATGCATTCTTCCAATGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >RG207865 representing NM_153713
 Red=Cloning site Green=Tags(s)

METMRAQLQPGVGTSGRGTLRALRPGVTGAAAATATPPAGPPPAPPPAPPPPLLLSGAPGLPLPPGA
 AGSPAVLREAVEAVVRSFAKHTQGYGRNVVEALQEFWQMKQSRGADLNKALVYEMVPSNSPPYVCYV
 TLPGGSCFGSFQFCPTKAEARRSAKIALMNSVFNEHPSRRITDEFIEKSVSEALASFNGNREEADNPNT
 GIGAFRFMLESNGKSMLEFQELMTVFQLLHWNGSLKAMRERQCSRQEVLAHYSHRALDDDIRHQMALDW
 VSREQSVPGALSRELASTERELDEARLAGKELRFHKEKKDILVLAAGQLGNMHSSNC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_153713

ORF Size: 1011 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_153713.1](#), [NP_714924.1](#)

RefSeq Size: 1467 bp


RefSeq ORF: 1014 bp

Locus ID: 128077

UniProt ID: [Q8IVB5](#)

Cytogenetics: 1q21.1

Product images:



Circular map for RG207865