

Product datasheet for **RG207903**

GTPBP4 (NM_012341) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	GTPBP4
Synonyms:	CRFG; NGB; NOG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG207903 representing NM_012341
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCACCATTACAACCTCAAGAAAATTACGGTGGTGCCGTCGCCCAAGGACTTCATAGACCTCACGTTGT
 CGAAGACTCAACGAAAGACTCCAACCGTTATTCATAAACATTACCAAATACATCGCATTAGACATTTTTC
 CATGAGAAAAGTCAAAATTTACTCAACAGAATTACCATGATAGACTTTTCACAAATTTCTAACAGATTTCCCC
 AAATTGGATGATATTCATCCGTTCTATGCTGATTTGATGAATATTCTCTACGACAAGGATCATTACAAGT
 TGGCTCTGGGGCAAATAAATATTGCCAAAAATTTAGTGACAATGTTGCTAAAGATTATGTGCGACTGAT
 GAAGTATGGCGACTCTCTCTACCGCTGCAAAACAGCTGAAGCGTGCAGCGCTGGGACGGATGTGCACAGTG
 ATCAAGAGGCGAGAAGCAGAGTTTGGAGTATTTGGAGCAAGTGCCTCAGCATTTATCCCGTTTGCCAAACCA
 TTGATCCGAATACCGAGCCCTGCTTTTGTGGGTACCCAAATGTTGGGAAGTCCAGCTTCATCAACAA
 GGTGACGAGAGCAGACGTGGATGTCCAGCCCTATGCGTTACAACCAAGTCTCTGTTTGTGGGCACATG
 GATTATAAGTATCTACGTTGGCAGGTTGTAGACACTCCTGGGATCCTGGACCACCTCTGGAGGATAGGA
 ACACCATCGAGATGCAGGCCATCACTGCCCTGGCCACCTCCGTGCTGCAGTCTGTATGTGATGGATTT
 GTCTGAGCAGTGTGGGCATGGGCTGAGGGAGCAGCTAGAACTCTTCCAGAACATCAGACCTCTCTTCATC
 AACAAAGCTCTCATAGTTGTAGCCAAACAAATGTGATGTGAAGAGAATAGCTGAACCTTTCTGAAGATGATC
 AGAAAATATTTACAGATTTGCAAGTCTGAAGGATTCCTGTAAATAGAGACCAGCACCTGACTGAGGAAGG
 TGTTATTAAGTTAAAACAGAGGCTTGCATAGGCTTTTGGCTCATCGAGTGGAACCAAAATGAAGGGA
 AATAAGTGAATGAGGTGCTGAATAGACTGCACCTGGCTATCCCAACCAGGAGGACGATAAGGAGAGGC
 CCCCTTTTATCCCTGAAGGAGTGGTGGCTCGCAGGAAGAGGATGGAACCTGAGGAGTCCAGGAAGAAGAG
 GGAACGAGATCTTGAGCTGGAATGGGAGATGATTATATTTTGGATCTTCAGAAAGTACTGGGATTTAATG
 AATTTGTCTGAAAAACATGATAAGATACCAGAAATCTGGGAAGGCCATAATATAGCTGATTATATTGATC
 CAGCCATCATGAAGAAATTGGAAGAATTAGAAAAAGAGAAGAGCTGAGAACAGCTGCTGGAGAGTATGA
 CAGTGTATCTGAGAGTGAAGACGAAGAGATGCTGGAATCCGACAGCTGGCAAAGCAAATTCGAGAGAAA
 AAGAAGTTGAAAAATCTGGAGTCCAAAGAAAAGAATACACAGGGACCCAGGATGCCGCGAACTGCTAAGA
 AGGTTCCAGAGGACAGTTTGGAGAAGGAGATGCGTAGTCTTGGTGTGACATGGACGATAAAGACGATGC
 CCATTACGAGTCCAGGCAAGAAGATCCCGGAGCATCACTAGGAAAAGAAAGCGGGAAGACTCTGCTCCC
 CCGTCTCTGTGGCCCGAGTGGGAGTTGCTCTGAACTCCACGTGACGTTTCTGGTCTTAGGGATGTCA
 AGATGGTGAAGAAAGCCAGACTATGATGAAGAATGCTCAGAAGAAGATGAATCGGTTGGGGAAGAAAGG
 GGAGGCGGATAGACACGTGTTTGATATGAAGCCCAAGCACTTGCTGTCTGGGAAGAGAAAGCTGGTAA
 AAGGACAGGAGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG207903 representing NM_012341
 Red=Cloning site Green=Tags(s)

MAHYNFKKITVVPKAKDFIDLTLSTQKRTPTVIHKHYQIHRIRHFYMRKVKFTQQNYHDRLSQILTDFF
 KLDDIHPFYADLMNILYDKDHYKLALGQINIAKNLVDNVAKYVRLMKYGDSLVRCKQLKRAALGRMCTV
 IKRQKQSLEYLEQVRQHL SRLPTIDPNTRTL LCGYPNVGKSSFINKVTRADVDPYAFATKSLFVGHM
 DYKYLRLWQVVDTPGILDHPLDRNTIEMQAITALHLRAAVLYVMDLSEQCGHGLREQLFLQNIPLFI
 NKPLIVVANKCDVKRIAE SEDDQKIFDQLQSEGFPVIETSTL TEEGVIVKTEACDRLLAHRVETKMKG
 NKVNEVLNRLHLAIPTRRDDKERPPFIPGTVARRKRMETESRKRERDLELEMDDYILDQLKYWDLM
 NLSEKHDKIPEIWEHNIADYIDPAIMKKLEELEKEELRTAAGEYDSVSESEDEMLEIRQLAKQIREK
 KKLKILESEKENTQGRMPRTAKKVQRTVLEKEMRSLGVDMDKDDAHYAVQARRSRSTIRKREDSAP
 PSSVARSGSCSRTPRDVSGLRDVKMVKKAKTMMKNAQKKMNRLGKKGEADRHVDFMKPKHLLSGKRKAGK
 KDRR

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_012341

ORF Size: 1902 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_012341.2</u> , <u>NP_036473.2</u>
RefSeq Size:	2537 bp
RefSeq ORF:	1905 bp
Locus ID:	23560
UniProt ID:	<u>Q9BZE4</u>
Cytogenetics:	10p15.3
Protein Families:	Druggable Genome
Gene Summary:	GTP-binding proteins are GTPases and function as molecular switches that can flip between two states: active, when GTP is bound, and inactive, when GDP is bound. 'Active' in this context usually means that the molecule acts as a signal to trigger other events in the cell. When an extracellular ligand binds to a G-protein-linked receptor, the receptor changes its conformation and switches on the trimeric G proteins that associate with it by causing them to eject their GDP and replace it with GTP. The switch is turned off when the G protein hydrolyzes its own bound GTP, converting it back to GDP. But before that occurs, the active protein has an opportunity to diffuse away from the receptor and deliver its message for a prolonged period to its downstream target. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG207903