

## Product datasheet for **RG208035**

### **TMEM9B (NM\_020644) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TMEM9B (NM\_020644) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** TMEM9B  
**Synonyms:** C11orf15  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG208035 representing NM\_020644  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCGACCCGTGGGGAGGCCTTCTTCGGCTTGGCTCCTTGCTCAGCCTGTCGTGCCTGGCGCTTCCG  
 TGCTGCTGCTGGCGCAGCTGTCAGACGCCGCCAAGAATTCGAGGATGTCAGATGTAATGTATCTGCC  
 TCCCTATAAGAAAATTCTGGGCATATTTATAATAAGAACATATCTCAGAAAGATTGTGATTGCCTTCAT  
 GTTGTGGAGCCCATGCCTGTGCGGGGCTGATGTAGAAGCATACTGTCTACGCTGTGAATGCAAAATATG  
 AAGAAAGAAGCTCTGTACAATCAAGGTTACCATTATAATTTATCTCCTATTTGGCCTTCTACTTCT  
 GTACATGGTATATCTTACTCTGGTTGAGCCCACTGAAGAGGCGCTCTTTGGACATGCACAGTTGATA  
 CAGAGTGATGATATTGGGGATCACCAGCCTTTTGAATGCACACGATGTGCTAGCCCGCTCCCGCA  
 GTCGAGCCAACGTGCTGAACAAGGTAGAATATGCACAGCAGCGCTGGAAGCTTCAAGTCCAAGAGCAGCG  
 AAAGTCTGTCTTTGACCGGCATGTTGTCCTCAGC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >RG208035 representing NM\_020644  
 Red=Cloning site Green=Tags(s)

MATLWGGLLRLGSLLSLSCLALSVLLLAQLSDAAKNFEDVRCKCICPPYKENSNGHIYNKNISQKDCDCLH  
 VVEPMPVRGPDVEAYCLRCECKYEERSSVTIKVTIIYLSILGLLLLVMVYLTLEPILKRRLFHQAQLI  
 QSDDDIGDHQPFANAHDVLRARSRRANVLNKVEYAQRWKLQVQEQRKSVFDRHVLS

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_020644

**ORF Size:** 594 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_020644.3](#)

RefSeq Size:	1753 bp
RefSeq ORF:	597 bp
Locus ID:	56674
UniProt ID:	<a href="#">Q9NQ34</a>
Cytogenetics:	11p15.4
Domains:	CXCXC
Protein Families:	Transmembrane
Gene Summary:	Enhances production of proinflammatory cytokines induced by TNF, IL1B, and TLR ligands. Has a role in TNF activation of both the NF-kappaB and MAPK pathways.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for RG208035