

Product datasheet for **RG208187**

MSRB3 (NM_001031679) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: MSRB3 (NM_001031679) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: MSRB3
Synonyms: DFNB74
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG208187 representing NM_001031679
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCTGCATTCAACCTGCTGCATTTGGTGACAAAGAGCCAGCCAGTAGCCCTTCGAGCCTGTGGGCTTC
 CCTCAGGGTCGTGTAGGGATAAAAAGAAGTGAAGGTGGTCTTTCCAGCAGGAACTGAGGAAGCGGCT
 AACACCCCTGCAGTACCATGTCACTCAGGAGAAAGGGACCGAAAGTGCCTTTGAAGGAGAATACACACAT
 CACAAAGATCCTGGAATATATAAATGTGTGTTTGTGGAAGTCCATTGTTAAGTCAGAAACCAATTTG
 ACTCCGGTTCAGGTTGGCCTTCATTCACGATGTGATCAATTCTGAGGCAATCACATTCACAGATGACTT
 TTCTATGGGATGCACAGGGTGGAAACAAGCTGCTCTCAGTGTGGTGCTCACCTTGGGCACATTTTGTAT
 GATGGGCCTCGTCCAAGTGGGAAAAGATACTGCATAAATTCGGCTGCCTTGTCTTTACACCTGCGGATA
 GCAGTGGCACCGCCGAGGGAGGCAGTGGGGTCGCCAGCCCGCCAGGCAGACAAAGCGGAGCTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG208187 representing NM_001031679
 Red=Cloning site Green=Tags(s)

MSAFNLLHLVTKSQPVVALRACGLPSGSCRDKKNCKVVFSSQQLRKRLTPLQYHVTQEKGTESAFEGEYTH
 HKDPGIYKCVVCGTPLFKSETKFDSGSGWPSFHDVINSEAITFTDDFSYGMHRVETSCSQCGAHLGHIFD
 DGPRPTGKRYCINSAALSFTPADSSGTAEGGSGVASPAQADKAEL

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001031679

ORF Size: 555 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001031679.3](#)

RefSeq Size: 4369 bp


RefSeq ORF: 558 bp

Locus ID: 253827

UniProt ID: [Q8IXL7](#)
 Cytogenetics: 12q14.3

Gene Summary: The protein encoded by this gene catalyzes the reduction of methionine sulfoxide to methionine. This enzyme acts as a monomer and requires zinc as a cofactor. Several transcript variants encoding two different isoforms have been found for this gene. One of the isoforms localizes to mitochondria while the other localizes to endoplasmic reticula. [provided by RefSeq, Jul 2010]

Product images:



Circular map for RG208187