

Product datasheet for **RG208345**

CCDC103 (NM_213607) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC103 (NM_213607) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CCDC103
Synonyms:	CILD17; PR46b; SMH
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG208345 representing NM_213607 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGAAGGAATGACATCATCAACTTCAAGGCTTTGGAGAAAGAGCTGCAGGCTGCACTCACTGCTGATG
AGAAGTACAAACGGGAGAATGCTGCCAAGTTACGGGCAGTGGAACAGAGGGTGGCTTCTATGAGGAGTT
CAGGGGTATTGTCCTTGATCACAATCTGAAGCCACTGGAGCGGAAGGATAAGATGGGAGGAAAGAGAACT
GTGCCCTGGAAGTGTCACTATTCCAGGAAGGACCTCCAGGATGTGGCCACTGAAATCTCCCCGAGA
AAGCCCCCTCCAGCCGAGACGTCTGCTGACTTCTATCGTGATTGGCGACGACACTTGCCAAGTGGGCC
AGAGCGCTACCAGGCTCTACTGCAGCTTGGGGTCCAAGGCTCGGCTGCCTCTCCAGACAGATGTGGGA
TTTGGACTTCTTGGGAGCTGCTGGTGGCACTGGCTGATCACGTGGGGCCGGCTGACCGGCAGCGGTGC
TGGGGATCCTATGCAGCCTGGCGAGCACTGGGCGCTTACCCTGAACCTAAGCCTGCTGAGCCGGGCAGA
GAGAGAGAGCTGCAAGGGCTTGTTCAGAAGCTGCAAGCCATGGGCAACCCAGATCCGTGAAGGAGGGG
CTCAGCTGGGAGGAGCAGGGTCTGGAGGAGCAGTCTGGTGGGCTCCAGGAGGAGGAGAGGGCTCTGCAGG
AGCTGCTGGAGCTGTACCAGTTGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG208345 representing NM_213607
Red=Cloning site Green=Tags(s)

MERNDIINFKALEKELQAALTADEKYKRENAAKLRAVEQQRVASYYEEFRGIVLASHLKPLERKDKMGGKRT
 VPWNCHTIQGRTFQDVATEISPEKAPLQPETSADFYRDWRRHLPSGPERYQALLQLGGPRLGCLFQTDVG
 FGLLGELLVALADHVGVPADRAAVLGILCSLASTGRFTLNL.SLLSRAERESCKGLFQKQLQAMGNPRSVKEG
 LSWEEQGLEEQSGGLQEEERLLQELLELYQVD

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



ACCN: NM_213607

ORF Size: 726 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_213607.3](#)

RefSeq Size: 1272 bp

RefSeq ORF: 729 bp


Locus ID: 388389

UniProt ID: [Q8IWI40](#)

Cytogenetics: 17q21.31

Gene Summary: This gene encodes a protein that contains a coiled-coil domain. [provided by RefSeq, Apr 2012]

Product images:



Circular map for RG208345