

Product datasheet for **RG208475**

OR2H1 (NM_030883) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR2H1 (NM_030883) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	OR2H1
Synonyms:	6M1-16; dj994E9.4; HS6M1-16; OLFR42A-9004-14; OLFR42A-9004.14/9026.2; OR2H6; OR2H8; OR6-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG208475 representing NM_030883 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTTAACCAAAGCTCCCCATGGGCTTCTCCTTCTGGGCTTCTCTGAACACCCAGCACTGGAAAGGA
CTCTCTTTGTGGTTGTCTTCACTTCTACCTCTTGACCCTGGTGGGCAACACACTCATCATCCTGTCTGTC
TGACTGTACCCAGGCTCCACTCTCCAATGACTTTTTCTCTGACCTCTCCTTCTGGACCTCTGC
TTTACCACAAGTTGTGTCCCCAGATGCTGGTCAACCTCTGGGGCCAAAGAAGACCATCAGCTTCTGG
GATGCTCTGTCCAGCTCTTCATCTTCTGTCCCTGGGGACCAGTGCATCCTCTGACAGTATGTC
CTTTGACCGATACGTGGCTGTCTGCCAGCCCTCCACTATGCCACCATCATCCACCCCGCTGTGCTGG
CAGCTGGCATCTGTGGCCTGGGTTATGAGTCTGGTTCAATCGATAGTCCAGACCATCCACCCCTCACT
TGCCCTTCTGTCCCACAGCAGATAGATGACTTTTTATGTGAGGTCCCATCTCTGATTCGACTCTCCTG
TGGAGATACCTCTACAATGAAATCCAGTTGGCTGTGTCCAGTGTTCATCTTCGTGGTTGTGCCTCTCAGC
CTCATCCTTGCTCTTATGGAGCCACTGCCAGGCACTGCTGAGGATTAAGTCTGCCACAGCATGGAGAA
AGGCCTTTGGGACCTGCTCCTCCATCTCACTGTGGTCAACCTCTTCTACAGCTCAGTCATTGCTGTCTA
CCTCCAGCCAAAATCCGATAGCCCAAGGAGGGCAAGTTCTTTGGTCTTTCTATGCAGTGGGCACT
CCTTCACTAACCCCTCTCGTATACCCCTGAGGAACAAGGAGATAAAGCGAGCACTCAGGAGGTTACTAG
GGAAGGAAAGAGACTCCAGGGAAGCTGGAGAGCTGCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG208475 representing NM_030883
 Red=Cloning site Green=Tags(s)

MVNQSSPMGFLLLGFSEHPALERTLFVVVFTSYLLTLVGNTLIILLSVLYPRLHSPMYFFLSDLSFLDLC
 FTTSCVPQMLVNLWGPKKTISFLGCSVQLFIFLSLGTTECILLTVMAFDRYVAVCQPLHYATIIHPRLCW
 QLASVAWVMSLVQSIIVQTPSTLHLPFCPHQQIDDFLCEVPSLIRLSCGDTSYNEIQLAVSSVIFVVVPLS
 LILASYGATAQAVLRINSATAWRKAFGTCSHLLTVVTLFYSSVIAVYLQPKNPYAQGRGKFFGLFYAVGT
 PSLNPLVYTLRNKEIKRALRRLLGKERDSRESWRAA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_030883

ORF Size: 948 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030883.4](#), [NP_112145.1](#)

RefSeq Size: 3009 bp

RefSeq ORF: 951 bp

Locus ID: 26716

UniProt ID: [Q9GZK4](#)


Cytogenetics: 6p22.1

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Olfactory transduction

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG208475