

Product datasheet for **RG208684**

RPAIN (NM_001033002) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RPAIN (NM_001033002) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	RPAIN
Synonyms:	HRIP; RIP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG208684 representing NM_001033002 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGAGTCGTTGAGGTCTCCGCGCCGCTCCCTGTACAACTGGTGGGCTCGCCGCTTGAAAGAGG
CTTCCGGCAGAGATGCCTGGAGAGAATGAGAAACAGCCGGGACAGGCTCCTAACAGGTACCGCCAGGC
TGAAGCAGTGGCCAGGAATTCTCAGAACAGCTTTCTAGTTCAAGAGGTGATGGAAGAAGAGTGAAT
GCTTTCAGTCAGTGGAGAATTGTCCAGAAGACTTGGCTCAGTTGGAGGAGCTGATAGACATGGCTGTGC
TGGAGGAAATTCAACAGGAGCTGATCAAGCAAGAGCAGTCCATCATCAGCGAGTATGAGAAGAGCTTGCA
GTTTGATGAAAAGTGTCTCAGCATCATGCTGGCTGAGTGGGAGGCAAACCACTCATCTGTCTGTATGT
ACAAAGTACAACCTGAGAATCACAAGCGGTGTGGTGGTGTGTCAGTGTGGCCTGTCCATCCCATCTCATT
CTTCTGAGTTGACAGAGCAGAAGCTTCGTGCCTGTTTAGAGGGTAGTATAAATGAGCACAGTGCACATTG
TCCCCACACCTGAATTTTCAGTCACTGGAGGAACAGAAGAAAAGTCCAGTCTTCTCATGAGCTGTCTG
GCCTGTGATACTTGGGCTGTGATCCTC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG208684 representing NM_001033002
 Red=Cloning site Green=Tags(s)

MAESLRSPRRSL YKLVGSPWPKEAFRQRCLERMNRSDRLLNRYRQAGSSGPGNSQNSFLVQEVMEEEWN
 ALQSVENCPEDLAQLEELIDMAVLEEQELIKQEQSIISEYEKSLQFDEKCLSIMLAEWANPLICPVC
 TKYNLRITSGVVVCQGLSIPSHSELTEQKLRACLEGSINEHSAHCPTPEFSVTGGTEEKSSLLMSCL
 ACDTWAVIL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001033002

ORF Size: 657 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033002.1](#), [NP_001028174.1](#)

RefSeq Size: 1536 bp

RefSeq ORF: 660 bp

Locus ID: 84268

UniProt ID: [Q86UA6](#)

Cytogenetics: 17p13.2

Gene Summary: Mediates the import of RPA complex into the nucleus, possibly via some interaction with importin beta. Isoform 2 is sumoylated and mediates the localization of RPA complex into the PML body of the nucleus, thereby participating in RPA function in DNA metabolism. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG208684