

## Product datasheet for **RG208686**

### YPEL3 (NM\_031477) Human Tagged ORF Clone

#### Product data:

Product Type: Expression Plasmids  
Product Name: YPEL3 (NM\_031477) Human Tagged ORF Clone  
Tag: TurboGFP  
Symbol: YPEL3  
Mammalian Cell Selection: Neomycin  
Vector: pCMV6-AC-GFP (PS100010)  
E. coli Selection: Ampicillin (100 ug/mL)  
ORF Nucleotide Sequence: >RG208686 representing NM\_031477  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGTGTGTGGCCAGGTCTGACAGCCACCTACTCCCTCCCGGCAGGCACTGGGCTCCCTCTGCTCCC  
CGTGGGCCGCTCCCGCGTGGGGCCACTGCCCCGCCCCCGCCATGGTGCGGATTTCAAAGCCCAAGAC  
GTTTCAGGCCTACTTGGATGATTGTCACCGGAGGTATAGCTGTGCCCACTGCCGCGCTCACCTGGCCAAC  
CACGACGACCTCATCTCCAAGTCCTCCAGGGCAGTCAGGGGCGTGCCTACCTCTCAACTCAGTGGTGA  
ACGTGGGCTGCGGGCCAGCCGAGGAGCGGGTCTGCTGACCGGCTCCATGCTGTCGCCGACATCCACTG  
CGAGAAGTGAAGACCACCTTTGGGCTGGAATATGAACAGGCCTTTGAGAGCAGCCAGAAGTACAAGAG  
GGGAAGTACATCATTGAACTCAACCACATGATCAAAGACAACGGCTGGGAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG208686 representing NM\_031477  
Red=Cloning site Green=Tags(s)

MCVAQVLTAHLPPRQALGSLCSPWAAPRVGPLPPAPAMVRISKPKTFQAYLDDCHRRYSACRAHLAN  
HDDLISKSFQGSQGRAYLFNSVNVVCGPAEERVLLTGLHAVADIHCENCKTTLGWKYEQAFESSQKYKE  
GKYIIELNHMIKDNGWD

**TRTRPLE** - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_031477

**ORF Size:** 471 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_031477.5](#)


**RefSeq Size:** 1588 bp

**RefSeq ORF:** 474 bp

**Locus ID:** 83719

UniProt ID: [P61236](#)  
 Cytogenetics: 16p11.2  
 Domains: Yippee  
 Gene Summary: Involved in proliferation and apoptosis in myeloid precursor cells.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for RG208686