

Product datasheet for **RG208904**

RGS13 (NM_002927) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RGS13 (NM_002927) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: RGS13
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG208904 representing NM_002927
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCAGGCGGAATTGTTGGATTTGTAAGATGTGCAGAGATGAATCTAAGAGGCCCCCTCAAACCTTA
CTTTGGAGGAAGTATTACAGTGGGCCAGTCTTTTGAAAATTTAATGGCTACAAAATATGGTCCAGTAGT
CTATGCAGCATATTTAAAAATGGAGCACAGTGACGAGAATTTCAATTCTGGATGGCATGTGAAACCTAT
AAGAAAATTGCCTCACGGTGGAGCAGAATTTCTAGGGCAAAGAAGCTTTATAAGATTTACATCCAGCCAC
AGTCCCTAGAGAGATTAACATTGACAGTTCGACAAGAGAGACTATCATCAGGAACATTCAGGAACCCAC
TGAAACATGTTTTGAAGAAGCTCAGAAAATAGTCTATATGCATATGGAAAGGGATTCTACCCAGATTT
CTAAAGTCAGAAATGTACCAAAAATTTGAAAATATGCAGTCCAACAACAGTTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG208904 representing NM_002927
Red=Cloning site Green=Tags(s)

MSRRNCWICKMCRDESKRPPSNLTLEEVLQWAQSFENLMATKYGPVVYAAYLKMEHSDENIQFWMACETY
KKIASRWSRISRAKKLYKIYIQQSPREINIDSSTRETIIRNIQEPTETCFEEAQKIVYMHMERDSYPRF
LKSEMYQKLLKTMQSNNSF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_002927

ORF Size: 477 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_002927.5](#)

RefSeq Size: 1498 bp

RefSeq ORF: 480 bp

Locus ID: 6003


UniProt ID: [O14921](#)

Cytogenetics: 1q31.2

Protein Families: Druggable Genome

Gene Summary: The protein encoded by this gene is a member of the regulator of G protein signaling (RGS) family. RGS family members share similarity with *S. cerevisiae* SST2 and *C. elegans* egl-10 proteins, which contain a characteristic conserved RGS domain. RGS proteins accelerate GTPase activity of G protein alpha-subunits, thereby driving G protein into their inactive GDP-bound form, thus negatively regulating G protein signaling. RGS proteins have been implicated in the fine tuning of a variety of cellular events in response to G protein-coupled receptor activation. The biological function of this gene, however, is unknown. Two transcript variants encoding the same isoform exist. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG208904