

Product datasheet for **RG208921**

ELK1 (NM_005229) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: ELK1 (NM_005229) Human Tagged ORF Clone
 Tag: TurboGFP
 Symbol: ELK1
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-AC-GFP (PS100010)
 E. coli Selection: Ampicillin (100 ug/mL)
 ORF Nucleotide Sequence: >RG208921 representing NM_005229
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGACCCATCTGTGACGCTGTGGCAGTTTCTGCTGCAGCTGCTGAGAGAGCAAGGCAATGGCCACATCA
 TCTCCTGGACTTCACGGGATGGTGGTGAATTCAGCTGGTGGATGCAGAGGAGGTGGCCCGGCTGTGGGG
 GCTACGCAAGAACAAGACCAACATGAATTACGACAAGCTCAGCCGGGCTTGGCGTACTACTATGACAAG
 AACATCATCCGCAAGGTGAGCGGCCAGAAGTTCGCTACAAGTTGTGTCCTACCTGAGGTGCGAGGGT
 GCTCCACTGAGGACTGCCGCCAGCCAGCCAGAGGTGTCTGTTACCTCCACCATGCCAAATGTGGCCCTGC
 TGCTATACATGCCGCCAGGGGACACTGTCTCTGGAAAAGCCAGGCACACCCAAGGGTGCAGGAATGGCA
 GGCCAGGCGGTTTGGCACGCAGCAGCCGGAACGAGTACATGCGCTCGGGCTCTATTCCACCTTACCA
 TCCAGTCTCTGCAGCCGAGCCACCCCTCATCCTCGGCCTGCTGTGGTCTCCCAAGTGCAGCTCCTGC
 AGGGGCAGCAGCGCCCCCTCGGGGAGCAGGAGCACCAGTCCAAGCCCCTTGGAGGCTGTCTGGAGGCT
 GAAGAGGCGGCTTGCCTCTGCAGGTATCCTGACCCCGCCGAGGCCCAAACTGAAATCGGAAGAGC
 TTAATGTGGAGCCGGTTTGGGCCGGGCTTGGCCCCAGAAGTAAAAGTAGAAGGGCCCAAGGAAGAGTT
 GGAAGTTGCGGGGAGAGAGGGTTTGTGCCAGAAACCACCAAGGCCGAGCCAGAAGTCCCTCCACAGGAG
 GCGGTGCCAGCCGGCTGCCCGGTTGTTATGGACACCCGAGGGCAGGGGGCGGCCATGCGGCTTCCA
 GCCTGAGATCTCCAGCCGAGAAAGGGCCGGAAGCCCCGGGACCTAGAGCTTCCACTCAGCCCCAGCCT
 GCTAGGTGGGCCGGGACCCGAACGGACCCAGGATCGGGAAGTGCTCCGGCCTCCAGGCTCCGGGGCCG
 GCGCTGACCCCATCCCTGCTTCTACGCATACATTGACCCCGGTGCTGCTGACACCCAGCTCGTGCCTC
 CTAGCATTCACTTCTGGAGCACCTGAGTCCCATTGCGCCCCGTAGCCCGGCCAAGCTCTCCTTCCAGTT
 TCCATCCAGTGGCAGCGCCAGGTGCACATCCCTTCTATCAGCGTGGATGGCTCTCGACCCCGTGGTG
 CTCTCCCCAGGGCCCCAGAAGCCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG208921 representing NM_005229
 Red=Cloning site Green=Tags(s)


MDPSVTLWQFLQLLREQNGHIIISWTSRDGGFEKLVDAEEVARLWGLRKNKTNMNYDKLSRALRYYYDK
 NIIRKVSQKQFVYKFSYPEVAGCSTEDCPPQPEVSVTSTMPNVAPAAIHAAPGDTVSGKPGTPKGAGMA
 GPGGGLARSSRNEYMRSGLYSTFTIQSLQPQPPHPRPAVVLPSAAPAGAAAPPSGSRSTSPSPLEACLEA
 EEAGLPLQVILTPPEAPNLKSEELNVEPGLGRALPPEVKVEGPKKELEVAGERGFVPETTKAEPEVPPQE
 GVPARLPVAVMDTAGQAGGHAASSPEISQPQKGRKPRDLELPLSPSLLGGPGPERTPGSGSGSLQAPGP
 ALTPSLLPTHHTLTPVLLTPSSLPPSIHFSTLSPIAPRSPAKLSFQFPSSGSAQVHIPSISVDGLSTPVV
 LSPGPQKP

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja2066_e10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_005229


ORF Size: 1284 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_005229.4 , NP_005220.2
RefSeq Size:	2828 bp
RefSeq ORF:	1287 bp
Locus ID:	2002
UniProt ID:	P19419
Cytogenetics:	Xp11.23
Domains:	ETS
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Endometrial cancer, ErbB signaling pathway, Focal adhesion, GnRH signaling pathway, Insulin signaling pathway, MAPK signaling pathway, Prion diseases
Gene Summary:	This gene is a member of the Ets family of transcription factors and of the ternary complex factor (TCF) subfamily. Proteins of the TCF subfamily form a ternary complex by binding to the the serum response factor and the serum response element in the promoter of the c-fos proto-oncogene. The protein encoded by this gene is a nuclear target for the ras-raf-MAPK signaling cascade. This gene produces multiple isoforms by using alternative translational start codons and by alternative splicing. Related pseudogenes have been identified on chromosomes 7 and 14. [provided by RefSeq, Mar 2012]

Product images:



Circular map for RG208921