

Product datasheet for **RG208959**

LRRC57 (NM_153260) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: LRRC57 (NM_153260) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: LRRC57
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG208959 representing NM_153260
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGAAACAGTGCCTCCGCGCTCATGTGAAACGGCGCAGAAAACGGTGTCTTCAGCTTAAGGACC
GAGGGCTGACCGAGTTCCCCGCAGACTTGCAAGCTGACGAGCAATCTCAGGACCATCGACTTGCCAA
CAACAAGATCGAAAGCCTACCGCCTTTGCTGATAGGAAAGTTCACTCTGCTGAAGAGCCTCTCCCTGAAC
AACAACTAGTACTGTTCTGCTGATGAGATATGCAATCTGAAAACTAGAGACGCTAAGCCTAAACA
ACAATCACCTTAGAGAGCTGCCGTCTACCTTTGGGCAACTCTCTGCCCTCAAGACCTGAGCCTCTCTGG
GAACCACTGGGAGCATTACCTCCCAACTTTGTAGCCTACGGCACCTGGATGTGATGGATCTCTCTAAG
AACCAGATTCGAAGTATACCTGACTCAGTGGGAGAGCTGCAAGTCATCGAACTCAACCTCAACCAGAATC
AGATATCTCAGATCTCAGTGAAGATATCTTGTGTCCACGCCTTAAATCTTCGCCTGGAAGAGAATTG
TCTTGAGCTCAGCATGCTTCCCCAGAGCATCCTCAGTGATTCCCAGATCTGTCTGCTTGTGTGGAAGGC
AATCTTTTGAATAAAGAACTTCGAGAACTGGAAGGCTATGATAAGTACATGGAGAGTTACAGCCA
CCAAGAAGAAGTTTGCG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG208959 representing NM_153260
Red=Cloning site Green=Tags(s)

MGNSALRAHVETAQKTGVFQLKDRGLTEFPADLQKLTSNLRITDLSNNKIESLPPLIGKFTLLKSLSLN
NNKLTVPDEICNLKKLETLNHNHLRELSTFGQLSALKTSLSGNQLGALPPQLCSLRHLDVMDLSK
NQIRSPDSVGLQVIELNLNQNIQISVKISCCPRLKILRLEENCLLSMLPQSILSDSQICLLAVEG
NLFEIKKLRELEGYDKYMERFTATKKFA


TRTRPLE - GFP Tag - V



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_153260

ORF Size: 717 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_153260.3](#)

RefSeq Size: 2647 bp

RefSeq ORF: 720 bp


Locus ID: 255252

UniProt ID: [Q8N9N7](#)

Cytogenetics: 15q15.2

Protein Families: Druggable Genome

Product images:



Circular map for RG208959