

Product datasheet for **RG209111**

SAMD10 (NM_080621) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SAMD10 (NM_080621) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: SAMD10
Synonyms: C20orf136; dj591C20; dj591C20.7
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG209111 representing NM_080621
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAACTACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGTTACCCGAGCTGAGGTCCAAGCTGAGCCCCCGCGTGCCGCGCCGGGGCCGTGCGCGCGGGCTTCG
GGGAGCGCCGGGATGTGGACGCCACTGCCACTTCAGTTCTGCCGACCCCTCTGGAGCACACGGTGTC
AGCTGAGAGCATCCCCTGCCACTTGCTCGGACACCTGGCACCAGCCTCACGTGGCATGACTCCCGCAGC
CAGAGGGCGGCCAGCAGCAGGCCAATCAAGCTCCTGCAGCAGCCCGGCACAGACACCCCCAGGGCCGGC
TGTACTCTGACCACTATGGCCTGTACCATAAAGCCCTCGCTGGGTGGCCTGACCCGGCCCGTGGTCT
GTGGAGTCAGCAGGACGTCTGCAAGTGGCTCAAGAAGCACTGTCCCACAACCTACCTCGTCTACGTGGAG
GCCTTCTCCCAGCATGCCATACCGGCCGGGCACTGCTGCGGCTGAATGCGGAGAAGCTGCAGCGGATGG
GGCTGGCCAGGAGGCCAGAGGAGGTGCTGCAGCAGGTGCTCCGCCTGCAGGTGCGTGAGGAGGG
GCGGAGCTGCAGCTGCTCAGCCAAGCTTCCTTCGGGAAAATGTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG209111 representing NM_080621
Red=Cloning site Green=Tags(s)


MFTELRSKLSPPRGRAGAVRAGFGERRDVDATAHFSFCRTLLEHTVSAESIPCHLPRTPGTSLTWHDSRS
 QRAASSRPIKLLQPGTDTPQGRLYSDHYGLYHTSPSLGGLTRPVVLSQQDVCKWLKKHCPHNYLVYVE
 AFSQHAITGRALLRLNAEKLQRMGLAQEAQRQEVLLQVLRQLQVREEGRSLQLLSQASFGKMS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




View online »

Cloning Scheme:


ACCN:	NM_080621
ORF Size:	606 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_080621.5
RefSeq Size:	2204 bp
RefSeq ORF:	609 bp
Locus ID:	140700
UniProt ID:	Q9BYL1
Cytogenetics:	20q13.33

Product images:



Circular map for RG209111