

## Product datasheet for **RG209348**

### Cystatin D (CST5) (NM\_001900) Human Tagged ORF Clone

#### Product data:

Product Type: Expression Plasmids  
 Product Name: Cystatin D (CST5) (NM\_001900) Human Tagged ORF Clone  
 Tag: TurboGFP  
 Symbol: Cystatin D  
 Mammalian Cell Selection: Neomycin  
 Vector: pCMV6-AC-GFP (PS100010)  
 E. coli Selection: Ampicillin (100 ug/mL)  
 ORF Nucleotide Sequence: >RG209348 representing NM\_001900  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGATGTGGCCCATGCACCCCCACTGCTGCTGCTGACTGCCTTGATGGTGGCCGTGGCCGGGAGTGCCT  
 CGGCCCAATCTAGGACCTTGGCAGGTGGCATCCATGCCACAGACCTCAATGACAAGAGTGTGCAGTGTGC  
 CCTGGACTTTGCCATCAGCGAGTACAACAAGGTCATTAATAAGGATGAGTACTACAGCCGCCCTCTGCAG  
 GTGATGGCTGCCTACCAGCAGATCGTGGGTGGGTGAAGTACTACTTCAATGTGAAGTTCGGTCGAACCA  
 CATGCACCAAGTCCCAGCCCACTTGGACAAGTGTCCCTTCAATGACCAGCCAAAAGTAAAGAGGAAGA  
 GTTCTGCTCTTTCCAGATCAATGAAGTTCCTGGGAGGATAAAATTTCCATTCTGAAGTACAAGTGCCGG  
 AAAGTC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG209348 representing NM\_001900  
 Red=Cloning site Green=Tags(s)

MMWPMHTPLLLLTALMVAVAGSASQSRFLAGGIHATDLNDKSVQCALDFAISEYNKVINKDEYYSRPLQ  
 VMAAYQQIVGGVNYFNVKFGRTTCTKSQPNLDNCPFNDQPKLKEEFCSFQINEVPWEDKISILNYKCR  
 KV


**TR**TRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

## Cloning Scheme:



ACCN: NM\_001900

ORF Size: 426 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM\\_001900.5](#)

RefSeq Size: 790 bp

RefSeq ORF:	429 bp
Locus ID:	1473
UniProt ID:	<a href="#">P28325</a>
Cytogenetics:	20p11.21
Domains:	CY
Gene Summary:	The cystatin superfamily encompasses proteins that contain multiple cystatin-like sequences. Some of the members are active cysteine protease inhibitors, while others have lost or perhaps never acquired this inhibitory activity. There are three inhibitory families in the superfamily, including the type 1 cystatins (stefins), type 2 cystatins and the kininogens. The type 2 cystatin proteins are a class of cysteine proteinase inhibitors found in a variety of human fluids and secretions. The cystatin locus on chromosome 20 contains the majority of the type 2 cystatin genes and pseudogenes. This gene is located in the cystatin locus and encodes a protein found in saliva and tears. The encoded protein may play a protective role against proteinases present in the oral cavity. [provided by RefSeq, Jul 2008]

## Product images:



Circular map for RG209348