

Product datasheet for **RG209375**

RPL36A (NM_021029) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPL36A (NM_021029) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: RPL36A
Synonyms: L36A; L44L; MIG6; RPL44
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG209375 representing NM_021029
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTTAACGTCCCTAAAACCCGCGGACTTTCTGTAAGAAGTGTGGCAAGCACCAACCCATAAAGTGA
 CACAGTACAAGAAGGGCAAGGATTCTCTGTACGCCAGGGAAAGCGGCTTATGACAGGAAGCAGAGTGG
 CTATGGTGGGCAAATAAGCCGATTTTCCGAAAAAGGCTAAAATAACAAGAAGATTGTGCTAAGGCTT
 GAGTGCCTTGAAGCAACTGCAGATCTAAGAGAATGCTGGCTATTAAGAAGATGCAAGCATTTGAAGTGG
 GAGGAGATAAGAAGAGAAAGGGCAAGTATCCAGTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG209375 representing NM_021029
 Red=Cloning site Green=Tags(s)
 MVNVPKTRRTFCKKCGKHQPHKVTQYKKGKDSL YAQGKRRYDRKQSGYGGQTKPIFRKKAKTTKKIVLRL
 ECVEPNCRSKRMLAIKRCKHFELGGDKRKGQVIQF

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_021029

ORF Size: 318 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_021029.6](#)

RefSeq Size: 469 bp
RefSeq ORF: 321 bp
Locus ID: 6173
UniProt ID: [P83881](#)
Cytogenetics: Xq22.1
Domains: Ribosomal_L44
Protein Pathways: Ribosome

Gene Summary: Cytoplasmic ribosomes, organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein, which shares sequence similarity with yeast ribosomal protein L44, belongs to the L44E (L36AE) family of ribosomal proteins. Although this gene has been referred to as ribosomal protein L44 (RPL44), its official name is ribosomal protein L36a (RPL36A). This gene and the human gene officially named ribosomal protein L36a-like (RPL36AL) encode nearly identical proteins; however, they are distinct genes. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Naturally occurring read-through transcription occurs between this locus and the heterogeneous nuclear ribonucleoprotein H2 (H') gene. [provided by RefSeq, Jan 2011]

Product images:



Circular map for RG209375