

Product datasheet for **RG209499**

ARMC9 (NM_025139) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARMC9 (NM_025139) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ARMC9
Synonyms:	ARM; JBTS30; KU-MEL-1; NS21
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG209499 representing NM_025139
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGGGACATTCTGGCTCATGAATCTGAATTACTTGGACTAGTGAAAGAGTATTTAGATTTTGCTGAAT
 TTGAAGACACCTTGAAAACATTTTCAAAAGAATGCAAAATAAAAGGAAAACCACTGTGTAAACAGTAGG
 CGGATCTTTCAGAGACTCCAATCATTGACAATTCAGAAGGATCTTGTGCTGCAATTTGACAACGGAGAC
 CAGAAGGTGTTCTTCGATCTGTGGGAGGAGCACATTTCAAGTTCCATCCGAGATGGGGACTCCTTTGCC
 AGAAGCTGGAATTCTATCTCCACATCCATTTTGCATCTATTTTTGAAGTACTCTGTGGGAGACCGGA
 CAAAGAGGAGCTGGATGAAAAGATTTCTACTTCAAAACCTACCTGGAGACCAAAGGGGAGCCTTGAGC
 CAGACCACAGAGTTTCTCCTTTCTATGCCCTTCTTTTGTCCCAACCTATGGTGCACCCCTCATTTA
 AAGAACTCTTCCAGGATTCCTGGACTCCAGAGTTAAAGTTGAAGTTGAAAAGTTTCTAGCTTTAATATC
 TAAAGCCAGCAACACGCCAAAGCTTTTAACAATATATAAGGAGAATGGACAAAGTAACAAAGAAATCTTG
 CAGCAGCTCCACCAGCAGCTGGTTGAAGCTGAACGTAGGTCACTGACATACCTCAAACGGTACAATAAGA
 TCCAGGCCGACTACCACAATCTCATTGGAGTCACAGCAGAGCTGGTGGATTCTCTAGAGGCCACAGTCAG
 CGGCAAGATGATCACCCCTGAGTACCTCCAGAGCGTCTGTGTCCGCTGTTCAAGTAAACAGATGCGGCAG
 AGCCTGGCGCATAGTGTGGACTTCACGAGGCTGGGACGGCATCCACCATGTTACGAGCCTCCTTGGCAC
 CCGTGAATTTGAAGGATGTCCATTACTGCCCTCCTTGGATTATGAGAACTGAAGAAGGATTTGATTTT
 GGGGAGTGACCGCTTGAAGCCTTCTTGTGTCAGGCTCTGCGCTGGCGCTTGACCACATCCCATCCTGGA
 GAGCAGAGGGAGACCGTTCTGCAAGCCTACATCAGCAATGACCTCTTGGACTGTTATAGCCACAACCAGA
 GGAGTGTGCTTCAGTTGCTGCACTCCACGAGCGACGTGGTGGCAGTACATGGCCAGGCTCATCAATGC
 TTTTGCCTCACTGGCAGAAGGTCGCTCTACCTTGCCCAAGACAAAGGTGCTGCAGATGCTGGAGGGA
 AGGCTGAAGGAGGAGGACAAGGATATCATCACCAGGGAGAATGTTCTTGGGGCCCTGCAGAAGTTAGTC
 TCAGGCGCCCGCTGCAGACAGCGATGATTCAAGACGGCCTCATCTTCTGGCTGGTTGATGTTCTGAAGGA
 CCCTGACTGCCTGTCTGACTACACGCTGGAGTACTCGGTGGCTTTGCTCATGAACCTCTGCCTCCGAGC
 ACAGGGAAGAACATGTGTCCAAGGTGGCAGGCTCGTGCTCAAAGTCTTTCGGATCTTCTTGGCCATG
 AAAACCATGAGATACAGCCGTATGTGAATGGAGCTCTGTACAGCATCTTCTGTTCCATCCATTCTGTA
 GGAAGCAAGAGCAATGGGAATGGAAGACATCCTACGCTGCTTCATCAAAGAAGGCAATGCTGAAATGATC
 CGCCAGATAGAATTCATCATCAAGCAGCTAAATCCGAAGAGCTACCAGATGGTGTCTTGAATCTGATG
 ATGATGAAGATGAAGATGATGAAGAGGACCATGACATCATGGAAGCCGATCTGGACAAAGACGAATGAT
 CCAGCCCCAGCTCGGAGAACTCTCAGGAGAGAAGCTTCTGACCACGGAGTACCTGGGGATCATGACCAAC
 ACGGGGAAGACAAGGCGGAAGGGCTGGCTAATGTGCAGTGGAGCGGGATGAGCCCCTGCAAAGGCCCG
 TCACCCCGCGGCCACAGAAACGGGTACCCAGTG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG209499 representing NM_025139
Red=Cloning site Green=Tags(s)


MGDILAHESSELLGLVKEYLDFAEFEDTLKTF SKECKIKGKPLCKTVGGSFRDSKSLTIQKDLVAAFNDNGD
 QKVFFDLWEEHISSSIRDGDSFAQKLEFYLIHFAIYFLKYSVGRPDKEELDEKISYFKTYLETGAALS
 QTTEFLPFYALPFVPMVHPSFKELFQDSWTPELKLEKFLALISKASNTPKLLTIYKENGQSNKEIL
 QQLHQQQLVEAERRSVTYLKRYNKIQADYHNLIGVTAELVDSLEATVSGKMITPEYLQSVCVRLFSNQMRQ
 SLAHSVDFTRPGTASTMLRASLAPVKLDVPLLP SLDYEKLLKDLILGSDRLKAFLLQALRWRLTTSHPG
 EQRETVLQAYISNDLLDCYSHNQSVLQLLHSTSDVVRQYMARLINAFASLAEGRLYLAQNKTQVLQMLEG
 RLKEEDKDIITRENVLGALQKFSLRRLQTAMIQDGLIFWLVDVLKDPDCLSDYTLEYSVALLMNLCLRS
 TGKNMCAKVAGLVKVLSDLLGHENHEIQPYVNGALYSILSVPSIREEARAMGMEDILRCFIKEGNAEMI
 RQIEFIKQLNSEELPDGVLESDDDEDEDDEEDHDIMEADLDKDELIQPQLGELSGEKLTTTEYLGIMTN
 TGKTRRGLANVQWSGDEPLQRPVTPGGHRNGYPV

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_025139

ORF Size: 1995 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_025139.2](#), [NP_079415.2](#)

RefSeq Size: 2178 bp

RefSeq ORF: 1998 bp


Locus ID: 80210

UniProt ID: [Q7Z3E5](#)

Cytogenetics: 2q37.1

Gene Summary: Required for ciliogenesis.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG209499