

Product datasheet for **RG209528**

ZCCHC8 (NM_017612) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZCCHC8
Synonyms:	PFBMFT5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RG209528 representing NM_017612
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGCGGCCGGAATTCGTGCACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCCCGCAGAGGTGTATTTGGCGATCTAGAGCTCTTCGAGCCGTTCCGACCCAGAGGAGTCGATTC
 CGAAGCCCGTTCACTCGCTTCAAGGACGACGCGGACGAGGAGGACGAAAATGGGGTCGGCAGCG
 GGAGCTACGGGAGCGGCTTCGGCAGTGCGAGGAGACCTCGAGCAGCTCCGCGCCGAGAATCAAGAACTT
 AAACGAAAATTGAACATTCGACTCGACCGAGTGAATATTGGTGAACGATACTAAGTTAGATGGACCTA
 TATTACAGATTCTATTCATGAACAATGCTATTTCAAAGCAATATCATCAAGAAATAGAGGAATTTGTATC
 AAATTTAGTAAAAAGATTTGAGGAACAGCAGAAAAATGATGTGGAAAAGACTTCTTTAATCTTTTGCCC
 CAGCCATCCAGTATTGTGCTAGAGGAGGACCACAAAGTGAAGAGTCCTGTGCCATTAACAAACAAGG
 AAGCTTTCAGTGTGTAGGAAGTGTCTGTATTTTACTAATTTTGCCTTGATAAATGGGGCAACCGCT
 TCTAAATGAAAACCTCAGCTTTCGAAGGATGGGAAATACCAAGTACCATCAAGTCTTCAGCCACATT
 GTTTCCTCTAGAAGGCAAGAAATACAAGTAAAGGCAAAAAGGCCAAAGCCTCACTGTTTCAATTGTGGTT
 CTGAAGAACACCAATGAAAGATTGCCAATGCCTCGAATGCTGCTCGAATAAGTGAAAAGAGAAAAGA
 GTATATGGATGCCTGTGGAGAAGCAAAATCAGAAATTCAGCAGCGATACCACGCAAGAAAGTAGAA
 GAAAGATTTGGAAGATTCAAGCCAGGAGTTATTAGTGAAGAACTTCAAGATGCACTAGGTGTGACAGACA
 AGAGTCTTCCACCTTTTATCTATCGGATGCGCCAGCTAGGGTACCCACAGGGTGGCTCAAAGAGGTTGA
 ATTGGAGAATTGCGGGCTTGCACTCTATGATGGAAGATGGCACTGATGGGGAACAGAAAGTTGGAGAA
 ATACAACAGAATAAAAGTGTCACTTACGATCTCTCAAATTTGGTCAACTATCCTGGTTTTAATATCTA
 CTCCCAGAGGAATTCAGACGAATGGAGGATCTTGGTTCCATACCAATGCAGGCATGTGACGAGAAGGA
 TGTGTTTGCAATTACCTTACTTCTAACTTCCAAGCGCCAGGTGTGAAGTCTGGCAACAAGAGGTCTTCA
 TCTCACTCTAGCCAGGTAGTCCAAGAAGCAGAAGAATGAAAGCAACTCAGCGGGATCTCCCGCCGACA
 TGGAGCTCGATTGAGATATGGAGGTACCACATGGTTCTCAGAGCAGCGAAAGTTTTCAGTTTCAACCACC
 ATTACCTCTGACACTCCTCCACTCCCCGGGGAACCTCTCCACCCGCTTTCACCCCTCACTCCCAAAG
 GGCACCCCGCGCTGACTCCCAAGTGAATCAACCCAGACCAAGACAGCATCTGGAGCTGTGGATGAGGACG
 CACTGACTCTAGAAGAATTGAAGAACAGCAGAGGCGGATCTGGGCAGCTCTTGAGCAGGCCGAGAGCGT
 AAACAGCGACTCCGACGTTCTGTGGACACACCTTTAACTGGCAATTCGTTGCCTCATCACCTGTCCA
 AATGAGCTAGACCTCCCTGTCCCGAGGGAAAAACATCTGAAAAGCAGACGCTGGATGAGCTGAGGTAC
 CAGAGATTTTACAAAGAAATCAGAAGCTGGACATGCCTCCAGTCCAGACTCTGAGGTGACATCACTTTG
 TCAGAAGGAAAAAGCAGAGTTGGTCCGGTAAACACTGAAGGTGCCCTTCTTGATAATGGCAGTGTCTGA
 CCAAATGTGACATCAGCAATGGGGGCGCCAGAGCTCTTCTCAGACACCAAGTCTTCAACGGCCA
 CTAATTCATAGCCCTATACCTGACATGAGCAATTTGCAACTGGAATCACGCCATTTGAATTTGAGAA
 TATGGCAGAACTCTACTGGAATGTACCTCAGGATAAGAAGCTTGTAAAGAACTCACCCGAAACAGCAG
 AAAACAAAAAGGCTCTGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG209528 representing NM_017612
 Red=Cloning site Green=Tags(s)

MAAEVYFGDLELFEFPDHPEESIPKPVHTRFKDDDGDEEDENGVGDAELRERLRQCEETIEQLRAENQEL
 KRKLNILTRPSGILVNDTKLDGPILQILFMNNAISKQYHQEIEEFVSNLVKRFEEQQKNDVEKTSFNLLP
 QPSSIVLEEDHKVEESCAIKNNKEAFSVVGSVLYFTNFCCLKLGGQPLLNNENPQLSEGWEIPKYHQVFSHI
 VSLEGQEIQVKAIRPKPHCFNCGSEEHQMKDCPMRPAARISEKRKEYMDACGEANNQNFQORYHAEVE
 ERFGRFKPGVISEELQDALGVTDKSLPPFIYMRQLGYPPGWLKEVELENSGLALYDGKDGTDGETEVE
 IQQNKSVTYDL SKLVNYPGFNISTPRGIPDEWRIFGSIQMCAQKQKDFANYLTSNFQAPGVKSGNKRSS
 SHSSPGSPKKQKNESNSAGSPADMLDSDMEVPHGSQSSESFQFPPLPDTPPLPRGTPPPVFTPLPK
 GTPPLTPSDSPQTRTASGAVDEDALTLEELQQRRIWAALAEQSVNSDSVDPVDTPLTGNVSVSSPCP
 NELDLVPPEGKTSEKQTLDEPEVPEIFTKKSEAGHASSPDSEVTSQCQKEAELAPVNTGALLDNGSVV
 PNCDISNGGSQKLFPAADTSPSTATKIHSPIPDMSKFATGITPFEFENMAESTGMYLRIRSLKNSPRNQ
 KNKKASE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_017612

ORF Size: 2121 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_017612.2](#), [NP_060082.2](#)

RefSeq Size:	2874 bp
RefSeq ORF:	2124 bp
Locus ID:	55596
UniProt ID:	Q6NZY4
Cytogenetics:	12q24.31
Domains:	PSP

Gene Summary: This gene encodes a scaffold protein which serves as an assessory factor to the nuclear RNA exosome complex. The encoded protein forms a trimeric human nuclear exosome targeting (NEXT) complex, together with hMTR4 and the RNA-binding factor RBM7 which promotes the exosomal degradation of non-coding promoter-upstream transcripts, enhancer RNAs and 3'-extended products of histone- and small nuclear RNA transcription. This complex is also thought to recruit the exosome to degrade intronic RNAs via its interaction with both the exosome and the spliceosome. It contains both an N-terminal zinc-knuckle domain and a C-terminal proline-rich domain. [provided by RefSeq, Apr 2017]

Product images:



Circular map for RG209528