

Product datasheet for **RG209589**

ZNF784 (NM_203374) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF784 (NM_203374) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ZNF784
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG209589 representing NM_203374
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCGCTGCGCGCCAGAGGCCAGAGTCGGAGCTCACCAGCTCCGGAGTCGCGATCCCAGGAGCCAC
TGGACCTGGTCTGGTGCCTGATGACTGCCGGCTGGCACACCCCGAGTGACCTCATCGAGATCCAGGT
GGTGAAGGTGACGGACACACGCTGGTCCCTGAGCCCCGGAGCCAGGTTCTTTCCACTGTGCCTTGTGC
CTTGCTGCCTTCCGGCTGGTTTCCGAGCTGCTGTTCCACGAACATGGCCACTTGGCTGGGGCCGAGGGAG
GCGGGCAGGGTGGGGACCCGAGCCGGTGCCACGTGTGCGGCCACAGCTGCCGGGGCCCGCCAGCCTGCG
CGCGCACTACAGCTTGACACGGGGAGCGGCCCTACCGCTGCGCGCTCTGCCCGCGCCTTCAAGGCC
TTGGCGCCCTGCTCCGGCACACGACCGGCACGGGTGGAGCCGGGCACCTCTCGGAGGCCTCCGGACA
CAGCGGCCGTTGCAGAACAGAGGCCGGGGTGGCCCCGGAGAGGGCGGAGGTGGTGATGGCGGCGCGGC
GGCGGGCGCAGCGGTGGGAAGCCTTTGCCTGCAGGTTCTGCGCAAGCCCTTCCGCCGCTCCTCAGAC
ATGCGAGACCACGAGCGCTGCACACTGGCGAGCGGCCATACCATTGCGGCATCTGCGGAAGGGCTTCA
CCCAGTCTCGGTGCTTAGCGGCCACGCCGCATCCACACTGGCGAGCGCCGTTCCGCTGCACGCTCTG
CGACCGCACCTTCAACAACCTCTCCAATTCCGCAAGCACCGCACCACTTCCACGGGCCGGGGCCG
GGGCTGGGAGACTCTGGAGGCCAGCTGGGCTCGTCGGCGCTGAGGGGTGCGGGAGCGGGTGTGGGGTAG
GGGACCCTGCGGAGGAGGGCGGGGGAGACCGCGAAGGTGAAGGTGGAGGCCGACCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG209589 representing NM_203374
Red=Cloning site Green=Tags(s)



MAAARPEAQRSSPTPESRSQEPLDLVLVPDDCRPGTPPSDLIEIQVVKVTDITLVPEPPEPGSFHCALC
PAAFRLVSELLFHEHGLAGAEGGGQGGDPSRCHVCGHSCPGPASLRAHYSLHTGERPYRCALCPRAFKA
LAPLLRHQHRHGVPGTSRRPPDTAAVAEQRPVAPERAEEVVMAAAAAGAAVGKPFACRFCAKPFRRSSD
MRDHERVHTGERPYHCGICGKGFTQSSVLSGHARIHTGERPFRCTLCDRTFNSSNFRKHQRTHFHGPGP
GLGDSGGQLGSSAAEGSGSGCGVGDPAEEGRGETAKVKVEADQ

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_203374

ORF Size: 969 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_203374.2](#)

RefSeq Size: 2006 bp

RefSeq ORF: 972 bp


Locus ID: 147808

UniProt ID: [Q8NCA9](#)

Cytogenetics: 19q13.42

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG209589