

Product datasheet for **RG209591**

MAMDC2 (NM_153267) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAMDC2 (NM_153267) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MAMDC2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG209591 representing NM_153267
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGTTAAGGGCGTCTCTGGCGTTGCAAGCCCTGCAGCTCGCCGGTGCCCTCGACCTGCCCGCTG
 GGTCTGTGCCTTTGAAGAGAGCACTTGC GGCTTTGACTCCGTGTTGGCCTCTCTGCCGTGATTTTAAA
 TGAGGAAGGCCATTACATTTATGTGGATACCTCCTTTGGCAAGCAGGGGGAGAAAGCTGTGCTGCTAAGT
 CCTGACTACAGGCTGAGGAATGGAGCTGCCTCCGTTTGGTCTACCAGATAACCACATCTTCGGAGTCTC
 TGTGAGATCCCAGCCAGCTGAACCTCTACATGAGATTTGAAGATGAAAGCTTTGATCGCTTGCTTTGGTC
 AGCTAAGGAACCTTCAGACAGCTGGCTCATAGCCAGCTTGGATTTGCAAAAACAGTTCCAAGAAATCAAG
 ATTTTAAATAGAAGGTGTACTAGGACAGGAAACACAGCCAGCATCGCACTATTTGAAATCAAGATGACAA
 CCGGCTACTGTATTGAATGTGACTTTGAAGAAAATCATCTCTGTGGCTTTGTGAACCGCTGGAATCCCAA
 TGTGAAGTGGTTTGTGGAGGAGGAAGTATTCGGAATGTCCACTCCATTCTCCACAGGATCACACCTTC
 AAGAGTGAAGTGGCCACTACATGTACGTGGACTCAGTTTATGTGAAGCACTTCCAGGAGGTGGCACAGC
 TCATCTCCCCGTTGACCACGGCCCCATGGCTGGCTGCCTGTCTATTTTATTACCAGATCCAGCAGGGGAA
 TGACAATGTCTTTTCCCTTTACTCTCGGATGTGGCTGGCCTTTACGAGGAAATCTGAAAGCAGACAGG
 CCAGGGAATGCTGCCTGGAACCTTGC GGAGGTCGAGTTCAATGCTCCTTACCCCATGGAGGTTATTTTGT
 AAGTTGCTTTCAATGGTCCCAAGGGAGGTTATGTTGCCCTGGATGATATTTCACTCTCCTGTTCACTG
 CCAGAATCAGACAGAACTTCTGTTCAAGTCCGTTGGAAGCCAGCTGCAATTTGAGCAAGATCTCTGCAAC
 TTTTACCAAGATAAAGAAGGTCCAGGTTGACCCGAGTAAAAGTAAAACCAACATGTATCGGGCTGGAG
 ACCACACTACAGGCTTAGGGTATTACCTGCTAGCCAACACAAAGTTCACATCTCAGCCTGGCTACATTGG
 AAGGCTCTATGGGCCCTCCCTACCAGGAACTTGCAAGTATTGTCTGCGTTTTTCATTATGCCATCTATGGA
 TTTTTAAAAATGAGTGACACCCTAGCAGTTTACATCTTTGAAGAGAACCATGTGGTTCAAGAGAAGATCT
 GGTCTGTGTTGGAGTCCCAAGGGGTGTTGGATGCAAGCTGAAATCACCTTTAAGAAGCCCATGCCTAC
 CAAGGTGGTTTTTCATGAGCCTATGCAAAAGTTTCTGGGACTGTGGGCTTGTAGCCCTGGATGACATTACA
 ATACAATTGGGAAGCTGCTCATCTTCAGAGAACTTCCACCTCCACCTGGAGAGTGTACTTTGAGCAAG
 ATGAATGTACATTTACTCAGGAGAAAAGAAACCGAGCAGCTGGCACAGGAGGAGGGGAGAACTCCAC
 TTCTACACAGGACCAAGGGAGATCACACTACTGGGTAGGCTACTACATGTACATTGAGGCTCCCAT
 ATGGTGTATGGACAAAAGCAGCCTCTTGTCCAGGCCTCTGCGAGGAGTCTCTGAAAACACTGCTTGA
 CCTTTTCTACCACATGTATGGAGGGGCACTGGCCTGCTGAGTGTATCTGAAAAGGAAGAAGACAG
 TGAAGAGTCCCTTTATGGAGGAGAAGAGTGAACAGAGCATTTCTGGCTACGAGCACTGATTGAATAC
 AGCTGTGAGAGGCAACACCAGATAATTTTTGAAGCCATTTCGAGGAGTATCAATAAGAAGTATTTGCCA
 TTGATGATGTTAAATTTAGGCAGGACCCTGTGGAGAAATGGAAGATACAACTCAACAATCATCAGGATA
 TTCTGAGGACTTAAATGAAATTGAGTAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG209591 representing NM_153267
 Red=Cloning site Green=Tags(s)


MLLRGVLALQALQLAGALDLPAGSCAFEESTCGFDSVLASLPWILNEEGHYIYVDTSF GKQGEKAVLLS
 PDLQAEWWSCLRLVYQITTSSESLSDPSQLNLYMRFEDESFDRLLSAKEPSD SWLIASLDLQNSSKFK
 ILIEGVLGQGNTASIALFEIKMTTGYCIECDFEENHLCGFVNRWNPVNVFVGGG SIRNVHSILPQDHTF
 KSELGHYMYVDSVYVKHFQEVAQLISPLTTAPMAGCLSFYYQIQQGNDNVFSLYTRDVAGLYEEIWKADR
 PGNAAWNLAEEVFNAPYPMEVIFEVAFNGPKGGYVALDDISFSPVHCQNQT ELLFSAVEASCNFEQDL CN
 FYQDKEGPGWTRVKVKNMYRAGDHTTGLGYLLANTKFTSQPGYIGRLYGPSLPGNLQYCLRFHYAIYG
 FLKMSDTLAVYIFEENHVYQEKIWSVLESPRGVWMQAEITFKKPMPTKVVFMSLCKSFWDGLVALDDIT
 IQLGSCSSSEKLPPPGECTFEQDECTFTQEKRNRSWHRRTPTSYTGPKGDHTTGVGYMYIEASH
 MUYGQKARLLSRPLRGVSGKHCLTFFYHMYGGGTGLLSVYLKKEEDSEESLLWRRRGEQSISWLRALIEY
 SCERQHQIIFEAIRGVSIRSDIAIDDKVFQAGPCGEMEDTTQSSGYS EDLNEIEY

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_153267

ORF Size: 2058 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_153267.3](#), [NP_694999.2](#)

RefSeq Size: 3610 bp
RefSeq ORF: 2061 bp
Locus ID: 256691
UniProt ID: [Q7Z304](#)
Cytogenetics: 9q21.12
Domains: MAM
Protein Families: Secreted Protein