

Product datasheet for **RG209598**

Amisyn (STXBP6) (NM_014178) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Amisyn (STXBP6) (NM_014178) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: Amisyn
Synonyms: amisyn; HSPC156
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG209598 representing NM_014178
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGTGCCAAATCTGCTATCAGCAAGGAAATTTTGCACCTCTTGATGAAAGGATGCTGGGAGCTGTCC
AAGTCAAGAGGAGGACAAAGAAAAGATTCTTTCTTGGCAACTGGAGGTCAAGGCGAATATTTAACTTA
TATCTGCCTGTCAGTGACAAACAAGAAACCCACACAGGCGTCCATCACAAGGTCAAACAGTTTGAAGGC
TCCACATCATTTGTTCCGAGATCACAGTGGATGCTCGAGCAGCTTCGCCAGGTTAATGGTATCGATCCTA
ATGGGGATTTCGGCAGAGTTTGATTGTTGTTTGAAGTCTTTGACCAGTGGGTAGCCAGCACAGCGTC
AGAAAAATGCACCTTCTCCAGATCCTCCACCATACCTGCCAGAGGTACCTCACGGACAGGAAGCCAGAG
TTTATTAAGTCCAATCCAAAATTATGGGAGGAAACAGCATCCTCCATTAGCTGCTGACAGCGTGACCA
CGCAGTGCAGAAGGCAAGCCAGGCCTTGAATGAGCGTGGAGAGCGATTAGGCCGAGCAGAGGAGAAGAC
AGAAGACCTGAAGAACAGCGCCAGCAGTTTGCAGAACTGCGCACAAAGCTTGCCATGAAGCACAAATGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG209598 representing NM_014178
Red=Cloning site Green=Tags(s)

MSAKSAISKEIFAPLDERMLGAVQVKRRTKKKIPFLATGGQGEYLTYICLSVTNKKPTQASITKVKQFEG
STSFVRRSQWMLEQLRQVNGIDPNGDSEAFDILLFENAFDQWVASTASEKTFQILHHTCQRYLTDKPE
FINCQSKIMGGNSILHSAADSVTSVAVQKASQALNERGERLGRAEEKTEDLKNSAQQFAETAHKLAMKHKC

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_014178

ORF Size: 630 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_014178.7](#)

RefSeq Size: 1970 bp

RefSeq ORF: 633 bp

Locus ID: 29091


UniProt ID: [Q8NFX7](#)

Cytogenetics: 14q12

Gene Summary:

STXBP6 binds components of the SNARE complex (see MIM 603215) and may be involved in regulating SNARE complex formation (Scales et al., 2002 [PubMed 12145319]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RG209598