

Product datasheet for **RG209841**

Cripto1 (TDGF1) (NM_003212) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cripto1 (TDGF1) (NM_003212) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: Cripto1
Synonyms: CR; CR-1; CRGF; CRIPTO
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG209841 representing NM_003212
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGACTGCAGGAAGATGGCCCGCTTCTCTTACAGTGTGATTTGGATCATGGCCATTTCTAAAGCCTTTG
 AACTGGGATTAGTTGCCGGGCTGGGCCATCAGGAATTTGCTCGTCCATCTCGGGGATACCTGGCCTTCAG
 AGATGACAGCATTGGCCCCAGGAGGAGCCTGCAATTCGGCTCGGTCTCCAGCGTGTGCCGCCATG
 GGGATACAGCACAGTAAGGAGCTAAACAGAACCTGCTGCCTGAATGGGGAACTGCATGCTGGGTCCT
 TTTGTGCCTGCCCTCCCTCCTTCTACGGACGGAAGTGTGAGCAGCATGTGCGCAAAGAGAACTGTGGGTC
 TGTGCCCATGACACCTGGCTGCCAAGAAGTGTCCCTGTGAAATGCTGGCACGGTCAGCTCCGCTGC
 TTTCTCAGGCATTTCTACCCGGCTGTGATGGCCTTGTGATGGATGAGCACCTCGTGGCTCCAGGACTC
 CAGAACTACCACCGTCTGCACGTAACCACTTTATGCTAGTTGGCATCTGCCTTTCTATACAAAGCTA
 CTAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG209841 representing NM_003212
 Red=Cloning site Green=Tags(s)

MDCRKMARFSYSYIWIWIMAIKAFELGLVAGLGHQEFARPSRGLAFRDDSIWPQEPAIRPRSSQVRVPPM
 GIQHSKELNRTCCLNNGTCMLGSFCACPPSFYGRNCEHVRKENCOSVPHDTWLPKCKSLCKCWHGQLRC
 FPQAFPLPGCDGLVMDEHLVASRTPPELPPSARTTTFMLVGICLSIQSY

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




Cloning Scheme:


ACCN:	NM_003212
ORF Size:	564 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003212.1 , NP_003203.1
RefSeq Size:	2033 bp
RefSeq ORF:	567 bp
Locus ID:	6997
UniProt ID:	P13385
Cytogenetics:	3p21.31

Protein Families: Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Stem cell - Pluripotency, Transmembrane

Gene Summary: This gene encodes an epidermal growth factor-related protein that contains a cripto, FRL-1, and cryptic domain. The encoded protein is an extracellular, membrane-bound signaling protein that plays an essential role in embryonic development and tumor growth. Mutations in this gene are associated with forebrain defects. Pseudogenes of this gene are found on chromosomes 2, 3, 6, 8, 19 and X. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]

Product images:



Circular map for RG209841