

Product datasheet for **RG209873**

ZNF770 (NM_014106) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF770 (NM_014106) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ZNF770
Synonyms:	PRO1914
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG209873 representing NM_014106
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGATGGCTGAAAACAATTTAAAAATGCTAAAGATTCAACAGTGTGTGGTAGCCAACAACTACCTAGAA
 ACAGGCCATATGTTTGAATATTTGTTTTAAGCACTTTGAAACACCATCAAAATTAGCTAGGCACTATCT
 CATTCACTAGGTCAAAAGCCATTTGAATGTGATGTGTGTCATAAACCTTTAGACAACCTAGTTCATCTG
 GAGAGGCATCAACTAAGTCTGCTCTTTAAATGTAGTATTTGTCAGCGTCACTTTAAAAATCTGA
 AGACATTTGTGAAGCACAACCTTCAATGAAACCTATCAGAATAATGTTAAACAGGTCAGAAGATT
 GCTGGAGGCCAAGCAAGAAAAGTCAATGTATGGAGTGTATAACTTTTACCACAGAGGAAAGATGGGCA
 TTACACCCGTGCTCTAAGTCTGATCCCATGTATAGCATGAAAAGAAGAAAGATATTCATGCATGTACAA
 TCTGTGGCAAGATGTTTCCATCACAGTCAAACTTGATAGGCATGTACTTATTCATACTGGTCAGAGGCC
 TTTTAAATGTGTCTTGTGTACTAAATCTTTTCGACAGTCAACTCACTTAAAAATCCCAACTTACACAT
 TCAGAAGAAAGACCTTTTCAATGTTGTTTTGTCAAAAAGGATTTAAGATTCAAAGCAAACTTCTGAAGC
 ATAAACAAATCCATACTAGGAATAAGGCTTTTCGGGCTCTTTTATTAAAGAAGAGCGGTACAGAATCTCG
 CCCCTGCCTAATAAGTTAAATGCAATCAGGGTGGTTTTGAAAATGGTGAGATTGGTGAATCTGAGGAG
 AATAATCCACTTGATGTCCACTCAATTTATATTGTCCTTTTCAATGTCCAAAGTGTGAAAAGTGTTTTG
 AATCAGAGCAGATTCTCAATGAACACAGCTGTTTTGCTGCTAGAAGTGGCAAAATTTCAAGCAGGTTCAA
 AAGAAGCTACAATAAAACCATTGTTAAAAAATCTTGCCAAAGCTTAAGCGTGCTAGGAGTAAAAAA
 TTAGATAACTTTCAATCTGAGAAAAAGTTTTTAAAAAGAGTTTCTTGAGAAATTGTGATCTTATTTCTG
 GTGAGCAGAGCTCTGAACAAACCCAGAGAACATTTGTGGGTTCTCTTGCCAAACATGGAACATATAAAC
 AATTGGCAATAGAAAGAAGAAAACATTGACTTTGCCATTTTCTTGCCAAATATGGGAAAAAATTTGAAA
 GGCATCCTTACGACAGAAAACATATTAAGCATTGATAATTCAGTGAATAAGAAAGACTTGTCAATCTGTG
 GTTCATCAGGTGAGGAATCTTTAATAACTGTGAGGTACTTCAGTGTGGTTTTTCAGTTCCAAGGGAAAA
 CATACGTACTAGACATAAGATATGTCCTTGTGACAAATGTGAGAAGGTATTTCTTCTATATCCAACTA
 AAAAGACACTATTTAATTCATACTGGACAGAGGCCCTTTGGCTGTAATTTTGTGGGAAATCTTTAGAC
 AGTCAGCTCACTTAAAAAGACATGAACAGACTCATAATGAAAAGAGTCTTATGCATCTCTTTGCCAAGT
 AGAATTTGGAACTTCAACAATCTTTCTAATCATTGAGTAATAATGTTAACTATAATGCTTCCCAACAA
 TGTGAGGCTCCTGGTGTCAAAAATACGAGGTCTCAGAGTCAGATCAAATGTCAGGAGTTAAGGCAGAGT
 CACAGGATTTTATTCCTGGTAGCACCGGCAACCCTGTCTTCTAATGTACTTTTGAATCAGAGCAAAAG
 CAATCCTTTTTGCAGTTATTCAGAGCATCAGGAGAAAAATGATGTCTTCTGTACCGATGCAGTGTTTGT
 GCTAAAAGTTTCCGATCTCCATCTAAACTGGAAAGACACTACCTAATTCATGCAGGGCAGAAACCATTGT
 AATGCTCAGTTTGTGGCAAAACATTGACAGGCTCCTCACTGGAAGAGACATCAGCTTACTCACTTTAA
 AGAACGACCACAAGGAAAGTGTTGCCTTAGATTTCGTTATG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG209873 representing NM_014106
Red=Cloning site Green=Tags(s)


MMAENNLKMLKIQCCVANKLPNRNPYVCNICKHFETPSKLARHYLIHTGQKPFECDEVCHKTFRQLVHL
 ERHQLTHSLPFKCSICQRHFKNLKTFVKHQQLHNETYQNNVKQVRRLLLEAKQEKSMYGVYNTFTTEERWA
 LHPCSKSDPMYSMKRRKNIHACTICGKMFPSSQSKLDRHVL IHTGQRPFKCVLCTKSFRQSTHLKIHQLTH
 SEERPFQCCFCQKGFQKIQSKLLKHKQIHTRNKAFRALLLKKRRTESRPLPNKLNANQGGFENGEGESEE
 NNPLDVHSIYIVPFQCPKCEKCFESEQILNEHSCFAARSGKIPSRFKRSYNYKTIKKILAKLKRARSKK
 LDNFQSEKKVFKKSFLRNCDLISGEQSSEQTQRTFVGSGLKGHTYKTIGNRKKKTLTLPSWQNMGNLNLK
 GILTTENILSIDSVNKKDLSICGSSGEEFFNCEVLQCGFSVPRENIRTRHKICPCDKCEKVFPISKL
 KRHYLIHTGQRPFGCNICGKSFRQSAHLKRHEQTHNEKSPYASLCQVEFGNFNLSNHSGNNVNYNASQQ
 CQAPGVQKYEVSQDQMSGVKAESQDFIPGSTGQPCLPNVLLESEQSNPFCSEYSEHQEKNDVFLYRCSVC
 AKSFRSPSKLERHYLIHAGQKPFECVCGKTRQAPHWKRHQLTHFKERPQGGKVALDSVM

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_014106

ORF Size: 2073 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_014106.2](#), [NP_054825.2](#)

RefSeq Size: 5418 bp

RefSeq ORF: 2076 bp


Locus ID: 54989

UniProt ID: [Q6IQ21](#)

Cytogenetics: 15q14

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG209873