

Product datasheet for **RG209932**

Cezanne (OTUD7B) (NM_020205) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cezanne (OTUD7B) (NM_020205) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cezanne
Synonyms:	CEZANNE; ZA20D1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG209932 representing NM_020205
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGACCCTGGACATGGATGCTGTTCTGTACAGATTTTGTCCGTTCCACAGGAGCAGAGCCAGGGCTAGCGC
 GAGATCTCCTAGAAAGAAAGAATTGGGATGTGAATGCCGCCCTCAGTGATTTTGAACAGCTACGTCAAGT
 CCATGCTGGAAACCTACCCCATCCTTTAGTGAGGGGAGTGTTGGCTCCAGGACCCTGAAAAAGGGTTT
 TCTGACAGAGAGCCTACTCGCCCTCCCCAGCCATCCTCCAGCGCAGGATGACATCGTTCAAGAAAAAC
 GCCTGTCTAGGGGCATCTCCACGCCAGCTCCAGCATTGTTTCCCTGGCCCGTCCCATGTCTCCTCCAA
 TGGTGGGGTGGGGGAGCAATGAGCACCCCTGGAAATGCCATCTGTGCCTTCCAGTTCAGATCTC
 ACTGTATACAATGAAGACTTCCGCAGTTCATAGAGAGAGACCTATTGAGCAGTCCATGCTGGTTGCCT
 TGGAACAGGCAGGGCGTTTGAAGTGGTGGGTGAGTGTGGATCCCACCTCTCAGAGGCTGCTTCTTTGGC
 AACTACTGGAGATGGGAACTGCCTCCTGCATGCAGCCTCCCTTGAATGTGGGGTTTCCATGATCGGGAC
 TTGATGCTGCGGAAAGCTTTGTATGCACCTGATGGAGAAGGGAGTTGAGAAGGAAGCGTTGAAAAGGCCT
 GGAGGTGGCAGCAGACACAGCAGAATAAAGAGTCAGGGCTGGTATACACAGAAGATGAATGGCAGAAAGGA
 GTGGAATGAAGTATCAAGCTTGCCTCAAGTGAACCCGAATGCATCTAGGTACCAATGGAGCCAACTGT
 GGTGGGGTGGAGAGTTCTGAGGAGCCTGTATATGAGAGCCTTGAAGAGTTTACGCTTTTGTCTTGGCTC
 ATGTGCTTAGGAGGCCATAGTCGTCGTGGCAGACACCATGCTGAGGGACTCCGGAGGGGAAGCATTGTC
 CCCTATCCCTTTGGAGGAATCTATCTGCCTTTGGAGTCCCAGCCAGCCAGTGTACCCTCCCTCTG
 AAGTGTGATCCCACTTACAGATTCAGAGTAAAGCTGCTGCCCTTGCCTTGTGCTGAGAGGAGGAAACCAAGGAAC
 GGGCTGGGAGTGGGGCAAAGATGATAGTACAATGTCCGATTGGCCAGTGAATTCCTGTCCCTAGAGGTC
 AAATTGCATCTGCTGCATAGCTACATGAATGTGAAGTGGATCCCCTGTCTCTGATGCACAGGCTCCTC
 TGGCCAGCCTGAGTCCCCACCGCCTCAGCTGGAGATGAGCCCCGGTCCACTCCTGAGTCTGGAGACTC
 AGACAAGGAGTCAAGTGGCAGCAGTTCACCAGCAACGAGGGCGGCCGGAAGGAGAAGTCAAAGCGA
 GATCGGGAGAAGGACAAGAAGAGAGCAGATTCTGTGGCTAACAACTGGGCAGCTTTGGCAAAACCTTGG
 GCAGCAAGCTCAAGAAGAACATGGGGGGCTGATGCACAGCAAGGGTTCAAAGCCTGGAGGGTGGGGAC
 AGGTTGGGAGGAAGCAGCGCACTGAGACACTGGAGAAGAAGAAGAAAACTACTGAAGAGCTGGAAG
 GGTGGCAAGGAGGAGGCAGCTGGGGATGGGCCTGTGTCTGAGAAGCCCCAGCTGAGTCTGTTGGTAACG
 GAGGGAGCAAGTATAGCCAGGAGGTGATGCAGAGCCTGAGCATTCTGAGGACTGCCATGCAAGGGGAGGG
 GAAGTTTATTTTTGTTGGAACCTGAAGATGGGTACCCTCACCAGTATCAGGAGGAAATGATCCAGCGC
 TACCTTTCTGATGCTGAGGAGAGATTCTGGCAGAACAGAAGCAGAAGGAGGCAGAGAGGAAGATCATGA
 ATGGAGGAATAGGGGGTGGCCCTCCTCCAGCAAAAAGCCAGAGCCAGATGCTAGGGAAGAGCAGCCGAC
 CGGTCCCCCAGCAGAGTCCAGGGCAATGGCATTTCCTACTGGCTACCCTGGGACTTTACTATCCCTCGG
 CCGTCTGGGGCGGAGTCCACTGCCAGGAACCCCGAGGCAGTTGGCAGGGGTCCATGTGTCGGGGGCC
 TACCACCATATGCCACCTTCCCAGACAGTGCCTCCTGGGCGACCCTACCCCAACAGGACAGCATCCC
 TTCTCTGGAGCCAGGCAGCCACTAAGGATGGACTTACAGGGGTGCCTTGTACCACCCCTACCGA
 GTGGCTGATTCTATAGCAATGGCTACAGAGAGCCCCCTGAGCCAGATGGATGGGCTGGAGGTCTCCGGG
 GCCTTCCCCAACTCAGACCAAATGCAAAACACCGAACTGCAGCTTCTATGGACCCCTGAGACAAACAA
 CTTCTGTTCTGTTGTTACAGGGAAGAAGTGGAGGAGGGAGCGGAACCGGATGGGGAGCTCCTGGT
 CACAGTTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG209932 representing NM_020205
Red=Cloning site Green=Tags(s)

```
MTLDMDAVLSDFVVRSTGAEPGLARDLLEGGKNDVNAALSDFEQLRQVHAGNLPPSFSEGGSSRTPEKGF
SDREPTRPPRPILQRQDDIVQEKRLSRGISHASSIVSLARSHVSSNGGGGGSNEHPLEMPICAFQLPDL
TVYNEDFRSFIERDLIEQSMLVALEQAGRLNWWWVSDPTSQRLPLATTGDGNCLLHAASLGMWGFHDRD
LMLRKALYALMEKGVEKEALKRRWRWQQTQQNKESGLVYTEDEWQKEWNELIKLASSEPRMHLGTNGANC
GGVESSEEPVYESLEEFHVFVLAHVLRPPIVVVADTMLRDSGGEAFAPIPFGGIYLPLEVPASQCHRSP
VLAYDQAHFSALVSMEQKENTKEQAVIPLTDSEYKLLPLHFAVDPGKGWEWGKDDSDNVRLASVILSLEV
KLHLLHSYMNKWIPLSSDAQAPLAQPESPTASAGDEPRSTPESGSDKESVGSSTSNEGGRRKEKSKR
DREKDKKRADSVANKLGSFGKTLGSKLKKNMGGLMHSKGSKPGGVTGLGGSSGTETLEKKKNSLKS WK
GGKEEAAGDGPVSEKPPAESVNGGSKYSQEVMSLSILRTAMQGEKGFIFVGTLMGHRHQYQEEMIQR
YLSDAEERFLAEQKQKEAERKIMNGGIGGGPPAKKPEPDAREEQPTGPPAESRAMAFSTGYPGDFTIPR
PSGGGVHCQEP RRQLAGGPCVGG LPPYATFPRQCPPGRPYPHQDSIPSLEPGSHSKDGLHRGALLPPPYR
VADSYNGYREPPPEPDGWAGGLRGLPPTQTKCKQPNC SFYGH PETNNFCSCCYREELRRRERE PDGELLV
HRF
```

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-Mlul

Cloning Scheme:



ACCN: NM_020205

ORF Size: 2529 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_020205.4
RefSeq Size:	6413 bp
RefSeq ORF:	2532 bp
Locus ID:	56957
UniProt ID:	Q6GQQ9
Cytogenetics:	1q21.2
Gene Summary:	<p>Negative regulator of the non-canonical NF-kappa-B pathway that acts by mediating deubiquitination of TRAF3, an inhibitor of the NF-kappa-B pathway, thereby acting as a negative regulator of B-cell responses. In response to non-canonical NF-kappa-B stimuli, deubiquitinates 'Lys-48'-linked polyubiquitin chains of TRAF3, preventing TRAF3 proteolysis and over-activation of non-canonical NF-kappa-B. Negatively regulates mucosal immunity against infections (By similarity). Deubiquitinates ZAP70, and thereby regulates T cell receptor (TCR) signaling that leads to the activation of NF-kappa-B (PubMed:26903241). Plays a role in T cell homeostasis and is required for normal T cell responses, including production of IFNG and IL2 (By similarity). Mediates deubiquitination of EGFR (PubMed:22179831). Has deubiquitinating activity toward 'Lys-11', 'Lys-48' and 'Lys-63'-linked polyubiquitin chains (PubMed:27732584). Has a much higher catalytic rate with 'Lys-11'-linked polyubiquitin chains (in vitro); however the physiological significance of these data are unsure (PubMed:27732584). Hydrolyzes both linear and branched forms of polyubiquitin.[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for RG209932